



PRODUCT
CATALOG
VOLUME 13



BUSA[®]
SURGICAL



BUSA[®]
SURGICAL

How to order

For additional assistance contact Product Technical Support:
productinfo@BUSAMedical.com

Return Goods Policy

Contact your distributor regarding returned goods policy.

Warranty Policy

Contact your distributor for details regarding warranty.

Repairs

Contact your distributor for details regarding repairs.

Technical Support

For additional assistance contact Product Technical Support:
productinfo@BUSAMedical.com

Table of Contents

 BUSA®	 ConMed® Linvatec®	 Stryker®
 Zimmer Biomet®	 DePuy Synthes®	 Codman®, Sarns® & Medtronic® Midas Rex®


Twist Drills - Standard, Round Shank		
	BUSA® Twist Drills - Standard, Round Shank	7
Quick Connect Twist Drills - Solid & Cannulated		
	Compatible with DePuy Synthes® Quick Connect Coupling	8-10
K-Wires - Smooth/Plain, Extra Sharp, Fully Threaded		
	BUSA® K-Wires - Smooth/Plain, Extra Sharp, Fully Threaded	11-12
Steinmann Pins - Smooth/Plain, 12 Inch, Fully Threaded		
	BUSA® Steinmann Pins - Smooth/Plain, 12 Inch, Fully Threaded	13-14
Orthopaedic Pin/Drill Packs		
	Compatible with Stryker® Knee Instrumentation Systems	15
	Compatible with Stryker® Ortholock™ Navigation System	15
 	Compatible with DePuy Synthes® Knee Instrumentation Systems	16
	Compatible with Zimmer Biomet® Knee Instrumentation Systems	16
Guide Wires		
	Compatible with DePuy Synthes® 3.5mm & 4.0mm Cannulated Screws	17
Large Bone Saw Blades		
	Compatible with Stryker® Large Bone Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD	18-49

Table of Contents

 BUS A®	 ConMed® Linvatec®	 Stryker®
 Zimmer Biomet®	 DePuy Synthes®	 Codman®, Sarns® & Medtronic® Midas Rex®









Large Bone Saw Blades		
	Compatible with Stryker® Large Bone Reciprocating Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower	50-51
	Compatible with ConMed® Linvatec® Large Bone Oscillating Saw Systems: HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus	52-60
	Compatible with ConMed® Linvatec® Large Bone Oscillating Saw Systems: HALL® Titan®, VersiPower®, Series 4, Series 3 & Orthair	61
	Compatible with ConMed® Linvatec® Large Bone Reciprocating Saw Systems: HALL 50™, HALL® Titan®, MPower®, PowerPro®, VersiPower® Plus, VersiPower®, Series 4, Series 3 & Orthair; MicroAire® & Smith® & Nephew® Large Bone Reciprocating Saw Systems	62
	Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems: Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System (TRS) Sagittal Saw Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal	63-69
	Compatible with DePuy Synthes® Large Bone Reciprocating Saw Systems: Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System	70
	Compatible with Zimmer Biomet® Reciprocating Saw System: Universal	71
Small Bone Blades & Rasps		
	Compatible with Stryker® Oscillating & Sagittal Saw Systems: SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2	72-79

Table of Contents

 BUS A®	 ConMed® Linvatec®	 Stryker®
 Zimmer Biomet®	 DePuy Synthes®	 Codman®, Sarns® & Medtronic® Midas Rex®

















Small Bone Blades & Rasps		
	Compatible with Stryker® Reciprocating Saw Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2	80-81
	Compatible with Stryker® Small Bone Reciprocating Saw System: Micro	82-83
	Compatible with ConMed® Linvatec® Small Bone/Micro Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUS A®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems	84-86
	Compatible with ConMed® Linvatec® Small Bone/Micro Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUS A®, MicroAire® & Smith & Nephew® Small Bone/Micro Oscillating Saw Systems	87-89
	Compatible with ConMed® Linvatec® Small Bone/Micro Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUS A®, MicroAire® & Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems	90-91
	Compatible with ConMed® Linvatec® Small Bone/Micro Mini-Driver™ K100 & K200 Sagittal Saw Systems: K120 Attachment	92-93
	Compatible with ConMed® Linvatec® Small Bone/Micro Mini-Driver™ K100 & K200 Sagittal Saw Systems: K220 Attachment	94
Sternum Saw Blades		
	Compatible with Stryker® Sternum Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower	95-96
	Compatible with Stryker® Large Bone Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD	97
	Compatible with ConMed® Linvatec® Sternum Saw Systems	98

Table of Contents





 BUSA®	 ConMed® Linvatec®	 Stryker®
 Zimmer Biomet®	 DePuy Synthes®	 Codman®, Sarns® & Medtronic® Midas Rex®

Sternum Saw Blades		
	Compatible with ConMed® Linvatec® Large Bone Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair; ConMed® Linvatec® Small Bone/Micro Sagittal Saw System: Micro 100™; ConMed® Linvatec® Repeat Sternotomy Saw System; BUSA®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems	99
	Compatible with ConMed® Linvatec® Small Bone/Micro Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems	100
	Compatible with Codman® Sternum Saw Systems	101
	Compatible with Sarns® Sternal Saw II	102
Telescoping Burs		
	Compatible with Stryker® Signature Portfolio Elite Attachments	103-110
High Speed Burs		
	Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™; ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II	111-114
Micro Burs - Standard & Oral Surgery		
	Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2	115-124
	Compatible with Stryker® Impaction Drill Systems: 100K, CMD, CORE™ & Oral Max	125-126
	Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® - Small Bone/Micro Drill Systems	127-137
Cross Reference		138-175

ConMed®, Linvatec®, HALL®, HALL 50™, Micro100™, Micro E®, MicroChoice®, MicroPower®, Mini-Driver™, MPower®, PowerPro®, Surgairtome Two®, Titan®, and VersiPower® are trademarks of CONMED Corporation or its affiliates.
 DePuy Synthes® and Codman® are registered trademarks of DePuy Synthes, Inc., a Johnson & Johnson Company, or its affiliates.
 Medtronic® and Midas Rex™ are trademarks of Medtronic, Inc., or its affiliates.
 MicroAire® is a registered trademark of MicroAire Surgical Instruments, LLC, or its affiliates.
 Sarns® is a registered trademark of Sarns, Inc., or its affiliates.
 Stryker®, Command®, CORE™, Ortholock®, RemB®, TPS™, and TPX® are trademarks of Stryker Corporation or its affiliates.
 Surgairtome Two® is a registered trademark of American Sterilizer Company
 Zimmer Biomet® is a registered trademark of Zimmer Biomet, Inc., or its affiliates.
 BUSA® is a registered trademark of Peter Brasseler Holdings, LLC, or its affiliates.

BUSA® Twist Drills - Standard - Round Shank

Fabricated from autoclavable stainless steel.

BUSA®				
Order Number		Length	Diameter	
Non-Sterile (6/pk)	Sterile (1/pk)	in	mm	in
KM166-00-11	KM166-10-11	4	1.1	3/64
				
KM166-00-00	KM166-10-00	5	1.6	1/16
				
KM166-00-01	KM166-10-01	5	2.0	5/64
				
KM166-00-02	KM166-10-02	5	2.4	3/32
				
KM166-00-03	KM166-10-03	5	2.8	7/64
				
KM166-00-04	KM166-10-04	5	3.2	1/8
				
KM166-00-05	KM166-10-05	5	3.6	9/64
				
KM166-00-06	KM166-10-06	5	4.0	5/32
				
KM166-00-45	KM166-10-45	5	4.5	11/64
				
KM166-00-07	KM166-10-07	5	4.8	3/16
				
KM166-00-60	KM166-10-60	5	6.0	15/64
				
KM166-00-08	KM166-10-08	5	6.4	1/4
				

Sterilization Parameters


Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.

BUSA® Quick Connect Twist Drills - Solid Replacements - DePuy Synthes®

Compatible with DePuy Synthes® Quick Connect Coupling

Fabricated from autoclavable stainless steel. Features one piece construction.





DePuy Synthes®							
Order Number		Replacement		Diameter (mm)	Total Length (mm)	Effective Length (mm)	Number of Flutes
Non-Sterile (1/pk)	Sterile* (1/pk)	Non-Sterile	Sterile				
KM166-310-11	KM166-310-11S	310.11	310.11S	1.1	60	35	2
							
KM166-310-15	KM166-310-15S	310.15	310.15S	1.5	87	61	2
							
KM166-310-16	KM166-310-16S	310.16	310.16S	1.5	111	86	2
							
KM166-310-19	KM166-310-19S	310.19	310.19S	2.0	100	75	2
							
KM166-310-21	KM166-310-21S	310.21	310.21S	2.0	125	100	2
							
KM166-310-23	KM166-310-23S	310.23	310.23S	2.5	181	156	2
							
KM166-310-25	KM166-310-25S	310.25	310.25S	2.5	111	86	2
 TiN - Titanium Nitride Coating							
KM166-310-26	KM166-310-26S	310.26	310.26S	2.7	102	76	2
 TiN - Titanium Nitride Coating							
KM166-310-28	KM166-310-28S	310.28	310.28S	2.7	126	100	2
							
KM166-310-29	KM166-310-29S	310.29	310.29S	3.2	196	170	2
							
KM166-310-31	KM166-310-31S	310.31	310.31S	3.2	147	121	2
							
KM166-310-35	KM166-310-35S	310.35	310.35S	3.5	112	86	2
							

BUSA® Quick Connect Twist Drills - Solid Replacements - DePuy Synthes®



Compatible with DePuy Synthes® Quick Connect Coupling

Fabricated from autoclavable stainless steel. Features one piece construction.

DePuy Synthes®							
Order Number		Replacement		Diameter (mm)	Total Length (mm)	Effective Length (mm)	Number of Flutes
Non-Sterile (1/pk)	Sterile* (1/pk)	Non-Sterile	Sterile				
KM166-310-37	KM166-310-37S	310.37	310.37S	3.5	196	171	2
							
KM166-310-44	KM166-310-44S	310.44	310.44S	4.5	147	121	2
							
KM166-310-48	KM166-310-48S	310.48	310.48S	4.5	195	170	2
							
KM166-310-288	N/A	310.288	N/A	2.8	165	140	2
							
KM166-310-509	N/A	310.509	N/A	1.8	110	85	2
							
KM166-310-530	KM166-310-530S	310.530	310.530S	2.4	100	75	2
							
KM166-315-28	N/A	315.28	N/A	2.7	127	102	3
							
KM166-323-062	N/A	323.62	N/A	2.0	140	115	3
							
KM166-324-214	N/A	324.214	N/A	2.8	197	171	3
							
KM166-357-407	N/A	357.407	N/A	4.0	260	235	3
							

Sterilization Parameters





Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.


BUSA® Quick Connect Twist Drills - Cannulated Replacements - DePuy Synthes®

Compatible with DePuy Synthes® Quick Connect Coupling

Fabricated from autoclavable stainless steel. Features one piece construction.

DePuy Synthes®						
Order Number	Replaces (Non-Sterile)	Diameter (mm)	Cannulation (mm)	Total Length (mm)	Effective Length (mm)	Number of Flutes
Non-Sterile (1/pk)						
KM166-310-65	310.65	3.2	1.7	170	139	3
						
KM166-310-66	310.66	4.5	1.7	170	139	3
						
KM166-310-67	310.67	2.7	1.4	160	130	3
						
KM166-310-68	310.68	3.5	1.4	160	130	3
						
KM166-310-221	310.221	2.0	1.2	150	50	3
						

Note: Because of their hollow cross section, cannulated drill bits are more susceptible to breakage than solid drill bits of the same diameter. Use less axial force, avoid bending and advance the drill bit slowly to minimize the possibility of breakage.




Sterilization Parameters

Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.

BUSA® K-Wires - Smooth, Plain & Extra Sharp

Fabricated from precision ground 316L implant grade stainless steel (ASTM F 138).



Items not shown at actual size.

Diameter		Order Number								
		4" Length		5" Length		5.5" Length	6" Length		9" Length	
mm	in	Non-Sterile (6/pk)	Sterile (1/pk)	Non-Sterile (6/pk)	Sterile (1/pk)	Non-Sterile (6/pk)	Non-Sterile (6/pk)	Sterile (1/pk)	Non-Sterile (6/pk)	Sterile (1/pk)
Point Style 1. ST										
0.7	.028	KM172-14-28	KM172-14-28S	N/A	N/A	KM172-155-28	KM172-16-28	KM172-16-28S	N/A	N/A
0.9	.035	KM172-14-35	KM172-14-35S	KM172-15-35	N/A	N/A	KM172-16-35	KM172-16-35S	KM172-19-35	KM172-19-35S
1.1	.045	KM172-14-45	KM172-14-45S	KM172-15-45	N/A	N/A	KM172-16-45	KM172-16-45S	KM172-19-45	KM172-19-45S
1.4	.054	KM172-14-54	KM172-14-54S	KM172-15-54	N/A	KM172-155-54	KM172-16-54	KM172-16-54S	KM172-19-54	KM172-19-54S
1.6	.062	KM172-14-62	KM172-14-62S	KM172-15-62	N/A	KM172-155-62	KM172-16-62	KM172-16-62S	KM172-19-62	KM172-19-62S
Point Style 2. DT										
0.7	.028	KM172-24-28	KM172-24-28S	N/A	N/A	N/A	KM172-26-28	KM172-26-28S	KM172-29-28	KM172-29-28S
0.9	.035	KM172-24-35	KM172-24-35S	N/A	N/A	N/A	KM172-26-35	KM172-26-35S	KM172-29-35	KM172-29-35S
1.1	.045	KM172-24-45	KM172-24-45S	N/A	N/A	N/A	KM172-26-45	KM172-26-45S	KM172-29-45	KM172-29-45S
1.4	.054	KM172-24-54	KM172-24-54S	N/A	N/A	N/A	KM172-26-54	KM172-26-54S	KM172-29-54	KM172-29-54S
1.6	.062	KM172-24-62	KM172-24-62S	N/A	N/A	N/A	KM172-26-62	KM172-26-62S	KM172-29-62	KM172-29-62S
Point Style 3. SD										
0.7	.028	N/A	N/A	N/A	N/A	N/A	KM172-36-28	N/A	KM172-39-28	N/A
0.9	.035	N/A	N/A	N/A	N/A	N/A	KM172-36-35	N/A	KM172-39-35	N/A
1.1	.045	N/A	N/A	N/A	N/A	N/A	KM172-36-45	N/A	KM172-39-45	N/A
1.4	.054	N/A	N/A	N/A	N/A	N/A	N/A	N/A	KM172-39-54	N/A
1.6	.062	N/A	N/A	N/A	N/A	N/A	KM172-36-62	N/A	KM172-39-62	N/A
Point Style 4. DD										
0.7	.028	KM172-44-28	KM172-44-28S	KM172-45-28	KM172-45-28S	N/A	KM172-46-28	KM172-46-28S	KM172-49-28	KM172-49-28S
0.9	.035	KM172-44-35	KM172-44-35S	KM172-45-35	KM172-45-25S	N/A	KM172-46-35	KM172-46-35S	KM172-49-35	KM172-49-35S
1.1	.045	KM172-44-45	KM172-44-45S	KM172-45-45	KM172-45-45S	N/A	KM172-46-45	KM172-46-45S	KM172-49-45	KM172-49-45S
1.4	.054	KM172-44-54	KM172-44-54S	N/A	N/A	N/A	KM172-46-54	KM172-46-54S	KM172-49-54	KM172-49-54S
1.6	.062	KM172-44-62	KM172-44-62S	KM172-45-62	KM172-45-62S	N/A	KM172-46-62	KM172-46-62S	KM172-49-62	KM172-49-62S
Point Style 5. DD Extra Sharp										
0.7	.028	N/A	KM172-44-28ESS	N/A	N/A	N/A	N/A	N/A	N/A	N/A
0.9	.035	N/A	KM172-44-35ESS	N/A	N/A	N/A	N/A	N/A	N/A	N/A
1.1	.045	N/A	KM172-44-45ESS	N/A	N/A	N/A	N/A	KM172-46-45ESS	N/A	N/A
1.6	.062	N/A	KM172-44-62ESS	N/A	N/A	N/A	N/A	KM172-46-62ESS	N/A	N/A

BUSA® K-Wires - Fully Threaded

Fabricated from precision ground 316L implant grade stainless steel (ASTM F 138).



1. Single Trocar
3-Sided Point; Point on One End



2. Double Trocar
3-Sided Point; Point on Both Ends



3. Single Diamond
2-Sided Point; Point on One End



4. Double Diamond
2-Sided Point; Point on Both Ends



Items not shown at actual size.

BUSA®				
Point Style	Diameter		Order Number	
			9" Length	
	mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
1. ST	0.9	.035	KM173-19-35	KM173-19-35S
	1.1	.045	KM173-19-45	KM173-19-45S
	1.4	.054	N/A	N/A
	1.6	.062	KM173-19-62	KM173-19-62S
2. DT	0.9	.035	KM173-29-35	KM173-29-35S
	1.1	.045	KM173-29-45	KM173-29-45S
	1.4	.054	KM173-29-54	KM173-29-54S
	1.6	.062	KM173-29-62	KM173-29-62S
3. SD	0.9	.035	KM173-39-35	KM173-39-35S
	1.1	.045	KM173-39-45	KM173-39-45S
	1.4	.054	KM173-39-54	KM173-39-54S
	1.6	.062	KM173-39-62	KM173-39-62S
4. DD	0.9	.035	KM173-49-35	KM173-49-35S
	1.1	.045	KM173-49-45	KM173-49-45S
	1.4	.054	KM173-49-54	KM173-49-54S
	1.6	.062	KM173-49-62	KM173-49-62S

Sterilization Parameters

Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.

BUSA® Steinmann Pins - Smooth/Plain

Fabricated from precision ground 316L implant grade stainless steel (ASTM F 138).



1. Single Trocar
3-Sided Point; Point on One End



2. Double Trocar
3-Sided Point; Point on Both Ends



3. Single Diamond
2-Sided Point; Point on One End



4. Double Diamond
2-Sided Point; Point on Both Ends

Items not shown at actual size.

BUSA®				
Point Style	Diameter		Order Number	
			9" Length	
	mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
1. ST	2.0	5/64"	KM168-19-56	KM168-19-56S
	2.4	3/32"	KM168-19-33	KM168-19-33S
	2.8	7/64"	KM168-19-76	KM168-19-76S
	3.2	1/8"	KM168-19-18	KM168-19-18S
	3.6	9/64"	KM168-19-96	KM168-19-96S
	4.0	5/32"	KM168-19-53	KM168-19-53S
	4.8	3/16"	KM168-19-31	KM168-19-31S
2. DT	2.0	5/64"	KM168-29-56	KM168-29-56S
	2.4	3/32"	KM168-29-33	KM168-29-33S
	2.8	7/64"	KM168-29-76	KM168-29-76S
	3.2	1/8"	KM168-29-18	KM168-29-18S
	3.6	9/64"	KM168-29-96	KM168-29-96S
	4.0	5/32"	KM168-29-53	KM168-29-53S
	4.8	3/16"	KM168-29-31	KM168-29-31S

BUSA®				
Point Style	Diameter		Order Number	
			9" Length	
	mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
3. SD	2.0	5/64"	KM168-39-56	KM168-39-56S
	2.4	3/32"	KM168-39-33	KM168-39-33S
	2.8	7/64"	KM168-39-76	KM168-39-76S
	3.2	1/8"	KM168-39-18	KM168-39-18S
	3.6	9/64"	KM168-39-96	KM168-39-96S
	4.0	5/32"	KM168-39-53	KM168-39-53S
	4.8	3/16"	KM168-39-31	KM168-39-31S
4. DD	2.0	5/64"	KM168-49-56	KM168-49-56S
	2.4	3/32"	KM168-49-33	KM168-49-33S
	2.8	7/64"	KM168-49-76	KM168-49-76S
	3.2	1/8"	KM168-49-18	KM168-49-18S
	3.6	9/64"	KM168-49-96	KM168-49-96S
	4.0	5/32"	KM168-49-53	KM168-49-53S
	4.8	3/16"	KM168-49-31	KM168-49-31S

BUSA® 12 Inch Steinmann Pin

Fabricated from precision ground 316L Implant grade stainless steel (ASTM F 138).

BUSA®				
Point Style	Diameter		Order Number	OEM Catalog Number
	mm	in	Sterile (1/pk)	Zimmer Biomet®
4. DD	5/32	12	ZB27-362093	7-362093



12 INCH
4. Double Diamond
2-Sided Point; Point on Both Ends

BUSA® Steinmann Pins - Fully Threaded

Fabricated from precision ground 316L implant grade stainless steel (ASTM F 138).



1. Single Trocar
3-Sided Point; Point on One End



2. Double Trocar
3-Sided Point; Point on Both Ends



3. Single Diamond
2-Sided Point; Point on One End



4. Double Diamond
2-Sided Point; Point on Both Ends

Items not shown at actual size.

BUSA®				
Point Style	Diameter		Order Number	
			9" Length	
	mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
1. ST	2.0	5/64	KM169-19-56	KM169-19-56S
	2.4	3/32	KM169-19-33	KM169-19-33S
	2.8	7/64	KM169-19-76	KM169-19-76S
	3.2	1/8	KM169-19-18	KM169-19-18S
	3.6	9/64	KM169-19-96	KM169-19-96S
	4.0	5/32	KM169-19-53	KM169-19-53S
	4.8	3/16	KM169-19-31	KM169-19-31S
2. DT	2.0	5/64	KM169-29-56	KM169-29-56S
	2.4	3/32	KM169-29-33	KM169-29-33S
	2.8	7/64	KM169-29-76	KM169-29-76S
	3.2	1/8	KM169-29-18	KM169-29-18S
	3.6	9/64	KM169-29-96	KM169-29-96S
	4.0	5/32	KM169-29-53	KM169-29-53S
	4.8	3/16	KM169-29-31	KM169-29-31S

BUSA®				
Point Style	Diameter		Order Number	
			9" Length	
	mm	in	Non-Sterile (6/pk)	Sterile (1/pk)
3. SD	2.0	5/64	KM169-39-56	KM169-39-56S
	2.4	3/32	KM169-39-33	KM169-39-33S
	2.8	7/64	KM169-39-76	KM169-39-76S
	3.2	1/8	KM169-39-18	KM169-39-18S
	3.6	9/64	KM169-39-96	KM169-39-96S
	4.0	5/32	KM169-39-53	KM169-39-53S
	4.8	3/16	KM169-39-31	KM169-39-31S
4. DD	2.0	5/64	KM169-49-56	KM169-49-56S
	2.4	3/32	KM169-49-33	KM169-49-33S
	2.8	7/64	KM169-49-76	KM169-49-76S
	3.2	1/8	KM169-49-18	KM169-49-18S
	3.6	9/64	KM169-49-96	KM169-49-96S
	4.0	5/32	KM169-49-53	KM169-49-53S
	4.8	3/16	KM169-49-31	KM169-49-31S

Sterilization Parameters

Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes

- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.

BUSA® Orthopaedic Pin/Drill Packs Replacements

Compatible with Stryker® Knee Instrumentation Systems

Fabricated from 17-4 stainless steel.

Stryker®						
Order Number	Product Description	Diameter		Length		Packaging
		in	mm	in	mm	
BMA-S90008S	Headless Pin (Square End & Trocar Point)	1/8	3.1	3.5	88.9	4 Pins Per Pack - Sterile
Replaces 7650-2038, 7650-1038						
						
BMA-S90010S	Headless Pin, Fluted (Square End & Trocar Point)	1/8	3.1	3.2	88.9	4 Pins Per Pack - Sterile
Replaces 7650-2038A, 7650-1038A						
						
BMA-S90012S	Headless Pin (Square End & Trocar Point)	1/8	3.1	3.0	76.2	4 Pins Per Pack - Sterile
Replaces 6541-4-003, 6541-004-003						
						

Compatible with Stryker® Ortholock™ Navigation System

Fabricated from 17-4 stainless steel.

Stryker®						
Order Number	Product Description	Diameter		Length		Packaging
		in	mm	in	mm	
ST6007-003-100	Navigation Pin (Square End & Trocar Point)	1/8	3	3.9	100	2 Pins Per Pack - Sterile
Replaces 6007-003-100						
						
ST6007-003-150	Navigation Pin (Square End & Trocar Point)	1/8	3	5.9	150	2 Pins Per Pack - Sterile
Replaces 6007-003-150						
						

BUSA® Orthopaedic Pin/Drill Packs Replacements

Compatible with DePuy Synthes® Knee Instrumentation Systems

Fabricated from 17-4 stainless steel.

DePuy Synthes®						
Order Number	Product Description	Diameter		Length		Packaging
		in	mm	in	mm	
BMA-S90006S	Fixation Pin (Round End & Trocar Point)	1/8	3.1	3.0	76.2	4 Pins Per Pack - Sterile
Replaces 2490-95-000, 2490-95						
						
BMA-S90007S	Fixation Pin (Triangular End & Trocar Point) & Twist Drills (Round End)	1/8	3.1	3.0	76.2	4 Pins & 2 Drills Per Pack - Sterile
Replaces 86-4192, 86-4192-000						
						

Compatible with Zimmer Biomet® Knee Instrumentation Systems

Fabricated from 17-4 stainless steel.

Zimmer Biomet®						
Order Number	Product Description	Diameter		Length		Packaging
		in	mm	in	mm	
ZB00-5901-020-00	Headless Drill Pin (Hex End & Trocar Point)	1/8	3.2	2.95	75	4 Pins Per Pack - Sterile
Replaces 00-5901-020-00						
						

Sterilization Parameters

Steam Sterilization Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Pre-Vacuum (Wrapped)	270°F (132°C)	4 minutes	30 minutes
Gravity (Wrapped)	270°F (132°C)	15 minutes	30 minutes


- All cycles require a minimum of 30 minutes dry time.
- Observe AAMI and hospital guidelines for steam sterilization and the sterilizer manufacturer's written instructions for cycle parameters and load configurations.


BUSA® Guide Wire

Replacement - DePuy Synthes®

Compatible with DePuy Synthes® 3.5mm & 4.0mm Cannulated Screws





Fabricated from 17-4 stainless steel.

DePuy Synthes®						
Order Number	Product Description	Diameter		Length		Packaging
		in	mm	in	mm	
BR900-722	Partially Threaded Guide Wire with Diamond Point	0.049	1.25	5.9	150.0	1 Pin Per Pack - Non-Sterile
Replaces 900.722						
						



BUSA® Large Bone Saw Blades
Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1290-35	mm	90	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
Replaces Stryker® 6113-089-090, 6113-89-090, 2108-152-006					
					
BR1-1290-39	mm	90	12.5	0.99	0.99
	in	3.54	0.49	0.039	0.039
Replaces Stryker® 6113-097-090					
					
BR1-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
Replaces Stryker® 6113-119-090, 2108-352-000, 2108-352					
					
BR1-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
Replaces Stryker® 6113-127-090, 2108-356-000, 2108-356					
					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD




BUSA® Large Bone Saw Blades


Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)



STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1290-54	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
Replaces Stryker® 6113-137-090					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-1360-50	mm	60	13	1.27	1.27
	in	2.36	0.51	0.050	0.050
Replaces Stryker® 6113-127-060					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-1370-50F	mm	70	13	1.27	1.27
	in	2.76	0.51	0.050	0.050
Replaces Stryker® 6213-127-070					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-1375-50	mm	75	13	1.27	1.27
	in	2.95	0.51	0.050	0.050
Replaces Stryker® 6113-127-075					
 Aggressive Cut Diamond Cut No Tooth Offset					



BUSA® Large Bone Saw Blades
Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1390-42	mm	90	13	1.07	1.07
	in	3.54	0.51	0.042	0.042
Replaces Stryker® 6113-107-090					
					Aggressive Cut Diamond Cut No Tooth Offset
BR1-1875-35	mm	75	18	0.89	0.89
	in	2.95	0.71	0.035	0.035
Replaces Stryker® 6118-089-075					
					Aggressive Cut Diamond Cut No Tooth Offset
BR1-1875-50	mm	75	18	1.27	1.27
	in	2.95	0.71	0.050	0.050
Replaces Stryker® 6118-127-075					
					Aggressive Cut Diamond Cut No Tooth Offset
BR1-1990-35	mm	90	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
Replaces Stryker® 6118-089-090, 2108-109-001					
					Aggressive Cut Diamond Cut No Tooth Offset

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades

Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1290-54	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
Replaces Stryker® 6113-137-090					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-1990-39	mm	90	19	0.99	0.99
	in	3.54	0.75	0.039	0.039
Replaces Stryker® 6118-097-090					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-1990-47	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
Replaces Stryker® 6118-119-090					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-19100-47	mm	100	19	1.19	1.19
	in	3.94	0.75	0.047	0.047
Replaces Stryker® 6118-119-100					
 Aggressive Cut Diamond Cut No Tooth Offset					





STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD



BUSA® Large Bone Saw Blades

Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
Replaces Stryker® 6118-127-090, 2108-382-000, 2108-382					
					Aggressive Cut Diamond Cut No Tooth Offset
BR1-19100-50	mm	100	19	1.27	1.27
	in	3.94	0.75	0.050	0.050
Replaces Stryker® 6118-127-100, 2108-385-000, 2108-385					
					Aggressive Cut Diamond Cut No Tooth Offset
BR1-19110-50	mm	110	19	1.27	1.27
	in	4.33	0.75	0.050	0.050
Replaces Stryker® 6118-127-110					
					Aggressive Cut Diamond Cut No Tooth Offset
BR1-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
Replaces Stryker® 6118-137-090, 2108-351-02					
					Aggressive Cut Diamond Cut No Tooth Offset

BUSA® Large Bone Saw Blades

Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-19100-54	mm	100	19	1.37	1.37
	in	3.94	0.75	0.054	0.054
Replaces Stryker® 6118-137-100					
					
BR1-2190-39F	mm	90	21	0.99	0.99
	in	3.54	0.83	0.039	0.039
Replaces Stryker® 6221-100-090					
					
BR1-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
Replaces Stryker® 6221-119-090, 4221-120-085					
					
BR1-2190-50F	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
Replaces Stryker® 6221-127-090, 2108-382-001, 2108-382-1, 2108-382-003, 21080382-3					
					

STRYKER®





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD



BUSA® Large Bone Saw Blades

Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-2190-54F	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
Replaces Stryker® 6221-137-090, 2108-382-002, 2108-382-2					
				Aggressive Cut Diamond Cut No Tooth Offset Flared Blade	
BR1-2575-35	mm	75	25	0.89	0.89
	in	2.95	0.98	0.035	0.035
Replaces Stryker® 6125-089-075					
				Aggressive Cut Diamond Cut No Tooth Offset	
BR1-2590-35	mm	90	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
Replaces Stryker® 6125-089-090, 2108-302-000, 2108-302, 2108-102-003, 2108-102-008					
				Aggressive Cut Diamond Cut No Tooth Offset	
BR1-2590-39	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
Replaces Stryker® 6125-097-090, 2108-393-006, 2108-393-6					
				Aggressive Cut Diamond Cut No Tooth Offset	





STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

STRYKER®


BUSA® Large Bone Saw Blades

Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-2590-39F	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
Replaces Stryker® 2108-393-005, 2108-393-5					
				Aggressive Cut Diamond Cut No Tooth Offset Flared Blade	
BR1-2590-47	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
Replaces Stryker® 6125-119-090, 2108-351-000, 2108-351					
				Aggressive Cut Diamond Cut No Tooth Offset	
BR1-2590-47F	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
Replaces Stryker® 2108-351-002, 2108-351-2					
				Aggressive Cut Diamond Cut No Tooth Offset Flared Blade	
BR1-2590-50	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Replaces Stryker® 6125-127-090, 2108-376-002, 2108-386-2					
				Aggressive Cut Diamond Cut No Tooth Offset	

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD



BUSA® Large Bone Saw Blades
Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-25100-50	mm	100	25	1.27	1.27
	in	3.94	0.98	0.050	0.050
Replaces Stryker® 6125-127-100 & 2108-328-000, 2108-328					
					
BR1-2590-50F	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Replaces Stryker® 2108-383-000, 2108-383					
					
BR1-2590-54	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Replaces Stryker® 6125-137-090					
					
BR1-2590-54F	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Replaces Stryker® N/A					
					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD




BUSA® Large Bone Saw Blades


Sagittal/Oscillating BR2000™ Proprietary (BR1 Series - Stryker®)



STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-2590-58	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
Replaces Stryker® 6125-147-090, 2108-176-002, 2108-176-003					
 Aggressive Cut Diamond Cut No Tooth Offset					
BR1-2590-58F	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
Replaces Stryker® N/A					
 Aggressive Cut Diamond Cut No Tooth Offset Flared Blade					
BR1-3165-39F	mm	65	31	0.99	0.99
	in	2.56	1.22	0.039	0.039
Replaces Stryker® N/A					
 Aggressive Cut Diamond Cut No Tooth Offset Flared Blade					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4111-119-090	mm	90	11	0.99	1.19
	in	3.54	0.43	0.039	0.047
Replaces Stryker® 4111-119-090					
					
BR4111-127-090	mm	90	11	1.14	1.27
	in	3.54	0.43	0.045	0.050
Replaces Stryker® 4111-127-090					
					
BR4111-135-090	mm	90	11	1.19	1.35
	in	3.54	0.43	0.047	0.053
Replaces Stryker® 4111-135-090					
					
BR4111-147-075	mm	75	11	1.37	1.47
	in	2.95	0.43	0.054	0.058
Replaces Stryker® 4111-147-075					
					

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

STRYKER®





BUSA® Large Bone Saw Blades Sagittal Stryker® Replacements




STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4118-089-075	mm	75	18	0.79	0.89
	in	2.95	0.71	0.031	0.035
Replaces Stryker® 4118-089-075					
					
Aggressive Cut Tooth Offset					
BR4118-089-090	mm	90	18	0.79	0.89
	in	3.54	0.71	0.031	0.035
Replaces Stryker® 4118-089-090					
					
Aggressive Cut Tooth Offset					
BR4118-119-090	mm	90	18	0.99	1.19
	in	3.54	0.71	0.039	0.047
Replaces Stryker® 4118-119-090					
					
Aggressive Cut Tooth Offset					
BR4118-127-075	mm	75	18	1.13	1.27
	in	2.95	0.71	0.044	0.050
Replaces Stryker® 4118-127-075					
					
Aggressive Cut Tooth Offset					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4118-127-090	mm	90	18	1.14	1.27
	in	3.54	0.71	0.045	0.050
Replaces Stryker® 4118-127-090					
					
BR4118-127-100	mm	100	18	1.14	1.27
	in	3.94	0.71	0.045	0.050
Replaces Stryker® 4118-127-100					
					
BR4118-135-090	mm	90	18	1.19	1.35
	in	3.54	0.71	0.047	0.053
Replaces Stryker® 4118-135-090					
					
BR4118-135-105	mm	105	18	1.19	1.35
	in	4.13	0.71	0.047	0.053
Replaces Stryker® 4118-135-105					
					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades





Sagittal Stryker® Replacements




STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD




Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4118-137-090	mm	90	18	1.27	1.37
	in	3.54	0.71	0.050	0.054
Replaces Stryker® 4118-137-090					
					
BR4125-089-075	mm	75	25	0.79	0.89
	in	2.95	0.98	0.031	0.035
Replaces Stryker® 4125-089-075					
					
BR4125-089-090	mm	90	25	0.80	0.89
	in	3.54	0.98	0.031	0.035
Replaces Stryker® 4125-089-090					
					
BR4125-097-090	mm	90	25	0.88	0.97
	in	3.54	0.98	0.035	0.038
Replaces Stryker® 4125-097-090					
					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4125-107-090	mm	90	25	0.97	1.07
	in	3.54	0.98	0.038	0.042
Replaces Stryker® 4125-107-090					
					Aggressive Cut Tooth Offset
BR4125-119-090	mm	90	25	0.99	1.19
	in	3.54	0.98	0.039	0.047
Replaces Stryker® 4125-119-090					
					Aggressive Cut Tooth Offset
BR4125-127-090	mm	90	25	1.14	1.27
	in	3.54	0.98	0.045	0.050
Replaces Stryker® 4125-127-090					
					Aggressive Cut Tooth Offset
BR4125-127-100	mm	100	25	1.14	1.27
	in	3.94	0.98	0.045	0.050
Replaces Stryker® 4125-127-100					
					Aggressive Cut Tooth Offset

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD






STRYKER®


Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades

Sagittal Stryker® Replacements



Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4125-135-090	mm	90	25	1.19	1.35
	in	3.54	0.98	0.047	0.053
Replaces Stryker® 4125-135-090					
 Diagram of the Stryker BR4125-135-090 blade. The blade is 90 mm long and 25 mm wide. It features a 'FULL INSERT' label on the handle and a '10 mm' label on the blade body. The blade has a series of teeth along the right edge, labeled 'Aggressive Cut' and 'Tooth Offset'.					
BR4125-137-090	mm	90	25	1.27	1.37
	in	3.54	0.98	0.050	0.054
Replaces Stryker® 4125-137-090					
 Diagram of the Stryker BR4125-137-090 blade. The blade is 90 mm long and 25 mm wide. It features a 'FULL INSERT' label on the handle and a '10 mm' label on the blade body. The blade has a series of teeth along the right edge, labeled 'Aggressive Cut' and 'Tooth Offset'.					
BR4125-147-090	mm	90	25	1.37	1.47
	in	3.54	0.98	0.054	0.058
Replaces Stryker® 4125-147-090					
 Diagram of the Stryker BR4125-147-090 blade. The blade is 90 mm long and 25 mm wide. It features a 'FULL INSERT' label on the handle and a '10 mm' label on the blade body. The blade has a series of teeth along the right edge, labeled 'Aggressive Cut' and 'Tooth Offset'.					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR4135-089-064	mm	64	35	0.80	0.89
	in	2.52	1.38	0.031	0.035
Replaces Stryker® 4135-089-064					
 Diagram of the Stryker BR4135-089-064 blade. The blade is shown with a handle on the left and a cutting edge on the right. The handle has a 'FULL INSERT' label. The blade has a length scale from 0 to 60 mm. The cutting edge has a width of 35 mm. The blade is labeled 'Aggressive Cut' and 'Tooth Offset'.					
BR2108-107-4S1	mm	60	35	1.19	1.35
	in	2.36	1.38	0.047	0.053
Replaces Stryker® 2108-107-4S1 & 2108-107-4SI					
 Diagram of the Stryker BR2108-107-4S1 blade. The blade is shown with a handle on the left and a cutting edge on the right. The handle has a 'FULL INSERT' label. The blade has a length scale from 0 to 60 mm. The cutting edge has a width of 35 mm. The blade is labeled 'Aggressive Cut', 'Tooth Offset', and 'Specialty'.					





BUSA® Large Bone Saw Blades Sagittal Stryker® Replacements




STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD




Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-100	mm	70	25	0.64	0.99
	in	2.76	0.98	0.025	0.039
Replaces Stryker® 2108-100-000, 2108-100					
					
Tooth Offset					
BR2108-102	mm	85	25	0.89	0.94
	in	3.35	0.98	0.035	0.037
Replaces Stryker® 2108-102-000, 2108-102					
					
Tooth Offset					
BR2108-103	mm	90	12.5	0.89	0.99
	in	3.54	0.49	0.035	0.039
Replaces Stryker® 2108-103-000, 2108-103					
					
Tooth Offset					
BR2108-105	mm	70	25	0.89	0.94
	in	2.76	0.98	0.035	0.037
Replaces Stryker® 2108-105-000, 2108-105					
					
Tooth Offset					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-107-1	mm	73	29	0.89	0.97
	in	2.87	1.14	0.035	0.038
Replaces Stryker® 2108-107-001					
					
Tooth Offset					
BR2108-107-4	mm	64	36	1.19	1.55
	in	2.52	1.42	0.047	0.061
Replaces Stryker® 2108-107-004					
					
Tooth Offset					
BR2108-107-4S1	mm	60	35	1.19	1.35
	in	2.36	1.38	0.047	0.053
Replaces Stryker® 2108-107-4S1, 2108-107-4SI					
					
Aggressive Cut Tooth Offset					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD






STRYKER®


Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades

Sagittal Stryker® Replacements





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-107-9	mm	60	34.5	0.79	0.79
	in	2.36	1.36	0.031	0.031
Replaces Stryker® 2108-107-009, 2108-107-9					
					
BR2108-109	mm	80	21	0.89	0.94
	in	3.15	0.83	0.035	0.037
Replaces Stryker® 2108-109-000, 2108-109					
					
BR2108-109-000S9	mm	81	20	1.27	1.27
	in	3.19	0.79	0.05	0.05
Replaces Stryker® 2108-109-000S9					
					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-110	mm	70	25	1.27	1.47
	in	2.76	0.98	0.050	0.058
Replaces Stryker® 2108-110-000, 2108-110					
					
Tooth Offset					
BR2108-111	mm	75	25	1.27	1.47
	in	2.95	0.98	0.050	0.058
Replaces Stryker® 2108-111-000, 2108-111					
					
Tooth Offset					
BR2108-113	mm	30	16.5	1.27	1.65
	in	1.18	0.65	0.050	0.065
Replaces Stryker® 2108-113-000, 2108-113					
					
Tooth Offset					
BR2108-115	mm	70	18.5	0.79	0.91
	in	2.76	0.73	0.031	0.036
Replaces Stryker® 2108-115-000, 2108-115					
					
Tooth Offset					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





BUSA® Large Bone Saw Blades Sagittal Stryker® Replacements




STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD




Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-118	mm	75	25	1.27	1.27
	in	2.95	0.98	0.050	0.050
Replaces Stryker® 2108-118-000, 2108-118					
 No Tooth Offset Diamond Tooth Design					
BR2108-120	mm	60	18.5	0.64	0.97
	in	2.36	0.73	0.025	0.038
Replaces Stryker® 2108-120-000, 2108-120					
 Tooth Offset					
BR2108-120-6	mm	63	19	0.64	0.91
	in	2.48	0.73	0.025	0.036
Replaces Stryker® 2108-120-006					
 Tooth Offset					
BR2108-125	mm	45	42	0.64	1.02
	in	1.77	1.65	0.025	0.040
Replaces Stryker® 2108-125-000, 2108-125					
 Repeat Sternotomy Blade Tooth Offset					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-131	mm	47.5	18	0.64	1.04
	in	1.87	0.69	0.025	0.041
Replaces Stryker® 2108-131-000					
					
Tooth Offset					
BR2108-140	mm	55	29	0.64	1.02
	in	2.17	1.14	0.025	0.040
Replaces Stryker® 2108-140-000, 2108-140					
					
Tooth Offset					
BR2108-145	mm	20	9	0.64	0.79
	in	0.79	0.35	0.025	0.031
Replaces Stryker® 2108-145-000, 2108-145					
					
Tooth Offset					
BR2108-148	mm	30	9	0.64	0.79
	in	1.18	0.35	0.025	0.031
Replaces Stryker® 2108-148-000, 2108-148					
					
ACL Blade Tooth Offset					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades





Sagittal Stryker® Replacements




STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-150	mm	70	12.5	0.64	0.99
	in	2.76	0.49	0.025	0.039
Replaces Stryker® 2108-150-000, 2108-150					
					
Tooth Offset					
BR2108-151	mm	85	25	1.19	1.19
	in	3.35	0.98	0.047	0.047
Replaces Stryker® 2108-151-000, 2108-151					
					
No Tooth Offset Diamond Tooth Design					
BR2108-152	mm	80	12.5	1.19	1.19
	in	3.15	0.49	0.047	0.047
Replaces Stryker® 2108-152-000, 2108-152					
					
No Tooth Offset Diamond Tooth Design					
BR2108-155	mm	70	12.5	1.27	1.45
	in	2.76	0.49	0.050	0.057
Replaces Stryker® 2108-155-000, 2108-155					
					
Tooth Offset					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-156	mm	75	12.5	1.27	1.45
	in	2.95	0.49	0.050	0.057
Replaces Stryker® 2108-156-000, 2108-156					
					
Tooth Offset					
BR2108-158	mm	80	12.5	1.27	1.27
	in	3.15	0.49	0.050	0.050
Replaces Stryker® 2108-158-000, 2108-158					
					
No Tooth Offset Diamond Tooth Design					
BR2108-160	mm	70	25	0.89	1.27
	in	2.76	0.98	0.035	0.050
Replaces Stryker® 2108-160-000, 2108-160					
					
Tooth Offset					
BR2108-176	mm	75	25	1.47	1.47
	in	2.95	0.98	0.058	0.058
Replaces Stryker® 2108-176-000, 2108-176					
					
No Tooth Offset Diamond Tooth Design					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades





Sagittal Stryker® Replacements



STRYKER®





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-177	mm	85	25	1.32	1.32
	in	3.35	0.98	0.052	0.052
Replaces Stryker® 2108-177-000, 2108-177					
 No Tooth Offset Diamond Tooth Design					
BR2108-182	mm	80	18.5	1.27	1.27
	in	3.15	0.73	0.050	0.050
Replaces Stryker® 2108-182-000, 2108-182					
 No Tooth Offset Diamond Tooth Design					
BR2108-183	mm	80	29	1.27	1.27
	in	3.15	1.14	0.050	0.050
Replaces Stryker® 2108-183-000, 2108-183					
 No Tooth Offset Flared Blade Diamond Tooth Design					
BR2108-183-1	mm	84	29	0.97	0.97
	in	3.31	1.14	0.038	0.038
Replaces Stryker® 2108-183-001					
 No Tooth Offset Aggressive Cut					

BUSA® Large Bone Saw Blades Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-183-3	mm	84	29	1.32	1.32
	in	3.31	1.14	0.052	0.052
Replaces Stryker® 2108-183-003					
					
BR2108-185	mm	100	18.5	1.27	1.27
	in	3.94	0.73	0.050	0.050
Replaces Stryker® 2108-185-000, 2108-185					
					
BR2108-189	mm	80	25	1.27	1.27
	in	3.15	0.98	0.050	0.050
Replaces Stryker® 2108-189-000, 2108-189					
					
BR2108-193	mm	80	25	1.04	1.04
	in	3.15	0.98	0.041	0.041
Replaces Stryker® 2108-193-000, 2108-193					
					

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

STRYKER®







STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements



Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-195	mm	90	48.5	0.99	0.99
	in	3.54	1.91	0.039	0.039
Replaces Stryker® 2108-195-000, 2108-195					
					
BR2108-197	mm	60	48.5	0.64	0.99
	in	2.36	1.91	0.025	0.039
Replaces Stryker® 2108-197-000, 2108-197					
					
BR2108-218	mm	95	25	1.27	1.27
	in	3.74	0.98	0.050	0.050
Replaces Stryker® 2108-218-000, 2108-218					
					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-305	mm	73	25	0.89	0.89
	in	2.87	0.98	0.035	0.035
Replaces Stryker® 2108-305-000					
 Aggressive Cut No Tooth Offset					
BR2108-309	mm	84.5	21	0.89	0.89
	in	3.33	0.83	0.035	0.035
Replaces Stryker® 2108-309-000					
 Aggressive Cut No Tooth Offset					
BR2108-382-2	mm	89	19	1.37	1.37
	in	3.50	0.73	0.054	0.054
Replaces Stryker® 2108-382-002					
 Aggressive Cut No Tooth Offset					

STRYKER® | Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD






STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades Sagittal Stryker® Replacements





Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD
ConMed® Linvatec® Tooth Design

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR5072-101	mm	86	19.5	1.27	1.27
	in	3.39	0.77	0.050	0.050
Replaces ConMed® Linvatec® 5072-101, Stryker® 2108-109-000S15					
 Opposing Teeth No Tooth Offset Chromium Coating					
BR5072-181	mm	86	19.5	1.27	1.27
	in	3.39	0.77	0.050	0.050
Replaces ConMed® Linvatec® 5072-18, Stryker® 2108-109-000S11					
 Opposing Teeth No Tooth Offset					



BUSA® Large Bone Saw Blades
Sagittal Stryker® Replacements

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD
Zimmer Biomet® Tooth Design

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR11-2212	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Replaces Zimmer Biomet™ 11-2212, 25-0901-27S-G1, 25-0901-27T-G1, 25-0901-27Y-F1					
					
BR11-2214	mm	60	13	0.89	0.89
	in	2.36	0.51	0.035	0.035
Replaces Zimmer Biomet™ 1-2214, 13-0600-89S-F1, 13-0600-89T-F1, 13-0600-89Y-F1					
					
BR11-2220	mm	90	25	1.42	1.47
	in	3.54	0.98	0.056	0.058
Replaces Zimmer Biomet™ 11-2220, 25-0901-47S-B1, 25-0901-47T-B1, 25-0901-47Y-B1					
					
BR11-2222	mm	90	25	1.22	1.27
	in	3.54	0.98	0.048	0.050
Replaces Zimmer Biomet™ 11-2222, 25-0901-27S-B1, 25-0901-27T-B1, 25-0901-27Y-B1					
					

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

STRYKER®





STRYKER®

Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

BUSA® Large Bone Saw Blades

Sagittal Stryker® Replacements






Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD
Zimmer Biomet® Tooth Design

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR11-2325	mm	90	13	1.27	1.27
	in	3.54	0.51	0.050	0.050
Replaces Zimmer Biomet™ 11-2325, 13-0901-27S-F1, 13-0901-27T-F1, 13-0901-27Y-F1					
					
BR11-2608	mm	70	13	1.27	1.27
	in	2.76	0.51	0.050	0.050
Replaces Zimmer Biomet™ 11-2608, 13-0701-27S-L1, 13-0701-27T-L1, 13-0701-27Y-L1; Stryker® 6213-127-070					
					

BUSA® Large Bone Saw Blades

Sagittal Stryker® Replacements

Compatible with Stryker® Reciprocating Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-250R	mm	60	6.3	0.64	1.04
	in	2.36	0.25	0.025	0.041
Replaces Stryker® 277-096-250, 277-96-250					
					Tooth Offset Short Single Sided
KM-275R	mm	70	12.5	0.64	0.94
	in	2.76	0.49	0.025	0.037
Replaces Stryker® 277-096-275, 277-96-275					
					Tooth Offset Double Sided
KM-276R	mm	71.4	12.5	0.79	0.99
	in	2.81	0.49	0.031	0.039
Replaces Stryker® 277-096-276, 277-96-276					
					Tooth Offset Double Sided
KM-277R	mm	70.8	12.5	0.63	1.12
	in	2.79	0.49	0.025	0.044
Replaces Stryker® 277-096-277, 277-96-277					
					Tooth Offset Double Sided
KM-281R	mm	70.0	12.7	1.27	1.40
	in	2.76	0.50	0.050	0.055
Replaces Stryker® 277-096-281					
					Tooth Offset Double Sided

BUSA® Large Bone Saw Blades




Sagittal Stryker® Replacements



STRYKER®

Compatible with Stryker® Reciprocating Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower





Compatible with Stryker® Reciprocating Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-300R	mm	77.5	6.3	0.64	1.04
	in	3.05	0.25	0.025	0.041
Replaces Stryker® 277-096-300, 277-96-300					
					Long Tooth Offset Single Sided
KM-325R	mm	77.5	11.2	0.76	1.17
	in	3.05	0.44	0.030	0.046
Replaces Stryker® 277-096-325, 277-96-325					
					Tooth Offset Single Sided
KM-326R	mm	77.5	11.0	1.24	1.60
	in	3.05	0.43	0.049	0.063
Replaces Stryker® 277-096-326					
					Tooth Offset Single Sided
KM-2279-73	mm	76	12.7	0.90	0.90
	in	2.99	0.50	0.035	0.035
Replaces Stryker® "Replaces N/A - Similar to ConMed® Linvatec® KM-479"					
					No Tooth Offset Double Sided Chromium Coating, Ribbed



BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ Proprietary - ConMed® Linvatec® K7 Series

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-1290-35	mm	90	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
Replaces ConMed® Linvatec® 5071-340, T130-089-90, TN130-089-90, V130-089-90					
					Aggressive Cut No Tooth Offset
K7-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
Replaces ConMed® Linvatec® 5071-343, T130-119-90, TN130-119-90, V130-119-90					
					Aggressive Cut No Tooth Offset
K7-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
Replaces ConMed® Linvatec® 5071-344, T130-127-90, TN130-127-90, V130-127-90					
					Aggressive Cut No Tooth Offset
K7-1990-35	mm	90	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
Replaces ConMed® Linvatec® 5071-300, T190-089-90, TN190-089-90, V190-089-90					
					Aggressive Cut No Tooth Offset

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

CONMED® LINVATEC®





BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ Proprietary - ConMed® Linvatec® K7 Series



CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus





Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-1990-47	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
Replaces ConMed® Linvatec® 5071-303, T190-119-90, TN190-119-90, V190-119-90					
					
Aggressive Cut No Tooth Offset					
K7-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
Replaces ConMed® Linvatec® 5071-304, T190-127-90, TN190-127-90, V190-127-90					
					
Aggressive Cut No Tooth Offset					
K7-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
Replaces ConMed® Linvatec® 5071-305, T190-137-90, TN190-137-90, V190-137-90					
					
Aggressive Cut No Tooth Offset					
K7-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
Replaces ConMed® Linvatec® 5071-313, T193-119-90, TN193-119-90, V193-119-90					
					
Flared Blade Aggressive Cut No Tooth Offset					



BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ Proprietary - ConMed® Linvatec® K7 Series

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-2590-35	mm	90	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
Replaces ConMed® Linvatec® 5071-320, T250-089-90, TN250-089-90, V250-089-90					
					Aggressive Cut No Tooth Offset
K7-2590-39	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
Replaces ConMed® Linvatec® 5071-321, T250-100-90, TN250-100-90, V250-100-90					
					Aggressive Cut No Tooth Offset
K7-2590-47	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
Replaces ConMed® Linvatec® 5071-323, T250-119-90, TN250-119-90, V250-119-90					
					Aggressive Cut No Tooth Offset
K7-2590-50	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Replaces ConMed® Linvatec® 5071-324, T250-127-90, TN250-127-90, V250-127-90					
					Aggressive Cut No Tooth Offset

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

CONMED® LINVATEC®





CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ Proprietary - ConMed® Linvatec® K7 Series

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus





ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K7-2590-54	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Replaces ConMed® Linvatec® 5071-325, T250-137-90, TN250-137-90, V250-137-90					
					
K7-2590-58	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
Replaces ConMed® Linvatec® 5071-326, T250-147-90, TN250-147-90, V250-147-90					
					



BUSA® Large Bone Saw Blades - Oscillating
ConMed® Linatec® Replacements

Compatible with ConMed® Linatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-201	mm	88	19.5	1.27	1.27
	in	3.46	0.77	0.050	0.050
Replaces ConMed® Linatec® 5071-201					
 Opposing Teeth No Tooth Offset Chromium Coating					
KM71-233	mm	25	9.5	0.40	0.60
	in	0.98	0.37	0.016	0.024
Replaces ConMed® Linatec® 5071-133, 5071-233					
 Opposing Teeth Tooth Offset ACL Blade					
KM71-523	mm	97	19.5	1.27	1.27
	in	3.79	0.77	0.050	0.050
Replaces ConMed® Linatec® 5071-123, 5071-423, 5071-523					
 Opposing Teeth No Tooth Offset					
KM71-524	mm	97	25.4	1.27	1.27
	in	3.79	1.00	0.050	0.050
Replaces ConMed® Linatec® 5071-124, 5071-424, 5071-524					
 Opposing Teeth No Tooth Offset					

Compatible with ConMed® Linatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

CONMED® LINATEC®





BUSA® Large Bone Saw Blades - Oscillating ConMed® Linvatec® Replacements



CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus




Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-545	mm	70	13	1.20	1.50
	in	2.76	0.51	0.047	0.059
Replaces ConMed® Linvatec® 5071-145, 5071-445, 5071-545					
 No Opposing Teeth Tooth Offset					
KM71-551	mm	65	32.5	0.61	0.90
	in	2.56	1.28	0.024	0.035
Replaces ConMed® Linvatec® 5071-151, 5071-451, 5071-551					
 Opposing Teeth Tooth Offset Flared Blade					
KM71-552	mm	65	13.5	0.61	0.90
	in	2.56	0.53	0.024	0.035
Replaces ConMed® Linvatec® 5071-152, 5071-452, 5071-552					
 No Opposing Teeth Tooth Offset Flared Blade					
KM71-553	mm	65	19.5	0.61	0.90
	in	2.56	0.77	0.024	0.035
Replaces ConMed® Linvatec® 5071-153, 5071-453, 5071-553					
 Opposing Teeth Tooth Offset					



BUSA® Large Bone Saw Blades - Oscillating ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-556	mm	65	33	0.80	0.90
	in	2.56	1.30	0.031	0.035
Replaces ConMed® Linvatec® 5071-156, 5071-456, 5071-556					
 Opposing Teeth Tooth Offset					
KM71-561	mm	92	25.4	1.27	1.27
	in	3.62	1.00	0.050	0.050
Replaces ConMed® Linvatec® 5071-161, 5071-461, 5071-561					
 Opposing Teeth No Tooth Offset					
KM71-563	mm	92	25.4	1.37	1.37
	in	3.62	1.00	0.054	0.054
Replaces ConMed® Linvatec® 5071-163, 5071-463, 5071-563					
 Opposing Teeth No Tooth Offset					

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

CONMED® LINVATEC®

BUSA® Large Bone Saw Blades - Oscillating ConMed® Linvatec® Replacements



CONMED® LINVATEC®


Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-567	mm	92	19.5	1.20	1.20
	in	3.62	0.77	0.047	0.047


Replaces ConMed® Linvatec® 5071-167, 5071-467, 5071-567



Chromium Coating
Opposing Teeth
No Tooth Offset

KM71-572	mm	70	19.5	0.80	0.90
	in	2.76	0.77	0.031	0.035


Replaces ConMed® Linvatec® 5071-172, 5071-472, 5071-572



Opposing Teeth
Tooth Offset

KM71-580	mm	85	10	1.27	1.27
	in	3.35	0.39	0.050	0.050


Replaces ConMed® Linvatec® 5071-180, 5071-480, 5071-580




No Opposing Teeth
No Tooth Offset
Diamond Tooth Design

KM71-581	mm	88	19.5	1.27	1.27
	in	3.46	0.77	0.050	0.050

Replaces ConMed® Linvatec® 5071-181, 5071-481, 5071-581





Opposing Teeth
No Tooth Offset



BUSA® Large Bone Saw Blades - Oscillating ConMed® Linatec® Replacements

Compatible with ConMed® Linatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

ConMed® Linatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM71-585	mm	82	19.5	0.90	0.90
	in	3.23	0.77	0.035	0.035
Replaces ConMed® Linatec® 5071-185, 5071-485, 5071-585					
 Opposing Teeth No Tooth Offset					
KM71-586	mm	65	32	0.80	0.90
	in	2.56	1.26	0.031	0.035
Replaces ConMed® Linatec® 5071-186, 5071-486, 5071-586					
 Opposing Teeth Tooth Offset					

Compatible with ConMed® Linatec® Oscillating Saw Systems:
HALL 50™, HALL® Titan®, MPower®, PowerPro® & VersiPower® Plus

CONMED® LINATEC®

BUSA® Large Bone Saw Blades - Oscillating ConMed® Linvatec® Replacements



CONMED® LINVATEC®



Compatible with ConMed® Linvatec® Oscillating Saw Systems:
HALL® Titan®, VersiPower®, Series 4, Series 3 & Orthair

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-251	mm	60	32.5	0.61	0.90
	in	2.36	1.28	0.024	0.035
Replaces ConMed® Linvatec® 5052-51, 5052-051, 5052-251, 5052-451, 5052-551, Z251					
					Flared Blade Opposing Teeth Tooth Offset
KM-256	mm	60	33.4	0.80	0.90
	in	2.36	1.31	0.031	0.035
Replaces ConMed® Linvatec® “5052-56, 5052056, 5052-256, 5052-456, 5052-556, Z256					
					Opposing Teeth Tooth Offset
KM-451	mm	60	32.5	0.61	0.90
	in	2.36	1.28	0.024	0.035
Replaces ConMed® Linvatec® 5052-51, 5052-051, 5052-251, 5052-451, 5052-551					
					Flared Blade Tooth Offset
KM-453	mm	60	19.5	0.61	0.92
	in	2.36	0.77	0.024	0.036
Replaces ConMed® Linvatec® 5052-53, 5052-053, 5052-153, 5052-253, 5052-453, 5052-553					
					Tooth Offset Opposing Teeth

BUSA® Large Bone Saw Blades - Reciprocating Replacements - ConMed® Linvatec®





Compatible with ConMed® Linvatec® Reciprocating Saw Systems:

HALL 50™, HALL® Titan®, MPower®, PowerPro®, VersiPower® Plus, VersiPower®, Series 4, Series 3 &
Orthair; MicroAire® & Smith® & Nephew® Large Bone Reciprocating Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-458	mm	89	12.5	0.61	1.00
	in	3.50	0.49	0.024	0.039
Replaces ConMed® Linvatec® 5052-58, 5052-058, 5052-158, 5052-258, Z258					
					Tooth Offset
KM-459	mm	54.1	9.4	0.61	1.00
	in	2.13	0.37	0.024	0.039
Replaces ConMed® Linvatec® 5052-59, 5052-059, 5052-159, 5052-259, Z259					
					Tooth Offset
KM-461	mm	89	12.5	1.50	1.60
	in	3.50	0.49	0.059	0.063
Replaces ConMed® Linvatec® 5052-60, 5052-061, 5052-161, 5052-261, Z261					
					Tooth Offset
KM-473	mm	73	12.5	0.80	1.00
	in	2.87	0.49	0.031	0.039
Replaces ConMed® Linvatec® 5052-73, 5052-073, 5052-173, 5052-273, Z273					
					Tooth Offset Double Sided
KM-479	mm	77.5	12.7	0.90	0.90
	in	3.05	0.50	0.035	0.035
Replaces ConMed® Linvatec® 5052-79, 5052-079, 5052-179, 5052-279					
					Chromium Coating No Tooth Offset Double Sided


BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000 - K5 Series - DePuy Synthes® & Zimmer Biomet®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DePuy Synthes®					
Zimmer Biomet®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1290-35	mm	90	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
Replaces DePuy Synthes® 05.002.201S, 05.003.100S					
Replaces Zimmer Biomet® 13-0900-89A-F1					
 Aggressive Cut No Tooth Offset					
K5-1290-39	mm	90	12.5	0.99	1.00
	in	3.54	0.49	0.039	0.039
Replaces DePuy Synthes® 05.002.405S, 05.003.102S					
Replaces Zimmer Biomet® 13-0901-00A-F1					
 Aggressive Cut No Tooth Offset					
K5-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
Replaces DePuy Synthes® 05.002.202S, 05.003.104S					
Replaces Zimmer Biomet® 13-0901-19A-F1					
 Aggressive Cut No Tooth Offset					
K5-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
Replaces DePuy Synthes® 05.002.203S, 05.003.106S					
Replaces Zimmer Biomet® 13-0901-27A-F1					
 Aggressive Cut No Tooth Offset					

DEPUY SYNTHES®
ZIMMER BIOMET®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal



BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000 - K5 Series - DePuy Synthes® & Zimmer Biomet®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:

Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw

Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal


DePuy Synthes®

Zimmer Biomet®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1290-54	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054

Replaces DePuy Synthes® 05.002.204S, 05.003.107S

Replaces Zimmer Biomet® 13-0901-37A-F1




Aggressive Cut
No Tooth Offset

K5-1990-35

mm	90	19	0.89	0.89
in	3.54	0.75	0.035	0.035

Replaces DePuy Synthes® 05.002.001S, 05.003.110S

Replaces Zimmer Biomet® 19-0900-89A-F1




Aggressive Cut
No Tooth Offset

K5-1990-39

mm	90	19	0.99	0.99
in	3.54	0.75	0.039	0.039

Replaces DePuy Synthes® 05.002.007S, 05.003.112S

Replaces Zimmer Biomet® 19-0901-00A-F1




Aggressive Cut
No Tooth Offset

K5-1990-47

mm	90	19	1.19	1.19
in	3.54	0.75	0.047	0.047

Replaces DePuy Synthes® 05.002.002S, 05.003.114S

Replaces Zimmer Biomet® 19-0901-19A-G1



Aggressive Cut
No Tooth Offset





Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DEPUY SYNTHES®
ZIMMER BIOMET®

BUSA® Large Bone Saw Blades - Sagittal/Oscillating


BR2000 - K5 Series - DePuy Synthes® & Zimmer Biomet®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DePuy Synthes®					
Zimmer Biomet®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
Replaces DePuy Synthes® 05.002.003S, 05.003.116S					
Replaces Zimmer Biomet® 19-0901-27A-F1, 19-0901-27A-G1					
					Aggressive Cut No Tooth Offset
K5-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
Replaces DePuy Synthes® 05.002.004S, 05.003.117S					
Replaces Zimmer Biomet® 19-0901-37A-G1					
					Aggressive Cut No Tooth Offset
K5-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
Replaces DePuy Synthes® 05.002.302S, 05.003.124S					
Replaces Zimmer Biomet® 21-0901-19A-V1					
					Flared Blade Aggressive Cut No Tooth Offset
K5-2190-50F	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
Replaces DePuy Synthes® 05.002.303S, 05.003.126S					
Replaces Zimmer Biomet® 21-0901-27A-V1					
					Flared Blade Aggressive Cut No Tooth Offset

DEPUY SYNTHES®
ZIMMER BIOMET®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
 Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
 Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal



Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw

Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DEPUY SYNTHES®
ZIMMER BIOMET®

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000 - K5 Series - DePuy Synthes® & Zimmer Biomet®




Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:

Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw

Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DePuy Synthes®

Zimmer Biomet®




Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-2190-54F	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
Replaces DePuy Synthes® 05.002.304S, 05.003.127S					
Replaces Zimmer Biomet® 21-0901-37A-V1					
 Flared Blade Aggressive Cut No Tooth Offset					
K5-2590-35	mm	90	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
Replaces DePuy Synthes® 05.002.101S, 05.003.130S					
Replaces Zimmer Biomet® 25-0900-89A-F1					
 Aggressive Cut No Tooth Offset					
K5-2590-39	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
Replaces DePuy Synthes® 05.002.106S, 05.003.132S					
Replaces Zimmer Biomet® 25-0901-00A-F1					
 Aggressive Cut No Tooth Offset					

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000 - K5 Series - DePuy Synthes® & Zimmer Biomet®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:


Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw

Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DePuy Synthes®					
Zimmer Biomet®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-2590-47	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
Replaces DePuy Synthes® 05.002.102S, 05.003.134S					
Replaces Zimmer Biomet® 25-0901-19A-B1, 25-0901-19A-F1, 25-0901-19A-G1					
 Aggressive Cut No Tooth Offset					
K5-2590-50	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.047
Replaces DePuy Synthes® 05.002.102S, 05.003.134S					
Replaces Zimmer Biomet® 25-0901-19A-B1, 25-0901-19A-F1, 25-0901-19A-G1					
 Aggressive Cut No Tooth Offset					
K5-2590-58	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
Replaces DePuy Synthes® 05.002.105S, 05.003.139S					
Replaces Zimmer Biomet® 25-0901-47A-B1					
 Aggressive Cut No Tooth Offset					

DEPUY SYNTHES®
ZIMMER BIOMET®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal



Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DEPUY SYNTHES®

BUSA® Large Bone Saw Blades - Oscillating DePuy Synthes® Replacements

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:


Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw

Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DePuy Synthes®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM519-10	mm	50.0	27	0.40	0.66
	in	1.97	1.06	0.016	0.026


Replaces DePuy Synthes® 519.10, 519.100S



Tooth Offset

KM519-15	mm	50.0	14	0.40	0.66
	in	1.97	0.55	0.016	0.026


Replaces DePuy Synthes® 519.15, 519.150S



Tooth Offset

KM519-17	mm	70	27	0.64	0.81
	in	2.76	1.06	0.025	0.032

Replaces DePuy Synthes® 519.17, 519.170S





Tooth Offset

BUSA® Large Bone Saw Blades - Oscillating

DePuy Synthes® Replacements

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

DePuy Synthes®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM519-21	mm	70	27	0.99	1.19
	in	2.76	1.06	0.039	0.047
Replaces DePuy Synthes® 519.21, 519.210S					
 Diagram of the KM519-21 saw blade. It is a long, narrow blade with a curved handle on the left. The blade has a series of teeth along the right edge. A scale in millimeters is shown along the top edge, ranging from 0 to 70. The text 'mm' is written vertically next to the scale. The text 'Tooth Offset' is written at the bottom right of the diagram.					


KM-1201-84	mm	86	19.5	1.27	1.27
	in	3.39	0.77	0.050	0.050
Replaces DePuy Synthes® N/A (Similar to ConMed® Linvatec® 5071-201)					
 Diagram of the KM-1201-84 saw blade. It is a long, narrow blade with a curved handle on the left. The blade has a series of teeth along the right edge. The text 'Special', 'No Tooth Offset', 'Opposing Teeth', and 'Chromium Coating' are written at the bottom right of the diagram.					

DEPUY SYNTHES®

Compatible with DePuy Synthes® Large Bone Oscillating Saw Systems:
 Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System Sagittal Saw
 Compatible with Zimmer Biomet® Large Bone Oscillating Saw System: Universal

BUSA® Large Bone Saw Blades - Reciprocating DePuy Synthes® Replacements

Compatible with DePuy Synthes® Large Bone Reciprocating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System (TRS) Sagittal Saw

DePuy Synthes®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR511-905S	mm	80	10	1.19	1.05
	in	3.15	0.39	0.047	0.041
Replaces DePuy Synthes® 511.905S					
					Tooth Offset Single-Sided
BR511-907S	mm	55	10	1.19	1.05
	in	2.17	0.39	0.047	0.041
Replaces DePuy Synthes® 511.907S					
					Tooth Offset Single-Sided
BR511-909S	mm	55	10	1.19	0.85
	in	2.17	0.39	0.047	0.033
Replaces DePuy Synthes® 511.909S					
					Tooth Offset Single-Sided
BR511-912S	mm	68	12	1.19	1.09
	in	2.68	0.47	0.047	0.043
Replaces DePuy Synthes® 511.912S					
					Tooth Offset Single-Sided
BR511-915S	mm	40	10	1.19	1.10
	in	1.57	0.39	0.047	0.043
Replaces DePuy Synthes® 511-915S					
					Tooth Offset Single-Sided

Compatible with DePuy Synthes® Large Bone Reciprocating Saw Systems:
Air Power Line, Battery Power Line, Compact Air Drive & Trauma Recon System (TRS) Sagittal Saw


DEPUY SYNTHES®



BUSA® Large Bone Saw Blade - Reciprocating
Zimmer Biomet® Replacements

Compatible with Zimmer Biomet® Reciprocating Saw System: Universal

Zimmer Biomet®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
ZB12-0761-00U-R1	mm	70	12.5	0.61	0.90
	in	2.76	0.49	0.024	0.035
Replaces Zimmer Biomet® 12-0761-00U-R1					
					Tooth Offset Double-Sided Sterile







ZB12-0901-19U-R1	89.4	12.2	0.90	1.14	1.27
	3.52	0.48	0.035	0.045	0.050
Replaces Zimmer Biomet® 12-0901-19U-R1					
					Special Tooth Offset Single-Sided Sterile



BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-51	mm	18.5	6.5	0.38	0.43
	in	0.73	0.26	0.015	0.017
Replaces Stryker® 2296-3-51, 2296-003-051					
					Thin
KM3-101	mm	10	9	0.38	0.43
	in	0.39	0.35	0.015	0.017
Replaces Stryker® 2296-3-101, 2296-003-101					
					Thin
KM3-102	mm	18.5	16.5	0.38	0.43
	in	0.73	0.65	0.015	0.017
Replaces Stryker® 2296-3-102, 2296-003-102					
					Thin
KM3-103	mm	15	7	0.38	0.43
	in	0.59	0.28	0.015	0.017
Replaces Stryker® 2296-3-103, 2296-003-103					
					Thin
KM3-104	mm	34.5	16.5	0.38	0.43
	in	1.36	0.65	0.015	0.017
Replaces Stryker® 2296-3-104, 2296-003-104					
					Thin
KM3-105	mm	18.5	9	0.38	0.43
	in	0.73	0.35	0.015	0.017
Replaces Stryker® 2296-3-105, 2296-003-105					
					Thin

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

STRYKER®

BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements








STRYKER®

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2






Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-106	mm	34.5	13	0.38	0.43
	in	1.36	0.51	0.015	0.017
Replaces Stryker® 2296-3-106, 2296-003-106					
					Thin
KM3-108	mm	39	13	0.38	0.43
	in	1.54	0.51	0.015	0.017
Replaces Stryker® 2296-3-108, 2296-003-108					
					Thin
KM3-111	mm	25	9	0.38	0.43
	in	0.98	0.35	0.015	0.017
Replaces Stryker® 2296-3-111, 2296-003-111					
					Thin
KM3-114	mm	18.5	7	0.38	0.43
	in	0.73	0.28	0.015	0.017
Replaces Stryker® 2296-3-114, 2296-003-114					
					Thin
KM3-115	mm	29.5	7	0.38	0.43
	in	1.16	0.28	0.015	0.017
Replaces Stryker® 2296-3-115, 2296-003-115					
					Thin
KM3-125	mm	31	9	0.38	0.43
	in	1.22	0.35	0.015	0.017
Replaces Stryker® 2296-3-125, 2296-003-125					
					Thin



BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-206	mm	34.5	12	0.38	0.43
	in	1.36	0.47	0.015	0.017
Replaces Stryker® 2296-3-206, 2296-003-206					
					Thin
KM3-212	mm	18	5.5	0.38	0.43
	in	0.71	0.22	0.015	0.017
Replaces Stryker® 2296-3-212, 2296-003-212					
					Thin
KM3-225	mm	31	9	0.38	0.43
	in	1.22	0.35	0.015	0.017
Replaces Stryker® 2296-3-225, 2296-003-225					
					Thin
KM3-254	mm	41	19.5	0.38	0.68
	in	1.61	0.77	0.015	0.027
Replaces Stryker® 2296-3-254, 2296-003-254					
					Thin
KM3-410	mm	11.5	5.5	0.38	0.43
	in	0.45	0.22	0.015	0.017
Replaces Stryker® 2296-3-410, 2296-003-410					
					Thin

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

STRYKER®

BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements









STRYKER®

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2







Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-411	mm	11.5	9	0.38	0.43
	in	0.45	0.35	0.015	0.017
Replaces Stryker® 2296-3-411, 2296-003-411					
					Thin
KM3-412	mm	18	5.5	0.38	0.43
	in	0.71	0.22	0.015	0.017
Replaces Stryker® 2296-3-412, 2296-003-412					
					Thin
KM3-413	mm	18	2.5	0.38	0.43
	in	0.71	0.10	0.015	0.017
Replaces Stryker® 2296-3-413, 2296-003-413					
					Thin
KM3-414	mm	25	5.5	0.38	0.43
	in	0.98	0.22	0.015	0.017
Replaces Stryker® 2296-3-414, 2296-003-414					
					Thin
KM3-511	mm	25	9	0.38	0.64
	in	0.98	0.35	0.015	0.025
Replaces Stryker® 2296-3-511, 2296-003-511					
					Thin Aggressive Cut
KM3-525	mm	31	9	0.38	0.64
	in	1.22	0.35	0.015	0.025
Replaces Stryker® 2296-3-525, 2296-003-525					
					Thin Aggressive Cut



BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-101	mm	10	9	0.51	0.56
	in	0.39	0.35	0.020	0.022
Replaces Stryker® 2296-33-101, 2296-033-101					
					Thick
KM33-102	mm	18.5	16.5	0.51	0.56
	in	0.73	0.65	0.020	0.022
Replaces Stryker® 2296-33-102, 2296-033-102					
					Thick
KM33-103	mm	15	7	0.51	0.56
	in	0.59	0.28	0.020	0.022
Replaces Stryker® 2296-33-103, 2296-033-103					
					Thick
KM33-104	mm	34.5	16.5	0.51	0.56
	in	1.36	0.65	0.020	0.022
Replaces Stryker® 2296-33-104, 2296-033-104					
					Thick
KM33-105	mm	18.5	9	0.51	0.56
	in	0.73	0.35	0.020	0.022
Replaces Stryker® 2296-33-105, 2296-033-105					
					Thick
KM33-106	mm	34.5	13	0.51	0.56
	in	1.36	0.51	0.020	0.022
Replaces Stryker® 2296-33-106, 2296-033-106					
					Thick

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

STRYKER®

BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements








STRYKER®

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2







Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-108	mm	39	13	0.51	0.56
	in	1.54	0.51	0.020	0.022
Replaces Stryker® 2296-33-108, 2296-033-108					
					Thick
KM33-111	mm	25	9	0.51	0.56
	in	0.98	0.35	0.020	0.022
Replaces Stryker® 2296-33-111, 2296-033-111					
					Thick
KM33-114	mm	18.5	7	0.51	0.56
	in	0.73	0.28	0.020	0.022
Replaces Stryker® 2296-33-114, 2296-033-114					
					Thick
KM33-115	mm	29.5	7	0.51	0.56
	in	1.16	0.28	0.020	0.022
Replaces Stryker® 2296-33-115, 2296-033-115					
					Thick
KM33-125	mm	31	9	0.51	0.56
	in	1.22	0.35	0.020	0.022
Replaces Stryker® 2296-33-125, 2296-033-125					
					Thick
KM33-206	mm	34.5	12	0.51	0.56
	in	1.36	0.47	0.020	0.022
Replaces Stryker® 2296-33-206, 2296-033-206					
					Thick



BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:

SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-212	mm	18	5.5	0.51	0.56
	in	0.71	0.22	0.020	0.022
Replaces Stryker® 2296-33-212, 2296-033-212					
					Thick
KM33-225	mm	31	9	0.51	0.56
	in	1.22	0.35	0.020	0.022
Replaces Stryker® 2296-33-225, 2296-033-225					
					Thick
KM33-410	mm	11.5	5.5	0.51	0.56
	in	0.45	0.22	0.020	0.022
Replaces Stryker® 2296-33-410, 2296-033-410					
					Thick
KM33-411	mm	11.5	9	0.51	0.56
	in	0.45	0.35	0.020	0.022
Replaces Stryker® 2296-33-411, 2296-033-411					
					Thick
KM33-412	mm	18	5.5	0.51	0.56
	in	0.71	0.22	0.020	0.022
Replaces Stryker® 2296-33-412, 2296-033-412					
					Thick
KM33-413	mm	18	2.5	0.51	0.56
	in	0.71	0.10	0.020	0.022
Replaces Stryker® 2296-33-413, 2296-033-413					
					Thick

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

STRYKER®





STRYKER®

Compatible with Stryker® Oscillating & Sagittal Saw Systems:
SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:




SABO, Cordless Driver (Sagittal Saw Attachment), RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-414	mm	25	5.5	0.51	0.56
	in	0.98	0.22	0.020	0.022
Replaces Stryker® 2296-33-414, 2296-033-414					
					Thick
KM33-521	mm	25	9	0.51	0.74
	in	0.98	0.35	0.020	0.029
Replaces Stryker® 2296-33-521, 2296-033-521					
					Thick Aggressive Cut

BUSA® Small Bone Reciprocating Saw Blades

Stryker® Replacements





Compatible with Stryker® Reciprocating Saw Systems:
RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM5100-137-123	mm	27.0	6.4	0.38	0.51
	in	1.06	0.25	0.015	0.020
Replaces Stryker® 5100-137-123					
					Thin
KM5100-137-133	mm	21.8	6.4	0.38	0.51
	in	0.86	0.25	0.015	0.020
Replaces Stryker® 5100-137-133					
					Thin
KM5100-137-233	mm	27.4	6	0.38	0.51
	in	1.08	0.24	0.015	0.020
Replaces Stryker® 5100-137-233					
					Thin

BUSA® Small Bone Reciprocating Rasps


Stryker® Replacements

Compatible with Stryker® Reciprocating Saw Systems:
RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®			
Order Number	Unit	Length	Max. Cut Width
KM5100-37-113	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 5100-37-113, 5100-037-113			
			Small Tear Cross Cut
KM5100-37-114	mm	13	6.5
	in	0.51	0.26
Replaces Stryker® 5100-37-114, 5100-037-114			
			Large Tear Cross Cut
KM5100-37-115	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 5100-37-115, 5100-037-115			
			Small Tear
KM5100-37-116	mm	13	6.5
	in	0.51	0.26
Replaces Stryker® 5100-37-116, 5100-037-116			
			Large Tear

STRYKER®

Compatible with Stryker® Reciprocating Saw Systems:
RemB®, CORE™, TPS™, F1, TPX® & Command® 2



BUSA® Small Bone Reciprocating Saw Blades
Stryker® Replacements

Compatible with Stryker® Reciprocating Saw System: Micro




Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-062	mm	27.2	6	0.64	1.00
	in	1.07	0.24	0.025	0.039
Replaces Stryker® 1675-133					
					Thick

BUSA® Small Bone Reciprocating Rasps

Stryker® Replacements



Compatible with Stryker® Reciprocating Saw System: Micro

Stryker®			
Order Number	Unit	Length	Max. Cut Width
KM-800	mm	13	6.5
	in	0.51	0.26
Replaces Stryker® 1675-116			
			Large Tear
KM-801	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 1675-115			
			Small Tear
KM-802	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 1675-113			
			Small Tear Cross-Cut





STRYKER® | Compatible with Stryker® Reciprocating Saw System: Micro



BUSA® Small Bone Sagittal Saw Blades

ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Sagittal Saw Systems:
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
 & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3100C	mm	41.5	14	0.40	0.61
	in	1.63	0.55	0.016	0.024
Replaces ConMed® Linvatec® 5023-33, 5023-033, 5023-133, 5023-233, 5053-33, 5053-033, 5053-133, 5053-233					
					Coarse Straight Edge
KM-3101	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces ConMed® Linvatec® 5023-38, 5023-038, 5023-138, 5023-238, 5053-38, 5053-038, 5053-138, 5053-238, Z338					
					Fine Straight Edge
KM-3101C	mm	25.7	9.5	0.38	0.61
	in	1.01	0.37	0.015	0.024
Replaces ConMed® Linvatec® 5023-43, 5023-043, 5023-143, 5023-243, 5053-43, 5053-043, 5053-143, 5053-243					
					Coarse Straight Edge
KM-3102	mm	25.7	4.5	0.40	0.60
	in	1.01	0.18	0.016	0.024
Replaces ConMed® Linvatec® 5023-39, 5023-039, 5023-139, 5023-239, 5053-39, 5053-039, 5053-139, 5053-239, Z339					
					Fine Straight Edge
KM-3104	mm	16.1	9.5	0.40	0.63
	in	0.63	0.37	0.016	0.025
Replaces ConMed® Linvatec® 5044-145, Z332					
					Fine Straight Edge

Compatible with ConMed® Linvatec® Sagittal Saw Systems:
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
 & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

CONMED® LINVATEC®

BUSA® Small Bone Sagittal Saw Blades

ConMed® Linvatec® Replacements



CONMED® LINVATEC®





Compatible with ConMed® Linvatec® Sagittal Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

Compatible with ConMed® Linvatec® Sagittal Saw Systems:

MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®

& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3105	mm	16.1	4.5	0.40	0.63
	in	0.63	0.18	0.016	0.025
Replaces ConMed® Linvatec® 5023-044, 5023-144, 5023-244, 5053-044, 5053-144, 5053-244					
					Fine Straight Edge
KM-3106	mm	25.7	5.5	0.40	0.61
	in	1.01	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-32, 5023-032, 5023-132, 5023-232, 5053-32, 5053-032, 5053-132, 5053-232					
					Coarse Straight Edge
KM-3107	mm	34.3	6	0.40	0.60
	in	1.35	0.24	0.016	0.024
Replaces ConMed® Linvatec® 5023-41, 5023-041, 5023-141, 5023-241, 5053-41, 5053-041, 5053-141, 5053-241, Z341					
					Fine Straight Edge
KM-3183	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces ConMed® Linvatec® 5023-183					
					Fine Straight Edge
KM-3231	mm	41.5	19.5	0.40	0.61
	in	1.63	0.77	0.016	0.024
Replaces ConMed® Linvatec® 5023-31, 5023-031, 5023-131, 5023-231, 5053-31, 5053-031, 5053-131, 5053-231, Z331					
					Coarse Straight Edge





Compatible with ConMed® Linvatec® Sagittal Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

CONMED® LINVATEC®

BUSA® Small Bone Sagittal Saw Blades

ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Sagittal Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3237	mm	25.7	14	0.40	0.61
	in	1.01	0.55	0.016	0.024
Replaces ConMed® Linvatec® 5023-37, 5023-037, 5023-137, 5023-237, 5053-37, 5053-037, 5053-137, 5053-237					
					Coarse Straight Edge
KM-3240	mm	41.5	9.5	0.40	0.60
	in	1.63	0.37	0.016	0.024
Replaces ConMed® Linvatec® 5023-40, 5023-040, 5023-140, 5023-240, 5053-40, 5053-040, 5053-140, 5053-240, Z340					
					Fine Tapered Straight Edge




BUSA® Small Bone Oscillating Saw Blades

ConMed® Linvatec® Replacements



CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems




ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3000	mm	18.5	5.5	0.40	0.60
	in	0.73	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-62, 5023-062, 5023-162, 5023-262, 5053-62, 5053-062, 5053-162, 5053-262, Z362					
					Fine Straight Edge
KM-3001	mm	28.5	5.5	0.40	0.60
	in	1.12	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-64, 5023-064, 5023-164, 5023-264, 5053-64, 5053-064, 5053-164, 5053-264, Z364					
					Fine Straight Edge
KM-3002	mm	34.7	10	0.40	0.60
	in	1.37	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-69, 5023-069, 5023-169, 5023-269, 5053-73, 5053-073, 5053-173, 5053-273, Z373					
					Fine Straight Edge
KM-3003	mm	22	5.5	0.40	0.60
	in	0.87	0.22	0.016	0.024
Replaces ConMed® Linvatec® N/A					
					Fine Straight Edge
KM-3004	mm	27.9	9.5	0.40	0.60
	in	1.10	0.37	0.016	0.024
Replaces ConMed® Linvatec® Z365					
					Fine Straight Edge
KM-3005	mm	18.5	10	0.40	0.60
	in	0.73	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-63, 5023-063, 5023-163, 5023-263, 5053-63, 5053-063, 5053-163, 5053-263					
					Fine Straight Edge



BUSA® Small Bone Oscillating Saw Blades

ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3006	mm	12	5.5	0.40	0.60
	in	0.47	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-60, 5023-060, 5023-160, 5023-260, 5053-60, 5053-060, 5053-160, 5053-260, Z360					
					Fine Straight Edge
KM-3008	mm	28.5	18.5	0.40	0.60
	in	1.12	0.73	0.016	0.024
Replaces ConMed® Linvatec® 5023-66, 5023-066, 5023-166, 5023-266, 5053-66, 5053-066, 5053-166, 5053-266, Z366					
					Fine Radial Edge
KM-3261	mm	12	10	0.40	0.60
	in	0.47	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-61, 5023-061, 5023-161, 5023-261, 5053-61, 5053-061, 5053-161, 5053-261, Z361					
					Fine Straight Edge
KM-3265	mm	28.5	10	0.40	0.60
	in	1.12	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-65, 5023-065, 5023-165, 5023-265, 5053-65, 5053-065, 5053-165, 5053-265					
					Fine Radial Edge

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

BUSA® Small Bone Oscillating Saw Blades

ConMed® Linvatec® Replacements




CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Oscillating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

Compatible with ConMed® Linvatec® Oscillating Saw Systems:

MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3170	mm	12	9	0.40	0.60
	in	0.47	0.36	0.016	0.024
Replaces ConMed® Linvatec® 5023-70, 5023-170, 5023-270, 5053-70, 5053-170, 5053-270, Z370					
					Fine Crescentic
KM-3171	mm	25.5	13.8	0.40	0.60
	in	1.00	0.55	0.016	0.024
Replaces ConMed® Linvatec® 5023-71, 5023-171, 5023-271, 5053-71, 5053-171, 5053-271, Z371					
					Fine Crescentic
KM-3172	mm	32	18	0.40	0.60
	in	1.26	0.71	0.016	0.024
Replaces ConMed® Linvatec® 5023-72, 5023-172, 5023-272, 5053-72, 5053-172, 5053-272, Z372					
					Fine Crescentic
KM-3176	mm	32	13.8	0.40	0.60
	in	1.26	0.55	0.016	0.024
Replaces ConMed® Linvatec® 5023-76, 5023-176, 5023-276, 5053-76, 5053-176, 5053-276					
					Fine Crescentic

BUSA® Small Bone Reciprocating Saw Blades ConMed® Linvatec® Replacements





Compatible with ConMed® Linvatec® Reciprocating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-060	mm	25.5	7.6	0.40	0.64
	in	1.00	0.30	0.016	0.025
Replaces ConMed® Linvatec® 5053-20, 5053-020, 5053-120, 5053-220, Z320					
					Flat Shank
KM-070	mm	38	7.6	0.40	0.64
	in	1.50	0.30	0.016	0.025
Replaces ConMed® Linvatec® 5053-21, 5053-021, 5053-121, 5053-221, Z321					
					Flat Shank
KM-54	mm	25.0	7.4	0.40	0.60
	in	0.984	0.29	0.016	0.024
Replaces ConMed® Linvatec® 5053-54, 5053-054, 5053-154, 5053-254, Z354					
					Round Shank

BUSA® Small Bone Reciprocating Rasps


ConMed® Linatec® Replacements

Compatible with ConMed® Linatec® Reciprocating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems

ConMed® Linatec®			
Order Number	Unit	Length	Max. Cut Width
KM-80	mm	11	5.5
	in	0.43	0.22
Replaces ConMed® Linatec® 5053-080, 5053-80			
			
Small Tear Round Shank			
KM-81	mm	13	6.9
	in	0.51	0.27
Replaces ConMed® Linatec® 5053-081, 5053-81			
			
Large Tear Round Shank			
KM-82	mm	11	5.5
	in	0.43	0.22
Replaces ConMed® Linatec® 5053-082, 5053-82			
			
Small Tear Cross Cut Round Shank			
KM-83	mm	13	6.9
	in	0.51	0.27
Replaces ConMed® Linatec® 5053-083, 5053-83			
			
Large Tear Cross Cut Round Shank			





CONMED® LINATEC®

Compatible with ConMed® Linatec® Reciprocating Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems



BUSA® Small Bone Sagittal Saw Blades
ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Mini-Driver™
K100 & K200 Sagittal Saw Systems - K120 Attachment

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-21K	mm	15.1	6.4	0.71	0.89
	in	0.59	0.25	0.028	0.035
Replaces ConMed® Linvatec® K121					
 Fine					
KM-22K	mm	31	6.4	0.71	0.89
	in	1.22	0.25	0.028	0.035
Replaces ConMed® Linvatec® K122					
 Fine					
KM-23K	mm	15.1	10.4	0.71	0.89
	in	0.59	0.41	0.028	0.035
Replaces ConMed® Linvatec® K123					
 Fine					
KM-24K	mm	31	10.4	0.71	0.89
	in	1.22	0.41	0.028	0.035
Replaces ConMed® Linvatec® K124					
 Fine					

Compatible with ConMed® Linvatec® Mini-Driver™
K100 & K200 Sagittal Saw Systems - K120 Attachment

CONMED® LINVATEC®

BUSA® Small Bone Sagittal Saw Blades




ConMed® Linvatec® Replacements



CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Mini-Driver™
K100 & K200 Sagittal Saw Systems - K120 Attachment






**Compatible with ConMed® Linvatec® Mini-Driver™
K100 & K200 Sagittal Saw Systems - K120 Attachment**

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-31K	mm	15.1	6.4	0.38	0.56
	in	0.59	0.25	0.015	0.022
Replaces ConMed® Linvatec® K131					
					Fine
KM-33K	mm	15.1	10.4	0.38	0.56
	in	0.59	0.41	0.015	0.022
Replaces ConMed® Linvatec® K133					
					Fine
KM-34K	mm	31	10.4	0.38	0.56
	in	1.22	0.41	0.015	0.022
Replaces ConMed® Linvatec® K134					
					Fine



BUSA® Small Bone Sagittal Saw Blades
ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Mini-Driver™
K100 & K200 Sagittal Saw Systems - K220 Attachment

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-41K	mm	14.2	6.4	0.71	0.97
	in	0.56	0.25	0.028	0.038
Replaces ConMed® Linvatec® K141					
 Fine					
KM-43K	mm	14.2	10.2	0.71	0.97
	in	0.56	0.40	0.028	0.038
Replaces ConMed® Linvatec® K143					
 Fine					
KM-44K	mm	30	10.2	0.71	0.97
	in	1.18	0.40	0.028	0.038
Replaces ConMed® Linvatec® K144					
 Fine					
KM-48K	mm	27.4	9.5	0.38	0.65
	in	1.08	0.37	0.015	0.026
Replaces ConMed® Linvatec® K148					
 Fine					
KM-51K	mm	14.2	6.4	0.38	0.65
	in	0.56	0.25	0.015	0.026
Replaces ConMed® Linvatec® K151					
 Fine					

Compatible with ConMed® Linvatec® Mini-Driver™
K100 & K200 Sagittal Saw Systems - K220 Attachment

CONMED® LINVATEC®







BUSA® Sternum Saw Blades


Stryker® Replacements



STRYKER® | Compatible with Stryker® Sternum Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower



Compatible with Stryker® Sternum Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
ST6207-97-101	mm	30.5	6.5	0.64	1.02
	in	1.20	0.26	0.025	0.040
Replaces Stryker® 6207-97-101, 6207-097-101					
					1 Per Pack Sterile
ST6207-97-101N	mm	30.5	6.5	0.64	1.02
	in	1.20	0.26	0.025	0.040
Replaces Stryker® 6207-97-101, 6207-097-101					
					1 Per Pack Non-Sterile
ST296-97-102	mm	31.5	6.2	0.64	1.01
	in	1.24	0.24	0.025	0.04
Replaces Stryker® 296-91-102, 296-097-102					
					1 Per Pack Sterile
ST296-97-102N	mm	31.5	6.2	0.64	1.01
	in	1.24	0.24	0.025	0.04
Replaces Stryker® 296-91-102, 296-097-102					
					1 Per Pack Non-Sterile
KM-97	mm	32	6.4	0.79	1.07
	in	1.26	0.25	0.031	0.042
Replaces Stryker® 298-97-100, 298-097-100					
					1 Per Pack Sterile
KM-97N	mm	32	6.4	0.79	1.07
	in	1.26	0.25	0.031	0.042
Replaces Stryker® 2296-3-105, 2296-003-105					
					1 Per Pack Non-Sterile



BUSA® Sternum Saw Blades
Stryker® Replacements


Compatible with Stryker® Sternum Saw Systems: 9, 8, 7, 6, 5, 4, 2000, EHD & OrthoPower

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-98	mm	32	6.4	0.64	1.02
	in	1.26	0.25	0.025	0.040
Replaces Stryker® 298-97-101, 298-097-101					
					1 Per Pack Sterile
KM-98N	mm	32	6.4	0.64	1.02
	in	1.26	0.25	0.025	0.040
Replaces Stryker® 298-97-101, 298-097-101					
					1 Per Pack Non-Sterile



BUSA® Sternum Saw Blades
Stryker® Replacements



Compatible with Stryker® Sagittal Saw Systems: 9, 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-125	mm	45	42	0.64	1.02
	in	1.77	1.65	0.025	0.040
Replaces Stryker® 2108-125, 2108-125-000					
 1 Per Pack Sterile Tooth Offset Repeat Sternotomy Blade					

BUSA® Sternum Saw Blades

ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Sternum Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-31	mm	54	9.9	0.64	0.97
	in	2.13	0.39	0.025	0.038
Replaces ConMed® Linvatec® 5059-31, 5059-031, 5059-531, 5059-031-01					
 Sterile					
KM-31N	mm	54	9.9	0.64	0.97
	in	2.13	0.39	0.025	0.038
Replaces ConMed® Linvatec® 5059-31, 5059-031, 5059-531, 5059-031-01					
 Non-Sterile					
KM-32	mm	32.8	9.4	0.64	0.99
	in	1.29	0.37	0.025	0.039
Replaces ConMed® Linvatec® 5059-32, 5059-032, 5059-532, 5059-032-01					
 Sterile					
KM-32N	mm	32.8	9.4	0.64	0.99
	in	1.29	0.37	0.025	0.039
Replaces ConMed® Linvatec® 5059-32, 5059-032, 5059-532, 5059-032-01					
 Non-Sterile					





BUSA® Sternum Saw Blades

ConMed® Linvatec® Replacement & Proprietary

Compatible with ConMed® Linvatec® Oscillating Saw Systems: VersiPower®, Series 4 , Series 3 & Orthair;
ConMed® Linvatec® Micro 100™ Sagittal Saw System; ConMed® Linvatec® Repeat Sternotomy Saw System;
BUSA®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-279	mm	39.5	55	0.40	0.79
	in	1.56	2.17	0.016	0.031
Replaces ConMed® Linvatec® 5059-279, Stryker® 5301-10-279, DePuy® Synthes® 5608-13, MicroAire® Z0-079					
					
Sterile Repeat Sternotomy Blade Opposing Teeth					
KM-279PZ	mm	26	39	0.40	0.71
	in	1.02	1.54	0.016	0.028
Replaces ConMed® Linvatec® N/A					
					
Sterile Repeat Sternotomy Blade					

CONMED® LINVATEC®

Compatible with ConMed® Linvatec® Oscillating Saw Systems: VersiPower®, Series 4 , Series 3 & Orthair;
ConMed® Linvatec® Micro 100™ Sagittal Saw System; ConMed® Linvatec® Repeat Sternotomy Saw System;
BUSA®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems



BUSA® Sternum Saw Blades


Proprietary

Compatible with ConMed® Linvatec® Small Bone/Micro Sagittal Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3179SN	mm	42.5	54	0.40	0.69
	in	1.67	2.13	0.016	0.027


Replaces ConMed® Linvatec® N/A



KM-3231SN

mm	41.9	19	0.40	0.71
in	1.65	0.75	0.016	0.028


Replaces ConMed® Linvatec® N/A



KM-3231SS

mm	41.9	19	0.40	0.71
in	1.65	0.75	0.016	0.028

Replaces ConMed® Linvatec® N/A




Compatible with ConMed® Linvatec® Small Bone/Micro Sagittal Saw Systems:
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

CONMED® LINVATEC®

BUSA® Sternum Saw Blades

Codman® (Symmetry) Replacements

Compatible with Codman® Sternum Saw Systems

Codman®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-7010	mm	34	9.4	0.89	1.20
	in	1.34	0.37	0.035	0.047
Replaces Codman® 26-7010					
 1 Per Pack Sterile					

BUSA® Sternum Saw Blades

Replacements Terumo






Compatible with Sarns® Sternal Saw II

Codman®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-278	mm	32	6.4	0.79	1.02
	in	1.26	0.25	0.031	0.040
Replaces Terumo 5589					
 1 Per Pack Sterile					
KM-278N	mm	32	6.4	0.79	1.02
	in	1.26	0.25	0.031	0.040
Replaces Terumo 5756					
 1 Per Pack Non-Sterile					

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements

Compatible with Stryker® Signature Portfolio Elite Attachments




Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-010-120	Round Fluted Bur - Less Aggressive	mm	2.0	2.0
		in	0.08	0.08
Replaces Stryker® 5820-010-120				
				
BR5820-010-130	Round Fluted Bur - Less Aggressive	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5820-010-130				
				
BR5820-010-140	Round Fluted Bur - Less Aggressive	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5820-010-140				
				
BR5820-010-150	Round Fluted Bur - Less Aggressive	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5820-010-150				
				
BR5820-010-160	Round Fluted Bur - Less Aggressive	mm	6.0	6.0
		in	0.24	0.24
Replaces Stryker® 5820-010-160				
				
BR5820-010-170	Round Fluted Bur - Less Aggressive	mm	7.0	7.0
		in	0.28	0.28
Replaces Stryker® 5820-010-170				
				

STRYKER® | Compatible with Stryker® Signature Portfolio Elite Attachments

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements







Compatible with Stryker® Signature Portfolio Elite Attachments

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-010-230	Round Fluted Bur - Aggressive	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5820-010-230				
				
BR5820-010-240	Round Fluted Bur - Aggressive	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5820-010-240				
				
BR5820-010-250	Round Fluted Bur - Aggressive	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5820-010-250				
				
BR5820-107-020	Neuro Match Head Drill	mm	2.0	3.3
		in	0.08	0.13
Replaces Stryker® 5820-107-020				
				
BR5820-107-025	Neuro Match Head Drill	mm	2.5	3.7
		in	0.10	0.15
Replaces Stryker® 5820-107-025				
				
BR5820-107-030	Neuro Match Head Drill	mm	3.0	3.8
		in	0.12	0.15
Replaces Stryker® 5820-107-030				
				

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements

Compatible with Stryker® Signature Portfolio Elite Attachments







Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-107-040	Neuro Match Head Drill	mm	4.0	5.1
		in	0.16	0.20
Replaces Stryker® 5820-107-040				
				
BR5820-107-430	Neuro Match Head Drill - Less Aggressive	mm	3.0	3.8
		in	0.12	0.15
Replaces Stryker® 5820-107-430				
				
BR5820-107-430S1	Neuro Match Head Drill - Specialty	mm	3.0	3.8
		in	0.12	0.15
Replaces Stryker® 5820-107-430S1				
				
BR5820-107-530	Neuro Match Head Drill - Exact	mm	3.0	3.8
		in	0.12	0.15
Replaces Stryker® 5820-107-530				
				
BR5820-107-530C	Neuro Match Head Drill - Carbide	mm	3.0	3.8
		in	0.12	0.15
Replaces Stryker® 5820-107-530C				
				
BR5820-010-640	Cross Cut Round Bur	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5820-010-640				
				

STRYKER® | Compatible with Stryker® Signature Portfolio Elite Attachments

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements





Compatible with Stryker® Signature Portfolio Elite Attachments

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-010-650	Cross Cut Round Bur	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5820-010-650				
				
BR5820-010-660	Cross Cut Round Bur	mm	6.0	6.0
		in	0.24	0.24
Replaces Stryker® 5820-010-660				
				
BR5820-012-020	Round Diamond Bur	mm	2.0	2.0
		in	0.08	0.08
Replaces Stryker® 5820-012-020				
				
BR5820-012-030	Round Diamond Bur	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5820-012-030				
				
BR5820-012-040	Round Diamond Bur	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5820-012-040				
				
BR5820-012-050	Round Diamond Bur	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5820-012-050				
				

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements

Compatible with Stryker® Signature Portfolio Elite Attachments






Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-012-060	Round Diamond Bur	mm	6.0	6.0
		in	0.24	0.24
Replaces Stryker® 5820-012-060				
				
BR5820-013-020	Round Diamond Bur - Coarse	mm	2.0	2.0
		in	0.08	0.08
Replaces Stryker® 5820-013-020				
				
BR5820-013-030	Round Diamond Bur - Coarse	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5820-013-030				
				
BR5820-013-040	Round Diamond Bur - Coarse	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5820-013-040				
				
BR5820-013-050	Round Diamond Bur - Coarse	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5820-013-050				
				
BR5820-013-140	Round Diamond Bur - X-Coarse	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5820-013-140				
				

STRYKER® | Compatible with Stryker® Signature Portfolio Elite Attachments

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements





Compatible with Stryker® Signature Portfolio Elite Attachments

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-013-150	Round Diamond Bur - X-Coarse	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5820-013-150				
				
BR5820-012-310	Tapered Round Bur Diamond	mm	1.0	1.0
		in	0.04	0.04
Replaces Stryker® 5820-012-310				
				
BR5820-012-315	Tapered Round Bur Diamond	mm	1.5	1.5
		in	0.06	0.06
Replaces Stryker® 5820-012-315				
				
BR5820-012-320	Tapered Round Bur Diamond	mm	2.0	2
		in	0.08	0.08
Replaces Stryker® 5820-012-320				
				
BR5820-107-130	Neuro Diamond Drill	mm	3.0	3.8
		in	0.12	0.15
Replaces Stryker® 5820-107-130				
				
BR5820-015-040	Egg Bur	mm	4.0	9.9
		in	0.16	0.39
Replaces Stryker® 5820-015-040				
				

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements

Compatible with Stryker® Signature Portfolio Elite Attachments


Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-015-050	Egg Bur	mm	5.0	9.9
		in	0.20	0.39
Replaces Stryker® 5820-015-050				
				
BR5820-070-011	Straight Router	mm	1.1	6
		in	0.04	0.24
Replaces Stryker® 5820-070-011				
				
BR5820-080-021	Helocoidal Rasp	mm	2.1	19.6
		in	0.08	0.77
Replaces Stryker® 5820-080-021				
				
BR5820-103-010	Wire Pass Drill	mm	1.0	6.7
		in	0.04	0.26
Replaces Stryker® 5820-103-010				
				
BR5820-103-015	Wire Pass Drill	mm	1.5	6.7
		in	0.06	0.26
Replaces Stryker® 5820-103-015				
				
BR5820-103-020	Wire Pass Drill	mm	2.0	6.7
		in	0.08	0.26
Replaces Stryker® 5820-103-020				
				

STRYKER® | Compatible with Stryker® Signature Portfolio Elite Attachments

Brasseler USA® Telescoping High Speed Burs & Drills

Stryker® Replacements

Compatible with Stryker® Signature Portfolio Elite Attachments

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
BR5820-066-012	Cross Cut Fissure Bur	mm	1.2	4.4
		in	0.05	0.17
Replaces Stryker® 5820-066-012				
				
BR5820-066-016	Cross Cut Fissure Bur	mm	1.6	4.7
		in	0.06	0.19
Replaces Stryker® 5820-066-016				
				
BR5820-066-021	Cross Cut Fissure Bur	mm	2.1	2.1
		in	0.08	0.08
Replaces Stryker® 5820-066-021				
				
BR5820-066-216	Cross Cut Fissure Bur - Tungsten Carbide	mm	1.6	6.2
		in	0.06	0.24
Replaces Stryker® 5820-066-216				
				



BUSA® High Speed Burs

Stryker® 3.16mm TPS™ Shank Replacements



Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM-M8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Stryker® 5120-107-30, 5220-107-30, 5420-107-30, ConMed® Linvatec® 7021-237, Medtronic® Midas Rex® M-8				
				
KM-AM8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Stryker® 5130-107-30, 5230-107-30, 5430-107-30, ConMed® Linvatec® 7021-337, Medtronic® Midas Rex® AM-8				
				
KM-M51U	Metal Cutter (Diamond Wheel)	mm	25.4	0.8
		in	1.00	0.03
Replaces Stryker® 5120-130-250, ConMed® Linvatec® 7021-292, Medtronic® Midas Rex® M-51				
				
KM-M52U	Metal Cutter (Carbide)	mm	3.1	17.5
		in	0.12	0.69
Replaces Stryker® 5120-80-30, 5420-80-30, ConMed® Linvatec® 7021-294, Medtronic® Midas Rex® M-52				
				
KM-U295	Oval	mm	4.0	8.0
		in	0.16	0.31
Replaces Stryker® 5120-15-140, ConMed® Linvatec® 7021-295, Medtronic® Midas Rex® N/A				
				

STRYKER® | Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™


CONMED® LINVATEC® | Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

BUSA® High Speed Burs

Stryker® 3.16mm TPS™ Shank Replacements

Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM5110-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Stryker® 5110-15-40, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				
KM-M340DU	Round Diamond	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5120-12-40, ConMed® Linvatec® 7021-264, Medtronic® Midas Rex® M-340D				
				
KM-M330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5120-10-30, 5220-10-30, ConMed® Linvatec® 7021-226, Medtronic® Midas Rex® M-330				
				
KM5120-10-230	Round Fluted Extended Length	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5120-10-230, 5220-10-230, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				
KM-M340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5120-10-40, 5220-10-40, ConMed® Linvatec® 7021-227, Medtronic® Midas Rex® M-340				
				
KM5120-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5120-10-240, 5220-10-240, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				

BUSA® High Speed Burs



Stryker® 3.16mm TPS™ Shank Replacements



Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

STRYKER® | Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™
CONMED® LINVATEC® | Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II


Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM-M35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5120-10-50, 5220-10-50, 5420-10-50, ConMed® Linvatec® 7021-228, Medtronic®Midas Rex® M-35				
				
KM-M33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Stryker® 5120-10-60, 5220-10-60, ConMed® Linvatec® 7021-229, Medtronic® Midas Rex® M-33				
				
KM-AM330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5130-10-30, 5230-10-30, ConMed® Linvatec® 7021-326, Medtronic® Midas Rex® AM-330				
				
KM5130-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5130-10-240, 5230-10-240, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				
KM-AM35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5130-10-50, 5230-10-50, 5430-10-50, ConMed® Linvatec® 7021-328, Medtronic® Midas Rex® AM-35				
				
KM5130-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5130-10-250, 5230-10-250, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				

BUSA® High Speed Burs

Stryker® 3.16mm TPS™ Shank Replacements

Compatible with High Speed Drill Systems: Stryker® CORE™ & TPS™






Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM-AM32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Stryker® 5130-10-75, 5230-10-75, ConMed® Linvatec® 7021-331, Medtronic® Midas Rex® AM-32				
				

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2



Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0901	Oval - Carbide	mm	44.5	4.0	7.9	8
		in	1.75	0.16	0.31	
Replaces Stryker® 277-10-62, 277-010-062						
						
1 per pack, Sterile						
14-C0905	Oval - Carbide	mm	44.5	5.5	10.3	10
		in	1.75	0.22	0.41	
Replaces Stryker® N/A						
						
1 per pack, Sterile						
16-C0901	Oval - Carbide	mm	54.5	4.0	7.9	8
		in	2.15	0.16	0.31	
Replaces Stryker® N/A						
						
1 per pack, Sterile						
14-C0104	Round - Carbide	mm	44.5	1.8	1.8	6
		in	1.75	0.07	0.07	
Replaces Stryker® 277-10-418, 277-010-418						
						
1 per pack, Sterile						
14-C0105	Round - Carbide	mm	44.5	2.0	2.0	6
		in	1.75	0.08	0.08	
Replaces Stryker® 277-10-225, 277-010-225						
						
1 per pack, Sterile						
14-C0106	Round - Carbide	mm	44.5	2.3	2.3	8
		in	1.75	0.09	0.09	
Replaces Stryker® 277-10-423, 277-010-423						
						
8 Shannon 1 per pack, Sterile						

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements



Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0107	Round - Carbide	mm	44.5	3.0	3.0	6
		in	1.75	0.12	0.12	
Replaces Stryker® 277-10-230, 277-010-230						
						
1 per pack, Sterile						
14-C0108	Round - Carbide	mm	44.5	4.0	4.0	8
		in	1.75	0.16	0.16	
Replaces Stryker® 277-10-240, 277-010-240						
						
1 per pack, Sterile						
14-C0109	Round - Carbide	mm	44.5	5.0	5.0	8
		in	1.75	0.20	0.20	
Replaces Stryker® 277-10-250, 277-010-250						
						
1 per pack, Sterile						
14-C0110	Round - Carbide	mm	44.5	6.0	6.0	8
		in	1.75	0.24	0.24	
Replaces Stryker® N/A						
						
1 per pack, Sterile						
16-C0107	Round - Carbide	mm	54.5	3.0	3.0	6
		in	2.15	0.12	0.12	
Replaces Stryker® N/A						
						
1 per pack, Sterile						
16-C0108	Round - Carbide	mm	54.5	4.0	4.0	8
		in	2.15	0.16	0.16	
Replaces Stryker® N/A						
						
1 per pack, Sterile						

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0308	Straight Cross Cut - Carbide	mm	44.5	1.6	4.5	6
		in	1.75	0.06	0.17	
Replaces Stryker® 1607-2-103, 1607-002-103						
						560 1 per pack, Sterile
14-C0405	Straight Cross Cut - Carbide	mm	44.5	1.6	5.9	6
		in	1.75	0.06	0.23	
Replaces Stryker® N/A						
						1 per pack, Sterile
14-C0606	Straight Cross Cut - Carbide	mm	44.5	1.2	4.1	6
		in	1.75	0.05	0.16	
Replaces Stryker® 277-10-212, 277-010-212						
						1 per pack, Sterile
14-C0608	Tapered Cross Cut - Carbide	mm	44.5	1.6	4.3	6
		in	1.75	0.06	0.17	
Replaces Stryker® 277-10-216, 277-010-216						
						1 per pack, Sterile
14-C0616	Tapered Cross Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces Stryker® 277-10-221, 277-010-221						
						1 per pack, Sterile
16-D0130	Round - Diamond	mm	54.5	2.0	2.0	Fine
		in	2.15	0.08	0.08	
Replaces Stryker® 1608-6-87, 1608-006-087						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPST™, F1, TPX® & Command® 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-D0131	Round - Diamond	mm	54.5	2.4	2.4	Fine
		in	2.15	0.09	0.09	
Replaces Stryker® 1608-6-89, 1608-006-089						
						
1 per pack, Sterile						
16-D0132	Round - Diamond	mm	54.5	3.0	3.0	Medium
		in	2.15	0.11	0.11	
Replaces Stryker® 1608-6-91, 1608-006-091						
						
1 per pack, Sterile						
16-D0134	Round - Diamond	mm	54.5	3.5	3.5	Medium
		in	2.15	0.14	0.14	
Replaces Stryker® 1608-6-93, 1608-006-093						
						
1 per pack, Sterile						
16-D0135	Round - Diamond	mm	54.5	4.0	4.0	Medium
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-95, 1608-006-095						
						
1 per pack, Sterile						
16-D0137	Round - Diamond	mm	54.5	5.0	5.0	Medium
		in	2.15	0.20	0.20	
Replaces Stryker® 1608-6-99, 1608-006-099						
						
1 per pack, Sterile						
14-S0402	Cylinder Stainless Steel	mm	44.5	6.1	11.1	14
		in	1.75	0.24	0.44	
Replaces Stryker® 1607-2-1, 1607-002-001						
						
1 per pack, Sterile						

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®






Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0402	Cylinder Stainless Steel	mm	54.5	6.1	11.1	14
		in	2.15	0.24	0.44	
Replaces Stryker® 1608-2-1, 1608-002-001						
 1 per pack, Sterile						
15-S1101	Flame - Stainless Steel	mm	48.0	6.0	12.7	14
		in	1.89	0.24	0.50	
Replaces Stryker® 1607-2-47, 1607-002-047						
 1 per pack, Sterile						
14-S0201	Inverted Cone Stainless Steel	mm	44.5	3.1	3.3	10
		in	1.75	0.12	0.13	
Replaces Stryker® 1607-2-57, 1607-002-057						
 1 per pack, Sterile						
16-S1805	Lindemann Stainless Steel	mm	54.5	1.6	9.0	4
		in	2.15	0.06	0.35	
Replaces Stryker® N/A						
 1 per pack, Sterile						
16-S1801	Lindemann Stainless Steel	mm	54.5	2.3	9.5	4
		in	2.15	0.09	0.37	
Replaces Stryker® 1608-2-21, 1608-002-021						
 1 per pack, Sterile						
19-S1802	Lindemann Stainless Steel	mm	66.7	2.2	24.3	4
		in	2.63	0.09	0.96	
Replaces Stryker® 1608-2-43, 1608-002-043						
 1 per pack, Sterile						

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements






Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPST™, F1, TPX® & Command® 2

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
23-S1809	Lindemann Stainless Steel	mm	70.0	2.3	22.0	4
		in	2.76	0.09	0.87	
Replaces Stryker® N/A						
						
1 per pack, Sterile						
22-S1808	Lindemann Stainless Steel	mm	75.0	2.3	35.0	4
		in	2.95	0.09	1.38	
Replaces Stryker® N/A						
						
1 per pack, Sterile						
14-S0903	Oval - Stainless Steel	mm	44.5	4.0	9.5	10
		in	1.75	0.16	0.37	
Replaces Stryker® 1607-2-35, 1607-002-035						
						
1 per pack, Sterile						
14-S0909	Oval - Stainless Steel	mm	44.5	6.0	9.5	14
		in	1.75	0.24	0.37	
Replaces Stryker® 1607-2-3, 1607-002-003						
						
1 per pack, Sterile						
16-S0903	Oval - Stainless Steel	mm	54.5	4.0	9.5	10
		in	2.15	0.16	0.37	
Replaces Stryker® 1608-2-35, 1608-002-035						
						
1 per pack, Sterile						
16-S0909	Oval - Stainless Steel	mm	54.5	6.0	9.5	14
		in	2.15	0.24	0.37	
Replaces Stryker® 1608-2-3, 1608-002-003						
						
1 per pack, Sterile						

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-S1201	Pear - Stainless Steel	mm	44.5	6.1	10.0	14
		in	1.75	0.24	0.39	
Replaces Stryker® 1607-2-5, 1607-002-005						
						
1 per pack, Sterile						
16-S1201	Pear - Stainless Steel	mm	54.5	6.1	10.0	14
		in	2.15	0.24	0.39	
Replaces Stryker® 1608-2-5, 1608-002-005						
						
1 per pack, Sterile						
16-S0106	Round - Stainless Steel	mm	54.5	2.3	2.3	6
		in	2.15	0.09	0.09	
Replaces Stryker® 1608-6-147, 1608-006-147						
						
1 per pack, Sterile						
16-S0107	Round - Stainless Steel	mm	54.5	3.2	3.2	6
		in	2.15	0.126	0.126	
Replaces Stryker® 1608-6-149, 1608-006-149						
						
1 per pack, Sterile						
16-S0108	Round - Stainless Steel	mm	54.5	4.0	4.0	8
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-137, 1608-006-137						
						
1 per pack, Sterile						
16-S0109	Round - Stainless Steel	mm	54.5	5.1	5.1	8
		in	2.15	0.20	0.20	
Replaces Stryker® 1608-6-141, 1608-006-141						
						
1 per pack, Sterile						

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements





Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPST™, F1, TPX® & Command® 2

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0110	Round - Stainless Steel	mm	54.5	6.1	6.1	8
		in	2.15	0.24	0.24	
Replaces Stryker® 1608-6-143, 1608-006-143						
						1 per pack, Sterile
16-S0115	Round - Stainless Steel	mm	54.5	2.4	2.4	8
		in	2.15	0.095	0.095	
Replaces Stryker® 1608-2-9, 1608-002-009						
						1 per pack, Sterile
16-S0116	Round - Stainless Steel	mm	54.5	3.0	3.0	20
		in	2.15	0.12	0.12	
Replaces Stryker® 1608-6-163, 1608-006-163						
						1 per pack, Sterile
16-S0117	Round - Stainless Steel	mm	54.5	3.2	3.2	8
		in	2.15	0.126	0.126	
Replaces Stryker® 1608-2-11, 1608-002-011						
						1 per pack, Sterile
16-S0118	Round - Stainless Steel	mm	54.5	4.0	4.0	24
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-165, 1608-006-165						
						1 per pack, Sterile
16-S0119	Round - Stainless Steel	mm	54.5	4.0	4.0	6
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-151, 1608-006-151						
						1 per pack, Sterile

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0120	Round - Stainless Steel	mm	54.5	4.0	4.0	12
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-2-55, 1608-002-055						
						
1 per pack, Sterile						
16-S0121	Round - Stainless Steel	mm	54.5	4.6	4.6	6
		in	2.15	0.18	0.18	
Replaces Stryker® 1608-6-153, 1608-006-153						
						
1 per pack, Sterile						
16-S0122	Round - Stainless Steel	mm	54.5	4.6	4.6	8
		in	2.15	0.18	0.18	
Replaces Stryker® 1608-6-139, 1608-006-139						
						
1 per pack, Sterile						
16-S0123	Round - Stainless Steel	mm	54.5	4.5	4.5	12
		in	2.15	0.177	0.177	
Replaces Stryker® 1608-2-13, 1608-002-013						
						
1 per pack, Sterile						
16-S0126	Round - Stainless Steel	mm	54.5	7.1	7.1	8
		in	2.15	0.28	0.28	
Replaces Stryker® 1608-6-145, 1608-006-145						
						
1 per pack, Sterile						
16-S1401	Shannon Style Tapered - Stainless Steel	mm	54.5	2.1	11.1	3
		in	2.15	0.08	0.44	
Replaces Stryker® 1608-2-101, 1608-002-101						
						
1 per pack, Sterile						

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPS™, F1, TPX® & Command® 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements




Compatible with Stryker® Micro Drill Systems: RemB®, CORE™, TPST™, F1, TPX® & Command® 2

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S1601	Swanson Reamer Stainless Steel	mm	56.5	2.0	10.5	5
		in	2.23	0.08	0.41	
Replaces Stryker® 1608-2-131, 1608-002-131						
						
1 per pack, Sterile						
17-S1602	Swanson Reamer Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
Replaces Stryker® 1608-2-133, 1608-002-133						
						
1 per pack, Sterile						
17-S1603	Swanson Reamer Stainless Steel	mm	62.0	4.0	13.5	6
		in	2.44	0.16	0.53	
Replaces Stryker® 1608-2-135, 1608-002-135						
						
1 per pack, Sterile						
18-S1701	Wire Pass Drill Stainless Steel	mm	64.5	1.0	10.0	2
		in	2.54	0.039	0.39	
Replaces Stryker® 1608-2-63, 1608-002-063						
						
1 per pack, Sterile						
18-S1703	Wire Pass Drill Stainless Steel	mm	65.5	1.5	6.9	2
		in	2.57	0.059	0.27	
Replaces Stryker® 1608-2-59, 1608-002-059						
						
1 per pack, Sterile						
19-S1705	Wire Pass Drill Stainless Steel	mm	66.7	2.0	6.4	2
		in	2.63	0.079	0.25	
Replaces Stryker® 1608-2-61, 1608-002-061						
						
1 per pack, Sterile						

BUSA® Micro Burs - Oral Surgery

Stryker® Flat Shank Replacements

Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE™ & Oral Max

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS958N	Round - Carbide	mm	58.4	2.3	2.3	8
		in	2.30	0.09	0.09	
Replaces Stryker® 2296-101-423						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS94701N	Tapered Fissure Cross Cut - Carbide	mm	44.5	1.2	4.2	6
		in	1.75	0.05	0.17	
Replaces Stryker® 2296-100-212						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS94702N	Tapered Fissure Cross Cut - Carbide	mm	44.5	1.6	4.4	6
		in	1.75	0.06	0.17	
Replaces Stryker® 2296-100-216						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS94703N	Tapered Fissure Cross Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces Stryker® 2296-100-221						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS95701N	Tapered Fissure Cross Cut - Carbide	mm	58.4	1.2	4.2	6
		in	2.30	0.05	0.17	
Replaces Stryker® 2296-101-212						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS95702N	Tapered Fissure Cross Cut - Carbide	mm	58.4	1.6	4.4	6
		in	2.30	0.06	0.17	
Replaces Stryker® N/A						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	

STRYKER® | Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE™ & Oral Max

BUSA® Micro Burs - Oral Surgery

Stryker® Flat Shank Replacements






Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE™ & Oral Max

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS94703LN	Tapered Fissure Long Cross Cut - Carbide	mm	44.5	2.1	7.5	6
		in	1.75	0.08	0.30	
Replaces Stryker® 2296-100-321						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS95702LN	Tapered Fissure Long Cross Cut - Carbide	mm	58.4	1.6	6.0	6
		in	2.30	0.06	0.24	
Replaces Stryker® 2296-101-316						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS941703LN	Tapered Fluted Carbide	mm	44.5	2.1	7.0	8
		in	1.75	0.08	0.27	
Replaces Stryker® 2296-100-421						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KMS951703LN	Tapered Fluted Carbide	mm	58.4	2.1	7.0	8
		in	2.30	0.08	0.27	
Replaces Stryker® 2296-101-421						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	
KM94901N	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces Stryker® 2296-100-26, 2296-100-026						
					Non-Sterile	
KMS95901N	Oval - Carbide	mm	58.4	4.0	8.0	8
		in	2.30	0.16	0.31	
Replaces Stryker® 2296-101-26, 2296-101-026						
					5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)	

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
03-C0901	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces ConMed® Linvatec® 5091-36, 5091-136, 5091-236, M40B						
						
03-C0905	Oval - Carbide	mm	48.0	5.5	10.2	10
		in	1.89	0.22	0.40	
Replaces ConMed® Linvatec® 5091-38, 5091-138, 5091-238, M55B						
						
09-C0901	Oval - Carbide	mm	70.0	4.0	8.0	8
		in	2.76	0.16	0.31	
Replaces ConMed® Linvatec® 5092-36, 5092-136, 5092-236, L40B						
						
09-C0905	Oval - Carbide	mm	70.0	5.5	10.2	10
		in	2.76	0.22	0.40	
Replaces ConMed® Linvatec® 5092-38, 5092-138, 5092-238, L55B						
						
12-C0901	Oval - Carbide	mm	95.0	4.0	8.0	8
		in	3.74	0.16	0.31	
Replaces ConMed® Linvatec® 5093-36, 5093-136, 5093-236, X40B						
						
12-C0905	Oval - Carbide	mm	95.0	5.5	10.2	10
		in	3.74	0.22	0.40	
Replaces ConMed® Linvatec® 5093-38, 5093-138, 5093-238, X55B						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements


Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire®
Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0101	Round - Carbide	mm	51.0	1.0	1.0	6
		in	2.01	0.04	0.04	
Replaces ConMed® Linvatec® 5091-20, 5091-120, 5091-220						
						
04-C0103	Round - Carbide	mm	51.0	1.5	1.5	6
		in	2.01	0.06	0.06	
Replaces ConMed® Linvatec® 5091-22, 5091-122, 5091-222, M15R						
						
04-C0104	Round - Carbide	mm	51.0	1.8	1.8	6
		in	2.01	0.07	0.07	
Replaces ConMed® Linvatec® 5091-21, 5091-121, 5091-221, M18R						
						
04-C0105	Round - Carbide	mm	51.0	2.0	2.0	6
		in	2.01	0.08	0.08	
Replaces ConMed® Linvatec® 5091-24, 5091-124, 5091-224, M20R						
						
04-C0106	Round - Carbide	mm	51.0	2.3	2.3	8
		in	2.01	0.09	0.09	
Replaces ConMed® Linvatec® 5091-23, 5091-123, 5091-223, M23R, 5091-25, 5091-125, 5091-225, M24R						
						
04-C0107	Round - Carbide	mm	51.0	3.0	3.0	8
		in	2.01	0.12	0.12	
Replaces ConMed® Linvatec® 5091-26, 5091-126, 5091-226, M30R						
						

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems


ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0108	Round - Carbide	mm	51.0	4.0	4.0	8
		in	2.01	0.16	0.16	
Replaces ConMed® Linvatec® 5091-28, 5091-128, 5091-228, M40R						
						
04-C0109	Round - Carbide	mm	51.0	5.0	5.0	10
		in	2.01	0.20	0.20	
Replaces ConMed® Linvatec® 5091-30, 5091-130, 5091-230, M50R						
						
04-C0110	Round - Carbide	mm	51.0	6.0	6.0	8
		in	2.01	0.24	0.24	
Replaces ConMed® Linvatec® M60R						
						
04-C0170	Round - Carbide	mm	51.0	6.5	6.5	10
		in	2.01	0.26	0.26	
Replaces ConMed® Linvatec® 5091-32, 5091-132, 5091-232, M65R						
						
09-C0109	Round - Carbide	mm	70.0	5.0	5.0	10
		in	2.76	0.20	0.20	
Replaces ConMed® Linvatec® 5092-30, 5092-130, 5092-230, L50R						
						
09-C0110	Round - Carbide	mm	70.0	6.0	6.0	8
		in	2.76	0.24	0.24	
Replaces ConMed® Linvatec® L60R						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements


Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®,
MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire®
Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-C0170	Round - Carbide	mm	70.0	6.5	6.5	10
		in	2.76	0.26	0.26	
Replaces ConMed® Linvatec® 5092-32, 5092-132, 5092-232, L65R						
						
10-C0101	Round - Carbide	mm	71.5	1.0	1.0	6
		in	2.81	0.04	0.04	
Replaces ConMed® Linvatec® 5092-20, 5092-120, 5092-220, L10R						
						
10-C0103	Round - Carbide	mm	72.0	1.5	1.5	6
		in	2.83	0.06	0.06	
Replaces ConMed® Linvatec® 5092-22, 5092-122, 5092-222, L15R						
						
10-C0105	Round - Carbide	mm	72.5	2.0	2.0	6
		in	2.85	0.08	0.08	
Replaces ConMed® Linvatec® 5092-24, 5092-124, 5092-224, L20R						
						
10-C0107	Round - Carbide	mm	73.0	3.0	3.0	8
		in	2.87	0.12	0.12	
Replaces ConMed® Linvatec® 5092-26, 5092-126, 5092-226, L30R						
						
10-C0108	Round - Carbide	mm	73.0	4.0	4.0	8
		in	2.87	0.16	0.16	
Replaces ConMed® Linvatec® 5092-28, 5092-128, 5092-228, L40R						
						

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
12-C0105	Round - Carbide	mm	95.0	2.0	2.0	6
		in	3.74	0.08	0.08	
Replaces ConMed® Linvatec® 5093-24, 5093-124, 5093-224, X20R						
						
12-C0108	Round - Carbide	mm	95.0	4.0	4.0	8
		in	3.74	0.16	0.16	
Replaces ConMed® Linvatec® 5093-28, 5093-128, 5093-228, X40R						
						
12-C0109	Round - Carbide	mm	95.0	5.0	5.0	10
		in	3.74	0.20	0.20	
Replaces ConMed® Linvatec® 5093-30, 5093-130, 5093-230, X50R						
						
04-C0601	Tapered - Carbide	mm	51.0	0.9	5.2	6
		in	2.01	0.035	0.20	
Replaces ConMed® Linvatec® 5091-15, 5091-115, 5091-215						
						
04-C0606	Tapered - Carbide	mm	51.0	1.2	4.2	6
		in	2.01	0.05	0.17	
Replaces ConMed® Linvatec® M12T						
						
04-C0623	Tapered - Carbide	mm	51.0	1.2	6.0	6
		in	2.01	0.05	0.24	
Replaces ConMed® Linvatec® N/A						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements


Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®,
MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire®
Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0607	Tapered - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces ConMed® Linvatec® 5091-18, 5091-118, 5091-218						
						
04-C0609	Tapered - Carbide	mm	51.0	1.5	10.0	6
		in	2.01	0.06	0.39	
Replaces ConMed® Linvatec® 5091-02, 5091-102, 5091-202, M15T						
						
04-C0612	Tapered Cross Cut - Carbide	mm	51.0	2.1	4.9	6
		in	2.01	0.08	0.19	
Replaces ConMed® Linvatec® 5091-13, 5091-113, 5091-213, M18T						
						
04-C0615	Tapered - Carbide	mm	51.0	2.1	4.8	6
		in	2.01	0.08	0.19	
Replaces ConMed® Linvatec® N/A						
						
04-C0617	Tapered - Carbide	mm	51.0	2.0	7.5	6
		in	2.01	0.08	0.30	
Replaces ConMed® Linvatec® 5091-06, 5091-106, 5091-206, M27T						
						
06-C0619	Tapered - Carbide	mm	57.0	2.0	19.0	6
		in	2.24	0.08	0.75	
Replaces ConMed® Linvatec® 5091-14, 5091-114, 5091-214, M20T						
						

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems


ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0617	Tapered - Carbide	mm	73.0	2.0	7.5	6
		in	2.87	0.08	0.30	
Replaces ConMed® Linvatec® 5092-06, 5092-106, 5092-206, L27T						
						
11-C0619	Tapered - Carbide	mm	80.0	2.0	19.0	6
		in	3.15	0.08	0.75	
Replaces ConMed® Linvatec® 5092-14, 5092-114, 5092-214, L20T						
						
03-C0312	Straight Cross Cut - Carbide	mm	48.0	4.0	10.0	6
		in	1.89	0.16	0.39	
Replaces ConMed® Linvatec® 5091-54, 5091-154, 5091-254, M40S						
						
04-C0302	Straight Cross Cut - Carbide	mm	51.0	1.0	5.0	6
		in	2.01	0.04	0.20	
Replaces ConMed® Linvatec® 5091-01, 5091-101, 5091-201, M10S						
						
04-C0306	Straight Cross Cut - Carbide	mm	51.0	1.5	7.5	6
		in	2.01	0.06	0.30	
Replaces ConMed® Linvatec® 5091-04, 5091-104, 5091-204, M15S						
						
04-C0307	Straight - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces ConMed® Linvatec® 5091-17, 5091-117, 5091-217, M17S						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements





Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®,
MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire®
Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0306	Straight - Carbide	mm	73.0	1.5	7.5	6
		in	2.87	0.06	0.30	
Replaces ConMed® Linvatec® 5092-04, 5092-104, 5092-204, L13S, L17S						
						
10-C0310	Straight Cross Cut - Carbide	mm	73.0	2.1	10.0	8
		in	2.87	0.08	0.39	
Replaces ConMed® Linvatec® 5092-52, 5092-152, 5092-252, L20S						
						
03-D0135	Round - Diamond	mm	48.0	4.0	4.0	Medium
		in	1.89	0.16	0.16	
Replaces ConMed® Linvatec® 5091-264, 5091-164, 5091-264, M40D						
						
04-D0139	Round - Diamond	mm	51.0	2.0	2.0	Medium
		in	2.01	0.08	0.08	
Replaces ConMed® Linvatec® 5091-60, 5091-160, 5091-260, M20D						
						
04-D0133	Round - Diamond	mm	51.0	3.0	3.0	Medium
		in	2.01	0.12	0.12	
Replaces ConMed® Linvatec® 5091-62, 5091-162, 5091-262						
						
09-D0135	Round - Diamond	mm	70.0	4.0	4.0	Medium
		in	2.76	0.16	0.16	
Replaces ConMed® Linvatec® 5092-64, 5092-164, 5092-264, L40D, L60D						
						

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-D0141	Round - Diamond	mm	70.0	5.0	5.0	Coarse
		in	2.76	0.20	0.20	
Replaces ConMed® Linvatec® 5092-66, 5092-166, 5092-266, L50D						
						
09-D0138	Round - Diamond	mm	70.0	6.0	6.0	Coarse
		in	2.76	0.24	0.24	
Replaces ConMed® Linvatec® 5092-68, 5092-168, 5092-268						
						
10-D0139	Round - Diamond	mm	73.0	2.0	2.0	Medium
		in	2.87	0.08	0.08	
Replaces ConMed® Linvatec® 5092-60, 5092-160, 5092-260, L20D						
						
10-D0133	Round - Diamond	mm	73.0	3.0	3.0	Medium
		in	2.87	0.12	0.12	
Replaces ConMed® Linvatec® 5092-62, 5092-162, 5092-262, L30D						
						
08-S1802	Lindemann Stainless Steel	mm	67.0	2.3	25.0	4
		in	2.64	0.09	0.98	
Replaces ConMed® Linvatec® 5091-12, 5091-112, 5091-212, L17L						
						
09-S0115	Round - Stainless Steel	mm	70.0	2.4	2.4	8
		in	2.76	0.091	0.091	
Replaces ConMed® Linvatec® N/A						
						




CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®,
MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire®
Small Bone/Micro Drill Systems

ConMed® Linvatec®

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-S0120	Round - Stainless Steel	mm	70.0	4.0	4.0	12
		in	2.76	0.16	0.16	
Replaces ConMed® Linvatec® N/A						
						
12-S0153	Round - Stainless Steel	mm	95.0	3.1	3.1	10
		in	3.74	0.12	0.12	
Replaces ConMed® Linvatec® N/A						
						
04-S1501	Shannon Style - Tapered Stainless Steel	mm	51.0	2.3	12.0	4
		in	2.01	0.09	0.47	
Replaces ConMed® Linvatec® N/A						
						
06-S1601	Swanson Reamer Stainless Steel	mm	57.0	2.0	10.0	5
		in	2.24	0.08	0.39	
Replaces ConMed® Linvatec® 5091-42, 5091-142, 5091-242, M20F						
						
07-S1602	Swanson Reamer Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
Replaces ConMed® Linvatec® 5091-43, 5091-143, 5091-243, M30F						
						
07-S1605	Swanson Reamer Stainless Steel	mm	62.0	4.0	16.5	6
		in	2.44	0.16	0.65	
Replaces ConMed® Linvatec® 5091-44, 5091-144, 5091-244, M40F						
						

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-S1701	Wire Pass Drill Stainless Steel	mm	51.0	1.0	19.0	2
		in	2.01	0.04	0.75	
Replaces ConMed® Linvatec® 5091-47, 5091-147, 5091-247, M10W						
						
04-S1703	Wire Pass Drill Stainless Steel	mm	51.0	1.5	19.0	2
		in	2.01	0.059	0.75	
Replaces ConMed® Linvatec® 5091-48, 5091-148, 5091-248, M15W						
						
04-S1704	Wire Pass Drill Stainless Steel	mm	51.0	2.0	19.0	2
		in	2.01	0.079	0.75	
Replaces ConMed® Linvatec® 5091-49, 5091-149, 5091-249, M20W						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BMA-S90006S	N/A	N/A	2490-95 2490-95-000	N/A	N/A	Page 16
BMA-S90007S	N/A	N/A	86-4192 86-4192-000	N/A	N/A	Page 16
BMA-S90008S	7650-2038 7650-1038	N/A	N/A	N/A	N/A	Page 15
BMA-S90010S	7650-2038A 7650-1038A	N/A	N/A	N/A	N/A	Page 15
BMA-S90012S	6541-4-003 6541-004-003	N/A	N/A	N/A	N/A	Page 15
BR900-722	N/A	N/A	900.722	N/A	N/A	Page 17
KM166-00-00	N/A	N/A	1801-10 5604-60	00231801500 00231801506 00-2318-015-00 00-2318-015-06	8054-010NS	Page 7
KM166-00-01	277-082-085S7	N/A	1801-12 5604-61	00231802000 00231802006 00-2318-020-00 00-2318-020-06	8054-011NS 8054-020NS	Page 7
KM166-00-02	277-082-086S5	N/A	1801-14 5604-62	00231802500 00231802506 00-2318-025-00 00-2318-025-06	8054-012NS	Page 7
KM166-00-03	277-082-085S5	N/A	1801-16 5604-63	00231802700 00231802706 00-2318-027-00 00-2318-027-06	8054-013NS 8054-021NS	Page 7
KM166-00-04	277-082-087S6	N/A	1801-18 5604-64	00231803200 00231803206 00-2318-032-00 00-2318-032-06	8054-014NS 8054-022NS	Page 7
KM166-00-05	0277-082-086	N/A	1801-20 5604-66	00231803500 00231803506 00-2318-035-00 00-2318-035-06	8054-015NS 8054-023NS	Page 7
KM166-00-06	N/A	N/A	5604-67 1801-22	00231804000 00231804006 00-2318-040-00 00-2318-040-06	8054-016NS	Page 7 Page 11
KM166-00-07	N/A	N/A	1801-24 5604-69	00231804800 00231804806 00-2318-048-00 00-2318-048-06	8054-018NS	Page 7
KM166-00-08	N/A	N/A	N/A	00231806400 00231806406 00-2318-064-00 00-2318-064-06	8054-019NS	Page 7
KM166-00-11	N/A	N/A	N/A	00-2318-011-00 00-2318-011-06 00230801100 00231801106	N/A	Page 7
KM166-00-45	N/A	N/A	N/A	00-4805-130 2318-045	8054-017NS	Page 7
KM166-00-60	N/A	N/A	N/A	N/A	N/A	Page 7
KM166-10-00	277-082-085S3	N/A	1801-10 5604-60	00231801500 00-2318-015-00	8054-010	Page 7

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM166-10-01	277-082-085S4	N/A	1801-12 5604-61	00231802000 00-2318-020-00	8054-011 8054-020	Page 7
KM166-10-02	277-082-085S1	N/A	1801-14 5604-62	00231802500 00-2318-025-00	8054-012	Page 7
KM166-10-03	277-082-085S2	N/A	1801-16 5604-63	00231802700 00-2318-027-00	8054-021 8054-013	Page 7
KM166-10-04	277-082-087S1	N/A	1801-18 5604-64	00231803200 00-2318-032-00	8054-014 8054-022	Page 7
KM166-10-05	277-082-086S1	N/A	1801-20 5604-66	00231803500 00-2318-035-00	8054-015 8054-923	Page 7
KM166-10-06	277-082-086S2	N/A	1801-22 5604-67	00231804000 00-2318-040-00	8054-016	Page 7
KM166-10-07	277-082-086S4	N/A	N/A	00231804800 00-2318-048-00	8054-018	Page 7
KM166-10-08	277-082-087S7	N/A	N/A	00231806400 00-2318-064-00	8054-019	Page 7
KM166-10-11	N/A	N/A	N/A	N/A	N/A	Page 7
KM166-10-45	N/A	N/A	N/A	N/A	8054-017	Page 7
KM166-10-60	N/A	N/A	N/A	N/A	N/A	Page 7
KM166-310-11	N/A	N/A	310.11	00241001100 00-2410-011-00 00-4806-060-11 2410-11	8053-011NS	Page 8
KM166-310-11S	N/A	N/A	310.11S	N/A	8053-011	Page 8
KM166-310-15	N/A	N/A	310.15	00241001500 00-2410-015-00 00-4806-085-15 2410-15	8053-015NS	Page 8
KM166-310-15S	N/A	N/A	310.15S	N/A	8053-015	Page 8
KM166-310-16	N/A	N/A	310.16	00241001600 00-2410-016-00 2410-16	8053-016NS	Page 8
KM166-310-16S	N/A	N/A	310.16S	N/A	8053-016	Page 8
KM166-310-19	N/A	N/A	310.19	00241002000 00-2410-020-00 00-4806-100-20 2410-20	8053-019NS	Page 8
KM166-310-19S	N/A	N/A	310.19S	N/A	8053-019	Page 8
KM166-310-21	N/A	N/A	310.21	00241002100 00-2410-021-00 00-4806-125-20 2410-21	8053-021NS	Page 8
KM166-310-21S	N/A	N/A	310.21S	N/A	8053-021	Page 8
KM166-310-221	702453	N/A	310.221	N/A	N/A	Page 10
KM166-310-23	700351	N/A	310.23	00241002600 00-2410-026-00 2410-26	8053-023NS	Page 8
KM166-310-23S	N/A	N/A	310.23S	N/A	8053-023	Page 8

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM166-310-25	N/A	N/A	310.25	00241002500 00-2410-025-00 00-4806-110-25 2410-25	8053-025NS	Page 8
KM166-310-25S	N/A	N/A	310.25S	N/A	8053-025	Page 8
KM166-310-26	N/A	N/A	310.26	00241002700 00-2410-027-00 00-4806-100-27 2410-27	8053-026NS	Page 8
KM166-310-26S	N/A	N/A	310.26S	N/A	8053-026	Page 8
KM166-310-28	N/A	N/A	310.28	00241002800 00-2410-028-00 2410-28	8053-028NS	Page 8
KM166-310-288	N/A	N/A	310.288	N/A	N/A	Page 9
KM166-310-28S	N/A	N/A	310.28S	N/A	8053-028	Page 8
KM166-310-29	N/A	N/A	310.29	00241003300 00-2410-033-00 00-4806-195-32 2410-33	8053-029NS	Page 8
KM166-310-29S	N/A	N/A	310.29S	N/A	8053-029	Page 8
KM166-310-31	700358	N/A	310.31	00241003200 00-2410-032-00 00-4806-145-32 2410-32	8053-031NS	Page 8
KM166-310-31S	N/A	N/A	310.31S	N/A	8053-031	Page 8
KM166-310-35	N/A	N/A	310.35	00241003500 00-2410-035-00 00-4806-110-35 2410-35	8053-035NS	Page 8
KM166-310-35S	N/A	N/A	310.35S	N/A	8053-035	Page 8
KM166-310-37	N/A	N/A	310.37	00241003600 00-2410-036-00 00-4806-195-35 2410-36	8053-037NS	Page 9
KM166-310-37S	N/A	N/A	310.37S	N/A	8053-037	Page 9
KM166-310-44	700359	N/A	310.44	00241004500 00-2410-045-00 00-4806-145-45 2410-45	8053-044NS	Page 9
KM166-310-44S	N/A	N/A	310.44S	N/A	8053-044	Page 9
KM166-310-48	N/A	N/A	310.48	00241004600 00-2410-046-00 00-4806-195-45 2410-46	8053-048NS	Page 9
KM166-310-48S	N/A	N/A	310.48S	N/A	8053-048	Page 9
KM166-310-509	N/A	N/A	310.509 357.509	N/A	N/A	Page 9
KM166-310-530	N/A	N/A	310.530	N/A	N/A	Page 9
KM166-310-530S	N/A	N/A	310.530S	N/A	N/A	Page 9
KM166-310-65	N/A	N/A	310.65	N/A	N/A	Page 10
KM166-310-66	N/A	N/A	310.66	N/A	N/A	Page 10

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM166-310-67	702449	N/A	310.67	N/A	N/A	Page 10
KM166-310-68	702450	N/A	310.68	N/A	N/A	Page 10
KM166-315-28	N/A	N/A	315.28	N/A	N/A	Page 9
KM166-323-062	N/A	N/A	323.62	N/A	N/A	Page 9
KM166-324-214	N/A	N/A	324.214	N/A	N/A	Page 9
KM166-357-407	N/A	N/A	357.407	N/A	N/A	Page 9
KM168-19-18	3811-4-090	N/A	1627-18-000 1653-32	00-0187-004-59 187-04-59	1632-509NS	Page 13
KM168-19-18S	N/A	N/A	N/A	47-0187-004-59 187-07-59	1632-509	Page 13
KM168-19-31	3811-7-090	N/A	1630-18-000 1653-62		1648-509NS	Page 13
KM168-19-31S	N/A	N/A	N/A	47-0187-007-59	1648-509	Page 13
KM168-19-33	3811-2-090	N/A	1628-18-000 1653-12	00-0187-002-59 187-02-59	1624-509NS	Page 13
KM168-19-33S	N/A	N/A	N/A	47-0187-002-59	1624-509	Page 13
KM168-19-53	3811-6-090	N/A	1629-18-000 1653-53	00-0187-006-69 187-06-59	1640-509NS	Page 13
KM168-19-53S	N/A	N/A	N/A	47-0187-006-59	1640-509	Page 13
KM168-19-56	3811-1-090	N/A	1624-18-000 1653-02	00-0187-001-59 187-01-59	1620-509NS	Page 13
KM168-19-56S	N/A	N/A	N/A	47-0187-001-59	1620-509	Page 13
KM168-19-76	3811-3-090	N/A	1626-18-000 1653-22	00-0187-003-59 187-03-59	1628-509NS	Page 13
KM168-19-76S	N/A	N/A	N/A	47018700359	1628-509	Page 13
KM168-19-96	3811-5-090	N/A	1628-18-000 1653-42	00-0187-005-59 187-05-59	1636-509NS	Page 13
KM168-19-96S	N/A	N/A	N/A	N/A	1636-509	Page 13
KM168-29-18	3711-4-090	N/A	1627-38-000 1653-30	187-04-19	1632-109NS	Page 13
KM168-29-18S	5500-19-32	N/A	N/A	N/A	1632-109	Page 13
KM168-29-31	3711-7-090	N/A	1630-38-000 1653-60	187-07-19	1648-109NS	Page 13
KM168-29-31S	5500-19-48	N/A	N/A	N/A	1648-109	Page 13
KM168-29-33	3711-2-090	N/A	1625-38-000 1653-10	187-02-19	1624-109NS	Page 13
KM168-29-33S	5500-19-24	N/A	N/A	N/A	1624-109	Page 13
KM168-29-53	3711-6-090	N/A	1629-38-000 1653-50	187-06-19	1640-109NS	Page 13
KM168-29-53S	5500-19-40	N/A	N/A	N/A	1640-109	Page 13
KM168-29-56	3711-1-090	N/A	1624-38-000 1653-00	187-01-19	1620-109NS	Page 13
KM168-29-56S	5500-19-20	N/A	N/A	N/A	1620-109	Page 13
KM168-29-76	3711-3-090	N/A	1626-38-000 1653-20	187-03-19	1628-109NS	Page 13

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM168-29-76S	5500-19-28	N/A	N/A	N/A	1628-109	Page 13
KM168-29-96	3711-5-090	N/A	1628-38-000 1653-40	187-05-19	1636-109NS	Page 13
KM168-29-96S	5500-19-36	N/A	N/A	N/A	1636-109	Page 13
KM168-39-18	3810-4-090	N/A	1634-18-000 1653-35	00-0187-004-69 187-04-69 27-361692	1632-609NS	Page 13
KM168-39-18S	N/A	N/A	N/A	N/A	1632-609	Page 13
KM168-39-31	3810-7-090	N/A	1637-18-000 1653-64	00-0187-007-69 187-07-69	1648-609NS	Page 13
KM168-39-31S	N/A	N/A	N/A	N/A	1648-609	Page 13
KM168-39-33	3810-2-090	N/A	1632-18-000 1653-14	00-0187-002-69 187-02-69 27-361292	1624-609NS	Page 13
KM168-39-33S	N/A	N/A	N/A	N/A	1624-609	Page 13
KM168-39-53	3810-6-090	N/A	1636-18-000 1653-54	187-06-69 27-362092	1640-609NS	Page 13
KM168-39-53S	N/A	N/A	N/A	N/A	1640-609	Page 13
KM168-39-56	3810-1-090	N/A	1631-18-000 1653-04	00-0187-001-69 187-01-69 27-361092	1620-609NS	Page 13
KM168-39-56S	N/A	N/A	N/A	N/A	1620-609	Page 13
KM168-39-76	3810-3-090	N/A	1633-18-000 1653-24	00-0187-003-69 187-03-69 27-361492	1628-609NS	Page 13
KM168-39-76S	N/A	N/A	N/A	N/A	1628-609	Page 13
KM168-39-96	3810-5-090	N/A	1635-18-000 1653-44	00-0187-005-69 187-05-69 27-361892	1636-609NS	Page 13
KM168-39-96S	N/A	N/A	N/A	47-0187-005-69	1636-609	Page 13
KM168-49-18	3710-4-090	N/A	1631-23-000 1653-36	187-04-79	1632-709NS	Page 13
KM168-49-18S	5500-219-32	N/A	N/A	N/A	1632-709	Page 13
KM168-49-31	3710-7-090	N/A	1631-26-000 1653-66	187-07-79	1648-709NS	Page 13
KM168-49-31S	5500-219-48	N/A	N/A	N/A	1648-709	Page 13
KM168-49-33	3710-2-090	N/A	1631-21-000 1653-16	187-02-79	1624-709NS	Page 13
KM168-49-33S	5500-219-24	N/A	N/A	N/A	1624-709	Page 13
KM168-49-53	3710-6-090	N/A	1631-25-000 1653-56	187-06-79	1640-709NS	Page 13
KM168-49-53S	5500-219-40	N/A	N/A	N/A	1640-709	Page 13
KM168-49-56	3710-1-090	N/A	1631-20-000 1653-06	187-01-79	1620-709NS	Page 13
KM168-49-56S	5500-219-20	N/A	N/A	N/A	1620-709	Page 13
KM168-49-76	3710-3-090	N/A	1631-22-000 1653-26	187-03-79	1628-709NS	Page 13

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM168-49-76S	5500-219-28	N/A	N/A	N/A	1628-709	Page 13
KM168-49-96	3710-5-090	N/A	1631-24-000 1653-46	187-05-79	1636-709NS	Page 13
KM168-49-96S	5500-219-36	N/A	N/A	N/A	1636-709	Page 13
KM169-19-18	3593-4-090	N/A	1614-18-000 1653-33	2620-23	1632-509TNS	Page 14
KM169-19-18S	N/A	N/A	N/A	N/A	1632-509T	Page 14
KM169-19-31	3593-7-090	N/A	1611-17-000 1611-58-000 1617-18 1653-63	2620-13	1648-509TNS	Page 14
KM169-19-31S	N/A	N/A	N/A	47-2620-013-00	1648-509T	Page 14
KM169-19-33	3593-2-090	N/A	1612-18-000 1653-13	2620-29	1624-509TNS	Page 14
KM169-19-33S	N/A	N/A	N/A	N/A	1624-509T	Page 14
KM169-19-53	3593-6-090	N/A	1616-18-000 1653-53	2620-17	1640-509TNS	Page 14
KM169-19-53S	N/A	N/A	N/A	47-2620-017-00	1640-509T	Page 14
KM169-19-56	3593-1-090	N/A	1611-18-000 1653-03	00-2620-033-00 2620-33	1620-509TNS	Page 14
KM169-19-56S	N/A	N/A	N/A	47-2620-033-00	1620-509T	Page 14
KM169-19-76	3593-3-090	N/A	1613-18-000 1653-23	2620-25	1628-509TNS	Page 14
KM169-19-76S	N/A	N/A	N/A	47-2620-025-00	1628-509T	Page 14
KM169-19-96	3593-5-090	N/A	1615-18-000 1653-43	2620-19	1636-509TNS	Page 14
KM169-19-96S	N/A	N/A	N/A	N/A	1636-509T	Page 14
KM169-29-18	3598-4-090	N/A	1617-63-000 1653-31	2620-21	1632-109TNS	Page 14
KM169-29-18S	5500-119-32	N/A	N/A	N/A	1632-109T	Page 14
KM169-29-31	3598-7-090	N/A	1617-66-000 1653-61	2620-11	1648-109TNS	Page 14
KM169-29-31S	5500-119-48	N/A	N/A	N/A	1648-109T	Page 14
KM169-29-33	3598-2-090	N/A	1617-61-000 1653-11	2620-27	1624-109TNS	Page 14
KM169-29-33S	5500-119-24	N/A	N/A	N/A	1624-109T	Page 14
KM169-29-53	3598-6-090	N/A	1617-65-000 1653-51	2620-15	1640-109TNS	Page 14
KM169-29-53S	5500-119-40	N/A	N/A	N/A	1640-109T	Page 14
KM169-29-56	3598-1-090	N/A	1617-60-000 1653-01	2620-31	1620-109TNS	Page 14
KM169-29-56S	5500-119-20	N/A	N/A	N/A	1620-109T	Page 14
KM169-29-76	3598-3-090	N/A	1617-62-000 1653-21	2620-03	1628-109TNS	Page 14
KM169-29-76S	5500-119-28	N/A	N/A	N/A	1628-109T	Page 14

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM169-29-96	3598-5-090	N/A	1617-64-000 1653-41	2620-01	1636-109TNS	Page 14
KM169-29-96S	5500-119-36	N/A	N/A	N/A	1636-109T	Page 14
KM169-39-18	3594-4-090	N/A	1614-58-000 1617-73-000 1653-35	00-0262-004-04 262-04-04 27-363747	1632-609TNS	Page 14
KM169-39-18S	N/A	N/A	N/A	47-0262-004-04	1632-609T	Page 14
KM169-39-31	3594-7-090	N/A	1617-56-000 1617-76-000 1653-65	00-0262-001-04 262-01-04 27-363787	1648-609TNS	Page 14
KM169-39-31S	N/A	N/A	N/A	N/A	1648-609T	Page 14
KM169-39-33	3594-2-090	N/A	1612-58-000 1617-71-000 1653-15	00-0262-006-04 262-06-04 27-363727	1624-609TNS	Page 14
KM169-39-33S	N/A	N/A	N/A	47-0262-006-04	1624-609T	Page 14
KM169-39-53	3594-6-090	N/A	1616-56-000 1617-75-000 1653-55	00-0262-002-04 262-02-04	1640-609TNS	Page 14
KM169-39-53S	N/A	N/A	N/A	N/A	1640-609T	Page 14
KM169-39-56	3594-1-090	N/A	1611-58-000 1617-70-000 1653-05	00-0262-007-04 262-07-04 27-363717	1620-609TNS	Page 14
KM169-39-56S	N/A	N/A	N/A	N/A	1620-609T	Page 14
KM169-39-76	3594-3-090	N/A	1613-58-000 1617-72-000 1653-25	00-0262-005-04 262-05-04 27-363737	1628-609TNS	Page 14
KM169-39-76S	N/A	N/A	N/A	N/A	1628-609T	Page 14
KM169-39-96	3594-5-090	N/A	1615-58-000 1617-74-000 1653-45	00-0262-003-04 262-03-04 27-363757	1636-609TNS	Page 14
KM169-39-96S	N/A	N/A	N/A	47-0262-003-04	1636-609T	Page 14
KM169-49-18	3595-4-090	N/A	1617-73 1653-37	262-04-79	1632-709TNS	Page 14
KM169-49-18S	5500-319-32	N/A	N/A	N/A	1632-709T	Page 14
KM169-49-31	3595-7-090	N/A	1617-76 1653-67	262-01-79	1648-709TNS	Page 14
KM169-49-31S	5500-319-48	N/A	N/A	N/A	1648-709T	Page 14
KM169-49-33	3595-2-090	N/A	1617-71 1653-17	262-06-79	1624-709TNS	Page 14
KM169-49-33S	5500-319-24	N/A	N/A	N/A	1624-709T	Page 14
KM169-49-53	3595-6-090	N/A	1617-75 1653-57	262-02-79	1640-709TNS	Page 14
KM169-49-53S	5500-319-40	N/A	N/A	N/A	1640-709T	Page 14
KM169-49-56	3595-1-090	N/A	1617-70 1653-07	262-07-79	1620-709TNS	Page 14
KM169-49-56S	5500-319-20	N/A	N/A	N/A	1620-709T	Page 14
KM169-49-76	3595-3-090	N/A	1617-72 1653-27	262-05-79	1628-709TNS	Page 14

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM169-49-76S	5500-319-28	N/A	N/A	N/A	1628-709T	Page 14
KM169-49-96	3595-5-090	N/A	1617-74 1653-47	262-03-79	1636-709TNS	Page 14
KM169-49-96S	5500-319-36	N/A	N/A	N/A	1636-709T	Page 14
KM172-14-28	N/A	N/A	1650-00	186-04-54	1604-128NS	Page 11
KM172-14-28S	N/A	N/A	N/A	N/A	1604-128	Page 11
KM172-14-35	1346-001-000	N/A	1650-01	N/A	1604-135NS	Page 11
KM172-14-35S	N/A	N/A	N/A	N/A	1604-135	Page 11
KM172-14-45	1346-003-000	N/A	1650-02	N/A	1604-145NS	Page 11
KM172-14-45S	N/A	N/A	N/A	N/A	1604-145	Page 11
KM172-14-54	N/A	N/A	1650-03	N/A	1604-154NS	Page 11
KM172-14-54S	N/A	N/A	N/A	N/A	1604-154	Page 11
KM172-14-62	1346-005-000	N/A	1650-04	N/A	1604-162NS	Page 11
KM172-14-62S	N/A	N/A	N/A	1604-162	N/A	Page 11
KM172-15-35	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-15-45	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-155-28	N/A	N/A	5600-20	N/A	N/A	Page 11
KM172-15-54	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-155-54	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-155-62	N/A	N/A	5600-23	N/A	1600-004NS	Page 11
KM172-15-62	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-16-28	1346-026-000	N/A	N/A	N/A	N/A	Page 11
KM172-16-28S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-16-35	1346-002-000	N/A	N/A	N/A	N/A	Page 11
KM172-16-35S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-16-45	1346-004-000	N/A	292.10	186-60	N/A	Page 11
KM172-16-45S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-16-54	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-16-54S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-16-62	1346-006-000	N/A	292.16	N/A	N/A	Page 11
KM172-16-62S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-19-35	3813-1-090	N/A	1602-14	186-50 47-186-50	1600-9355NS	Page 11
KM172-19-35S	N/A	N/A	N/A	N/A	1600-9355	Page 11
KM172-19-45	3813-2-090	N/A	1604-14	186-66 47-186-66	1600-9455NS	Page 11
KM172-19-45S	N/A	N/A	1652-22-0000	186-66 47-186-66	1600-9455	Page 11
KM172-19-54	N/A	N/A	1606-14	N/A	1600-9545NS	Page 11
KM172-19-54S	N/A	N/A	N/A	N/A	1600-9545	Page 11

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM172-19-62	3813-3-090	N/A	1608-14	186-82	1600-9625NS	Page 11
KM172-19-62S	N/A	N/A	1652-42-0000	47-186-82	1600-9625	Page 11
KM172-24-28	1346-26	N/A	1642-00 1652-55	186-04-14	1600-428NS	Page 11
KM172-24-28S	5300-4-28	N/A	N/A	47-0186-004-14	1600-428	Page 11
KM172-24-35	1346-1 3715-1-040	N/A	1643-10 1652-65	N/A	1600-435NS	Page 11
KM172-24-35S	5300-4-35	N/A	N/A	N/A	1600-435	Page 11
KM172-24-45	1346-3 3715-2-040	N/A	1644-10 1652-75	N/A	1600-445NS	Page 11
KM172-24-45S	5300-4-45	N/A	N/A	N/A	1600-445	Page 11
KM172-24-54	N/A	N/A	1645-10 1652-85	N/A	1600-454NS	Page 11
KM172-24-54S	5300-4-54	N/A	N/A	N/A	1600-454	Page 11
KM172-24-62	1346-5 3715-3-040	N/A	1646-10 1652-95	N/A	1600-462NS	Page 11
KM172-24-62S	5300-4-62	N/A	N/A	N/A	1600-462	Page 11
KM172-26-28	N/A	N/A	1642-06-028 1652-59	N/A	1600-628NS	Page 11
KM172-26-28S	5300-6-28	5050-046 5050-246	N/A	N/A	1600-628	Page 11
KM172-26-35	1346-2	N/A	1642-06-035 1652-69	00-0186-044-00 186-44	1600-635NS	Page 11
KM172-26-35S	5300-6-35	5050-047 5050-247	N/A	47-0186-044-00 47-186-44	1600-635	Page 11
KM172-26-45	1346-4	N/A	292.50 1642-06-045 1652-79	00-0186-060-00	1600-645NS	Page 11
KM172-26-45S	5300-6-45	5050-048 5050-248	N/A	47-0186-060-00	1600-645	Page 11
KM172-26-54	N/A	N/A	1642-06-054 1652-89	N/A	1600-654NS	Page 11
KM172-26-54S	5300-6-54	N/A	N/A	N/A	1600-654	Page 12
KM172-26-62	1346-6	N/A	292.56 1642-06-062 1652-99	00-0186-076-00 186-76	1600-662NS	Page 11
KM172-26-62S	5300-6-62	5050-049 5050-249	N/A	47-0186-076-00 47-186-76	1600-662	Page 11
KM172-29-28	1346-9-1	N/A	N/A	1601-31-000 1638-32-000	N/A	Page 11
KM172-29-28S	5300-9-28	N/A	N/A	N/A	N/A	Page 11
KM172-29-35	1346-9-2 3713-1-090	N/A	1647-10 1652-12	00-0186-046-00	1600-935NS	Page 11
KM172-29-35S	5300-9-35	N/A	N/A	47-0186-046-00 47-186-46	1600-935	Page 11
KM172-29-45	1346-9-3 3713-2-090	N/A	1648-10-000 1652-22	00-0186-062-00 186-62	1600-945NS	Page 12

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM172-29-45S	5300-9-45	N/A	N/A	47-0186-062-00 47-186-62	1600-945	Page 11
KM172-29-54	N/A	N/A	1601-32-000	N/A	1600-954NS	Page 11
KM172-29-54S	5300-9-54	N/A	N/A	N/A	1600-954	Page 11
KM172-29-62	1346-9-4 3713-3-090	N/A	1649-10-000	00-0186-078-00 186-78	1600-962NS	Page 11
KM172-29-62S	5300-9-62	N/A	N/A	47-0186-078-00 47-186-78	1600-962	Page 11
KM172-36-28	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-36-35	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-36-45	N/A	N/A	N/A	00-0186-027-12	N/A	Page 11
KM172-36-62	N/A	N/A	N/A	00-0186-003-69 00-0186-027-16	N/A	Page 11
KM172-39-28	N/A	N/A	1601-36-000	N/A	N/A	Page 11
KM172-39-35	3813-4-090	N/A	1603-14-000	00-0186-001-69 26-352518	N/A	Page 11
KM172-39-45	3813-5-090	N/A	1605-14-000	00-0186-002-69 26-352548	N/A	Page 11
KM172-39-54	N/A	N/A	1607-14-000	N/A	N/A	Page 11
KM172-39-62	3813-6-090	N/A	1609-14-000	186-03-69 26-351578	N/A	Page 11
KM172-44-28	3715-5-040	N/A	1650-20-000	N/A	N/A	Page 11
KM172-44-28ESS	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-44-28S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-44-35	3715-6-040	N/A	1650-21-000	00-0186-054-00 186-54	N/A	Page 11
KM172-44-35ESS	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-44-35S	N/A	N/A	N/A	47-0186-054-00 47-186-54	N/A	Page 11
KM172-44-45	3715-7-040	N/A	1650-22-000	00-0186-070-00 186-70	N/A	Page 11
KM172-44-45ESS	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-44-45S	N/A	N/A	N/A	47-0186-070-00 47-186-70	N/A	Page 11
KM172-44-54	N/A	N/A	1650-23-000	N/A	N/A	Page 11
KM172-44-54S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-44-62	3715-8-040	N/A	1650-24-000	00-0186-086-00 186-86	N/A	Page 11
KM172-44-62ESS	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-44-62S	N/A	N/A	N/A	47-0186-086-00 47-186-86	N/A	Page 11
KM172-45-28	N/A	5050-241	N/A	N/A	N/A	Page 11
KM172-45-28S	N/A	5050-041	N/A	N/A	N/A	Page 11
KM172-45-35	N/A	5050-242	N/A	186-56	N/A	Page 11

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM172-45-35S	N/A	5050-042	N/A	N/A	N/A	Page 11
KM172-45-45	N/A	5050-243	N/A	186-72	N/A	Page 11
KM172-45-45S	N/A	5050-043	N/A	N/A	N/A	Page 11
KM172-45-62	N/A	5050-244	N/A	186-88	N/A	Page 11
KM172-45-62S	N/A	5050-044	N/A	N/A	N/A	Page 11
KM172-46-28	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-46-28S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-46-35	N/A	N/A	N/A	00-0186-058-00 186-58	N/A	Page 11
KM172-46-35S	N/A	N/A	N/A	47-0186-058-00 47-186-58	N/A	Page 11
KM172-46-45	N/A	N/A	N/A	00-0186-074-00 186-74	N/A	Page 11
KM172-46-45ESS	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-46-45S	N/A	N/A	N/A	47-0186-074-00 47-186-74	N/A	Page 11
KM172-46-54	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-46-54S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-46-62	N/A	N/A	N/A	00-0186-003-79 00-0186-090-00 186-90	N/A	Page 11
KM172-46-62ESS	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-46-62S	N/A	N/A	N/A	47-0186-090-00 47-186-90	N/A	Page 11
KM172-49-28	N/A	N/A	1601-37-000	N/A	N/A	Page 11
KM172-49-28S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-49-35	3713-4-090	N/A	1601-38-000	00-0186-001-79 186-01-79	N/A	Page 11
KM172-49-35S	N/A	N/A	N/A	47-0186-001-79 47-186-01-79	N/A	Page 11
KM172-49-45	3713-5-090	N/A	1601-39-000	00-0186-002-79 186-02-79	N/A	Page 11
KM172-49-45S	N/A	N/A	N/A	47-0186-002-79 47-186-02-79	N/A	Page 11
KM172-49-54	N/A	N/A	1601-40-000	N/A	N/A	Page 11
KM172-49-54S	N/A	N/A	N/A	N/A	N/A	Page 11
KM172-49-62	3713-6-090	N/A	1601-41-000	186-03-79	N/A	Page 11
KM172-49-62S	N/A	N/A	N/A	47-0186-003-79 47-186-03-79	N/A	Page 11
KM173-19-35	3596-1-090	N/A	1639-18-000	2620-09	N/A	Page 12
KM173-19-35S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-19-45	3596-2-090	N/A	1640-18-000	2620-07	N/A	Page 12
KM173-19-45S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-19-62	3596-3-090	N/A	1610-18	2620-05	N/A	Page 12

Cross Reference

Accessories

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM173-19-62S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-29-35	3597-1-090	N/A	1638-34-000	00-2620-039-00 2620-39	1600-935TNS	Page 12
KM173-29-35S	5500-9-35	N/A	N/A	N/A	N/A	Page 12
KM173-29-45	3597-2-090	N/A	1638-36-000	00-2620-037-00 2620-37	1600-945TNS	Page 12
KM173-29-45S	5500-9-45	N/A	N/A	N/A	N/A	Page 12
KM173-29-54	N/A	N/A	1638-38-000	N/A	1600-954TNS	Page 12
KM173-29-54S	5500-9-54	N/A	N/A	N/A	N/A	Page 12
KM173-29-62	3597-3-090	N/A	1638-40-000	00-2620-035-00 2620-35	1600-962TNS	Page 12
KM173-29-62S	5500-9-62	N/A	N/A	N/A	N/A	Page 12
KM173-39-35	3596-4-090	N/A	1638-24-000	262-10-04	N/A	Page 12
KM173-39-35S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-39-45	3596-5-090	N/A	1638-26-000	262-09-04	N/A	Page 12
KM173-39-45S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-39-54	N/A	N/A	1638-28-000	N/A	N/A	Page 12
KM173-39-54S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-39-62	3596-6-090	N/A	1610-58-000	00-0262-008-04 262-08-04	N/A	Page 12
KM173-39-62S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-49-35	3597-4-090	N/A	1638-44-000	00-0262-011-00 262-11-00	N/A	Page 12
KM173-49-35S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-49-45	3597-5-090	N/A	1638-46-000	00-0262-012-00 262-12-00	N/A	Page 12
KM173-49-45S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-49-54	N/A	N/A	1638-48-000	N/A	N/A	Page 12
KM173-49-54S	N/A	N/A	N/A	N/A	N/A	Page 12
KM173-49-62	3597-6-090	N/A	1638-50-000	00-0262-013-00 262-13-00	N/A	Page 12
KM173-49-62S	N/A	N/A	N/A	N/A	N/A	Page 12
ST6007-003-100	6007-003-100	N/A	N/A	N/A	N/A	Page 15
ST6007-003-150	6007-003-150	N/A	N/A	N/A	N/A	Page 15
ZB00-5901-020-00	N/A	N/A	N/A	00-5901-020-00	N/A	Page 16
ZB27-362093	N/A	N/A	N/A	7-362093	N/A	Page 13

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR11-2212	N/A	N/A	5627-82, 5627-82-000	11-2212 25-0901-27S-G1 25-0901-27T-G1 25-0901-27Y-F1	N/A	Page 48
BR11-2214	N/A	N/A	5627-44-00	1-2214 13-0600-89S-F1 13-0600-89T-F1 13-0600-89Y-F1	N/A	Page 48
BR11-2220	N/A	N/A	N/A	11-2220 25-0901-47S-B1 25-0901-47T-B1 25-0901-47Y-B1	N/A	Page 48
BR11-2222	N/A	N/A	N/A	11-2222 25-0901-27S-B1 25-0901-27T-B1 25-0901-27Y-B1	N/A	Page 48
BR11-2325	N/A	N/A	N/A	11-2325 13-0901-27S-F1 13-0901-27T-F1 13-0901-27Y-F1	N/A	Page 49
BR11-2608	6213-127-70 6213-127-070	N/A	5627-87	11-2608 13-0701-27S-L1 13-0701-27T-L1 13-0701-27Y-L1	N/A	Page 49
BR1-1290-35	6113-089-090 6113-89-090 2108-152-006	STK130-089-90	5627-44	11-3663 13-0900-89S-F1 13-0900-89T-F1 13-0900-89Y-F1	N/A	Page 18
BR1-1290-39	6113-097-090	STK130-100-90	N/A	11-3706 13-0901-00S-F1 13-0901-00T-F1 13-0901-00Y-F1	N/A	Page 18
BR1-1290-47	6113-119-090 2108-352-000 2108-352	STK130-119-90	5627-50 5627-11-50	11-3704 13-0901-19S-T1 13-0901-19T-F1 13-0901-19Y-F1	SO-6352	Page 18
BR1-1290-50	6113-127-090 2108-356-000 2108-356	STK130-127-90	N/A	11-3662 13-0901-27S-F1 13-0901-27T-F1 13-0901-19Y-F1 13-0901-27Y-G1	SO-6356	Page 18
BR1-1290-54	6113-137-090	STK130-137-90	N/A	13-0901-37S-F1 13-0901-37T-F1 13-0901-37Y-F1	SO-6580	Page 21
BR1-1360-50	6113-127-060	N/A	N/A	N/A	N/A	Page 19
BR1-1370-50F	6213-127-070	N/A	N/A	N/A	N/A	Page 19
BR1-1375-50	6113-127-075	N/A	N/A	N/A	N/A	Page 19
BR1-1390-42	6113-107-090	N/A	N/A	N/A	N/A	Page 20
BR1-1875-35	6118-089-075	N/A	N/A	N/A	N/A	Page 20
BR1-1875-50	6118-127-075	N/A	N/A	N/A	N/A	Page 20
BR1-19100-47	6118-119-100	STK190-119-105	N/A	19-1001-19Y-G1	N/A	Page 21
BR1-19100-50	6118-127-100 2108-385-000 2108-385	STK190-127-105	N/A	19-1001-27Y-G1	SO-685	Page 22
BR1-19100-54	6118-137-100	N/A	N/A	516061	N/A	Page 23

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR1-19110-50	6118-127-110	N/A	N/A	11-2216 11-3205 19-1101-27Y-G1	N/A	Page 22
BR1-1990-35	6118-089-090 2108-109-001	STK190-089-90	N/A	11-3453 11-3697 19-0900-89S-F1 19-0900-89T-F1 19-0900-89Y-F1	SO-6309	Page 20
BR1-1990-39	6118-097-090	STK190-100-90	N/A	11-2865 11-3696 19-0901-00S-F1 19-0901-00T-F1 19-0901-00Y-F1	N/A	Page 21
BR1-1990-47	6118-119-090	STK190-119-90	N/A	11-3110 11-3255 11-3693 11-3694 19-0901-19S-G1 19-0901-19T-G1 19-0901-19Y-G1	N/A	Page 21
BR1-1990-50	6118-127-090 2108-382-000, 2108-382	STK190-127-90	N/A	11-2210 11-3254 11-3676 11-3692 11-3965 19-0901-27S-G1 19-0901-27T-G1 19-0901-27Y-G1	SO-682	Page 22
BR1-1990-54	6118-137-090 2108-351-02	STK190-137-90	N/A	11-3690 11-3691 19-0901-37S-G1 19-0901-37T-G1 19-0901-37Y-G1	N/A	Page 22
BR1-2190-39F	6221-100-090	STK193-100-90	N/A	11-2787 11-3702 21-0901-00S-L1 21-0901-00T-L1 21-0901-00Y-L1	N/A	Page 23
BR1-2190-47F	6221-119-090 4221-120-085	STK193-119-90	5627-10-55 5627-10-651	11-2714 11-3291 11-3671 11-3701 21-0901-19S-V1 21-0901-19T-V1 21-0901-19Y-V1	N/A	Page 23
BR1-2190-50F	6221-127-090 2108-382-001 2108-382-1 2108-382-003 21080382-3	STK193-127-90	N/A	11-2721 11-3290 11-3658 11-3700 21-0901-27S-V1 21-0901-27T-V1 21-0901-27Y-V1	N/A	Page 23
BR1-2190-54F	6221-137-090 2108-382-002 2108-382-2	STK193-137-90	N/A	11-2735 11-3289 11-3672 11-3699 21-0901-37S-V1 21-0901-37T-V1 21-0901-37Y-V1	SO-6380	Page 24

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR1-25100-50	6125-127-100 2108-328-000 2108-328	N/A	N/A	N/A	SO-718	Page 26
BR1-2575-35	6125-089-075	N/A	N/A	N/A	N/A	Page 24
BR1-2590-35	6125-089-090 2108-302-000 2108-302 2108-102-003 2108-102-008	STK250-089-90	5627-43-60 5627-43-00	11-2213 11-3688 25-0900-89S-F1 25-0900-89T-F1 25-0900-89Y-F1	SO-6302	Page 24
BR1-2590-39	6125-097-090 2108-393-006 2108-393-6	STK250-100-90	N/A	11-3659 25-0901-00S-F1 25-0901-00T-F1 25-0901-00Y-F1	SO-696	Page 24
BR1-2590-39F	2108-393-005 2108-393-5	N/A	N/A	25-0901-00S-L1 25-0901-00T-L1 25-0901-00Y-L1	N/A	Page 25
BR1-2590-47	6125-119-090 2108-351-000 2108-351	STK250-119-90	5627-10-50 6527-10-60 5627-10-601	11-2586 11-3237 11-3669 11-3686 25-0901-19A-GC 25-0901-19S-B1 25-0901-19S-F1 25-0901-19S-G1 25-0901-19T-B1 25-0901-19T-F1 25-0901-19T-G1 25-0901-19Y-B1 25-0901-19Y-G1	SO-699	Page 25
BR1-2590-47F	2108-351-002 2108-351-2	N/A	5627-10-55 6527-10-65	25-0901-19S-V1 25-0901-19T-V1 25-0901-19Y-V1	N/A	Page 25
BR1-2590-50	6125-127-090 2108-376-002 2108-386-2	STK250-127-90	5627-82-60	11-3208 11-3657 11-3685 25-0901-27A-GC 25-0901-27S-F1 25-0901-27S-G1 25-0901-27T-F1 25-0901-27T-G1 25-0901-27Y-B1 25-0901-27Y-F1 25-0901-27Y-G1	SO-661	Page 25
BR1-2590-50F	2108-383-000 2108-383	N/A	N/A	25-0901-27S-V1 25-0901-27T-V1 25-0901-27Y-V1	N/A	Page 26
BR1-2590-54	6125-137-090 6125-137-90	STK250-137-90	N/A	11-3156 11-3221 11-3670 25-0901-37A-GC 25-0901-37S-B1 25-0901-37S-G1 25-0901-37T-B1 25-0901-37T-G1 25-0901-37Y-B1 25-0901-37Y-G1	SO-663	Page 26
BR1-2590-54F	N/A	N/A	N/A	25-0901-37S-V1 25-0901-37T-V1 25-0901-37Y-V1	N/A	Page 26

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR1-2590-58	6125-147-090 2108-176-002 2108-176-003	STK250-147-90	5627-11-60 5627-11-601 5627-53-00	11-3162 11-3213 11-3684 25-0901-47R-B1 25-0901-47S-B1 25-0901-47T-B1 25-0901-47Y-B1 254090147SS	N/A	Page 27
BR1-2590-58F	N/A	N/A	N/A	25-0901-47S-V1 25-0901-47T-V1 25-0901-47Y-V1 254090147SS	N/A	Page 27
BR1-3165-39F	N/A	N/A	N/A	310065099SS	N/A	Page 27
BR2108-100	2108-100-000 2108-100	5072-100	N/A	N/A	S0-600 S06-100	Page 35
BR2108-102	2108-102-000 2108-102	5072-102	N/A	N/A	S0-602 S06-102	Page 35
BR2108-103	2108-103-000 2108-103	5072-103	N/A	N/A	S0-603 S06-103	Page 35
BR2108-105	2108-105-000 2108-105	5072-105	N/A	N/A	S0-605 S06-105	Page 35
BR2108-107-1	2108-107-001	N/A	N/A	N/A	N/A	Page 36
BR2108-107-4	2108-107-004	N/A	N/A	N/A	N/A	Page 36
BR2108-107-4S1	2108-107-4S1 2108-107-4SI	N/A	N/A	N/A	N/A	Page 34 & 36
BR2108-107-9	2108-107-009 2108-107-9	N/A	N/A	N/A	N/A	Page 37
BR2108-109	2108-109-000 2108-109	5072-109	N/A	N/A	S0-609 S06-109	Page 37
BR2108-109-000S9	2108-109-000S9	N/A	N/A	N/A	N/A	Page 37
BR2108-110	2108-110-000 2108-110	5072-110	N/A	N/A	S0-610 S06-110	Page 38
BR2108-111	2108-111-000 2108-111	5072-111	N/A	N/A	S0-611 S06-111	Page 38
BR2108-113	2108-113-000 2108-113	5072-113	N/A	N/A	S0-613 S06-113	Page 38
BR2108-115	2108-115-000 2108-115	5072-115	N/A	N/A	S0-615 S06-115	Page 38
BR2108-118	2108-118-000 2108-118	5072-118	N/A	N/A	S0-618 S06-118	Page 39
BR2108-120	2108-120-000 2108-120	5072-120	N/A	N/A	S0-620 S06-120	Page 39
BR2108-120-6	2108-120-006	N/A	N/A	N/A	N/A	Page 39
BR2108-125	2108-125-000 2108-125	N/A	N/A	N/A	N/A	Page 39
BR2108-131	2108-131-000	N/A	N/A	N/A	N/A	Page 40
BR2108-140	2108-140-000 2108-140	5072-140	N/A	N/A	S0-640 S06-140	Page 40
BR2108-145	2108-145-000 2108-145	5072-145	N/A	N/A	S0-645 S06-145	Page 40

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR2108-148	2108-148-000 2108-148	5072-148	N/A	N/A	S0-648 S06-148	Page 40
BR2108-150	2108-150-000 2108-150	5072-150	N/A	N/A	S0-650 S06-150	Page 41
BR2108-151	2108-151-000 2108-151	5072-151	N/A	N/A	S0-651 S06-151	Page 41
BR2108-152	2108-152-000 2108-152	5072-152	N/A	N/A	S0-652 S06-152	Page 41
BR2108-155	2108-155-000 2108-155	5072-155	N/A	N/A	S0-655 S06-155	Page 41
BR2108-156	2108-156-000 2108-156	5072-156	N/A	N/A	S0-656 S06-156	Page 42
BR2108-158	2108-158-000 2108-158	5072-158	N/A	N/A	S0-658 S06-158	Page 42
BR2108-160	2108-160-000 2108-160	5072-160	N/A	N/A	S0-660 S06-160	Page 42
BR2108-176	2108-176-000 2108-176	5072-176	N/A	N/A	S0-676 S06-176	Page 42
BR2108-177	2108-177-000 2108-177	5072-177	N/A	N/A	S0-677 S06-177	Page 43
BR2108-182	2108-182-000 2108-182	5072-182	N/A	N/A	S0-682 S06-182	Page 43
BR2108-183	2108-183-000 2108-183	5072-183	N/A	N/A	S0-683 S06-183	Page 43
BR2108-183-1	2108-183-001	N/A	N/A	N/A	N/A	Page 43
BR2108-183-3	2108-183-003	N/A	N/A	N/A	N/A	Page 44
BR2108-185	2108-185-000 2108-185	5072-185	N/A	N/A	S0-685 S06-185	Page 44
BR2108-189	2108-189-000 2108-189	5072-189	N/A	N/A	S0-689 S06-189	Page 44
BR2108-193	2108-193-000 2108-193	5072-193	N/A	N/A	S0-693 S06-193	Page 44
BR2108-195	2108-195-000 2108-195	5072-195	N/A	N/A	S0-695 S06-195	Page 45
BR2108-197	2108-197-000 2108-197	5072-197	N/A	N/A	S0-697 S06-197	Page 45
BR2108-218	2108-218-000 2108-218	5072-218	N/A	N/A	S06-218	Page 45
BR2108-305	2108-305-000	N/A	N/A	N/A	N/A	Page 46
BR2108-309	2108-309-000	N/A	N/A	N/A	N/A	Page 46
BR2108-382-2	2108-382-002	N/A	N/A	N/A	N/A	Page 46
BR4111-119-090	4111-119-090	N/A	N/A	N/A	S011-119-090 S0611-119-090	Page 28
BR4111-127-090	4111-127-090	N/A	N/A	N/A	S011-127-090 S0611-127-090	Page 28
BR4111-135-090	4111-135-090	N/A	N/A	N/A	N/A	Page 28
BR4111-147-075	4111-147-075	N/A	N/A	N/A	S011-147-075 S0611-147-075	Page 28

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR4118-089-075	4118-089-075	N/A	N/A	N/A	S018-089-075 S0618-089-075	Page 29
BR4118-089-090	4118-089-090	N/A	N/A	N/A	S018-089-090 S0618-089-090	Page 29
BR4118-119-090	4118-119-090	N/A	N/A	N/A	S018-119-090 S0618-119-090	Page 29
BR4118-127-075	4118-127-075	N/A	N/A	N/A	N/A	Page 29
BR4118-127-090	4118-127-090	N/A	N/A	N/A	S018-127-090 S0618-127-090	Page 30
BR4118-127-100	4118-127-100	N/A	N/A	N/A	S018-127-100 S0618-127-100	Page 30
BR4118-135-090	4118-135-090	N/A	N/A	N/A	S018-135-090 S0618-135-090	Page 30
BR4118-135-105	4118-135-105	N/A	N/A	N/A	N/A	Page 30
BR4118-137-090	4118-137-090	N/A	N/A	N/A	S018-137-090 S0618-137-090	Page 31
BR4125-089-075	4125-089-075	N/A	N/A	N/A	S025-089-075 S0625-089-075	Page 31
BR4125-089-090	4125-089-090	N/A	N/A	N/A	N/A	Page 31
BR4125-097-090	4125-097-090	N/A	N/A	N/A	N/A	Page 31
BR4125-107-090	4125-107-090	N/A	N/A	N/A	N/A	Page 32
BR4125-119-090	4125-119-090	N/A	N/A	N/A	S025-119-090 S0625-119-090	Page 32
BR4125-127-090	4125-127-090	N/A	N/A	N/A	S025-127-090 S0625-127-090	Page 32
BR4125-127-100	4125-127-100	N/A	N/A	N/A	S025-127-100 S0625-127-100	Page 32
BR4125-135-090	4125-135-090	N/A	N/A	N/A	S025-135-090 S0625-135-090	Page 33
BR4125-137-090	4125-137-090	N/A	N/A	N/A	S025-137-090 S0625-137-090	Page 33
BR4125-147-090	4125-147-090	N/A	N/A	N/A	S025-147-090 S0625-147-090	Page 33
BR4135-089-064	4135-089-064	N/A	N/A	N/A	N/A	Page 34
BR5072-101	2108-109-000S15	5072-101	N/A	N/A	S0-601	Page 47
BR5072-181	2108-109-000S11	5072-181	N/A	N/A	N/A	Page 47
BR511-905S	N/A	N/A	511.905S	N/A	N/A	Page 70
BR511-907S	N/A	N/A	511.907S	N/A	N/A	Page 70
BR511-909S	N/A	N/A	511.909S	N/A	N/A	Page 70
BR511-912S	N/A	N/A	511.912S	N/A	N/A	Page 70
BR511-915S	N/A	N/A	511-915S	N/A	N/A	Page 70
K5-1290-35	N/A	N/A	05.002.201S 05.003.100S	13-0900-89A-F1	N/A	Page 63
K5-1290-39	N/A	N/A	05.002.405S 05.003.102S	13-0901-00A-F1	N/A	Page 63

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
K5-1290-47	N/A	N/A	05.002.202S 05.003.104S	13-0901-19A-F1	N/A	Page 63
K5-1290-50	N/A	N/A	05.002.203S 05.003.106S	13-0901-27A-F1	N/A	Page 63
K5-1290-54	N/A	N/A	05.002.204S 05.003.107S	13-0901-37A-F1	N/A	Page 64
K5-1990-35	N/A	N/A	05.002.001S 05.003.110S	19-0900-89A-F1	N/A	Page 64
K5-1990-39	N/A	N/A	05.002.007S 05.003.112S	19-0901-00A-F1	N/A	Page 64
K5-1990-47	N/A	N/A	05.002.002S 05.003.114S	19-0901-19A-G1	N/A	Page 64
K5-1990-50	N/A	N/A	05.002.003S 05.003.116S	19-0901-27A-F1 19-0901-27A-G1	N/A	Page 65
K5-1990-54	N/A	N/A	05.002.004S 05.003.117S	19-0901-37A-G1	N/A	Page 65
K5-2190-47F	N/A	N/A	05.002.302S 05.003.124S	21-0901-19A-V1	N/A	Page 65
K5-2190-50F	N/A	N/A	05.002.303S 05.003.126S	21-0901-27A-V1	N/A	Page 65
K5-2190-54F	N/A	N/A	05.002.304S 05.003.127S	21-0901-37A-V1	N/A	Page 66
K5-2590-35	N/A	N/A	05.002.101S 05.003.130S	25-0900-89A-F1	N/A	Page 66
K5-2590-39	N/A	N/A	05.002.106S 05.003.132S	25-0901-00A-F1	N/A	Page 66
K5-2590-47	N/A	N/A	05.002.102S 05.003.134S	25-0901-19A-B1 25-0901-19A-F1 25-0901-19A-G1	N/A	Page 67
K5-2590-50	N/A	N/A	05.002.103S 05.003.136S	25-0901-27A-B1 25-0901-27A-F1 25-0901-27A-G1	N/A	Page 67
K5-2590-58	N/A	N/A	05.002.105S 05.003.139S	25-0901-47A-B1	N/A	Page 67
K7-1290-35	N/A	5071-340 T130-089-90 TN130-089-90 V130-089-90	N/A	13-0901-89L-F1	N/A	Page 52
K7-1290-47	N/A	5071-343 T130-119-90 TN130-119-90 V130-119-90	N/A	13-0901-19L-F1	N/A	Page 52
K7-1290-50	N/A	5071-344 T130-127-90 TN130-127-90 V130-127-90	N/A	13-0901-27L-F1	N/A	Page 52
K7-1990-35	N/A	5071-300 T190-089-90 TN190-089-90 V190-089-90	N/A	19-0900-89L-F1	N/A	Page 52
K7-1990-47	N/A	5071-303 T190-119-90 TN190-119-90 V190-119-90	N/A	19-0901-19L-G1	N/A	Page 53

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
K7-1990-50	N/A	5071-304 T190-127-90 TN190-127-90 V190-127-90	N/A	19-0901-27L-F1 19-0901-27L-G1	N/A	Page 53
K7-1990-54	N/A	5071-305 T190-137-90 TN190-137-90 V190-137-90	N/A	19-0901-37L-G1	N/A	Page 53
K7-2190-47F	N/A	5071-313 T193-119-90 TN193-119-90 V193-119-90	N/A	21-0901-19L-V1	N/A	Page 53
K7-2590-35	N/A	5071-320 T250-089-90 TN250-089-90 V250-089-90	N/A	25-0900-89L-F1	N/A	Page 54
K7-2590-39	N/A	5071-321 T250-100-90 TN250-100-90 V250-100-90	N/A	25-0901-00L-F1	N/A	Page 54
K7-2590-47	N/A	5071-323 T250-119-90 TN250-119-90 V250-119-90	N/A	25-0901-19L-B1 25-0901-19L-F1 25-0901-19L-G1	N/A	Page 54
K7-2590-50	N/A	5071-324 T250-127-90 TN250-127-90 V250-127-90	N/A	25-0901-27L-B1 25-0901-27L-F1 25-0901-27L-G1	N/A	Page 54
K7-2590-54	N/A	5071-325 T250-137-90 TN250-137-90 V250-137-90	N/A	25-0901-37L-B1 25-0901-37L-F1 25-0901-37L-G1	N/A	Page 55
K7-2590-58	N/A	5071-326 T250-147-90 TN250-147-90 V250-147-90	N/A	25-0901-47L-B1 25-0901-47L-F1 25-0901-47L-G1	N/A	Page 55
KM-1201-84	N/A	N/A	N/A	N/A	N/A	Page 69
KM-2279-73	N/A	KM-479	N/A	13-0760-89S-R1	N/A	Page 51
KM-250R	277-096-250 277-96-250	5052-330 S025	5606-59	N/A	SR-250	Page 50
KM-251	5301-10-51	5052-51 5052-051 5052-251 5052-451 5052-551 Z251	5607-89 5608-16	N/A	ZO-151	Page 61
KM-256	5301-10-56	5052-56 5052-256 5052-456 5052-556 Z256	5608-17	N/A	ZO-156	Page 61
KM-275R	277-096-275 277-96-275	5052-334	N/A	N/A	SR-275	Page 50
KM-276R	277-096-276 277-96-276	N/A	N/A	12-0761-00S-RT	N/A	Page 50
KM-277R	277-096-277 277-96-277	N/A	N/A	N/A	SR-277	Page 50

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM-281R	277-096-281	N/A	N/A	N/A	N/A	Page 50
KM-300R	277-096-300 277-96-300	5052-331 S030	5606-60	N/A	SR-300	Page 51
KM-325R	277-096-325 277-96-325	5052-332 S032	5606-61	10-0701-19S-R1	SR-325	Page 51
KM-326R	277-096-326	N/A	N/A	N/A	N/A	Page 51
KM-451	N/A	5052-51 5052-051 5052-251 5052-451 5052-551	N/A	N/A	ZO-051	Page 61
KM-453	5301-10-53	5052-53 5052-053 5052-153 5052-253 5052-453 5052-553	5607-91	N/A	ZO-053	Page 61
KM-458	5301-20-58	5052-58 5052-058 5052-158 5052-258 Z258	5608-36	N/A	ZR-058	Page 62
KM-459	5301-20-59	5052-59 5052-059 5052-159 5052-259 Z259	5608-37	N/A	ZR-059	Page 62
KM-461	5301-20-61	5052-60 5052-061 5052-161 5052-261 Z261	5608-39	N/A	ZR-061	Page 62
KM-473	5301-20-73	5052-73 5052-073 5052-173 5052-273 Z273	5608-40	N/A	ZR-073	Page 62
KM-479	N/A	5052-79 5052-079 5052-179 5052-279	N/A	13-0760-89L-R1	ZR-079	Page 62
KM519-10	N/A	N/A	519.10 519.100S	N/A	N/A	Page 68
KM519-15	N/A	N/A	519.15 519.150S	N/A	N/A	Page 68
KM519-17	N/A	N/A	519.17 519.170S	N/A	N/A	Page 68
KM519-21	N/A	N/A	519.21 519.210S	N/A	N/A	Page 69
KM71-201	N/A	5071-201	N/A	N/A	N/A	Page 56
KM71-233	N/A	5071-133 5071-233	N/A	N/A	N/A	Page 56
KM71-523	N/A	5071-123 5071-423 5071-523	N/A	N/A	N/A	Page 56

Cross Reference

Large Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM71-524	N/A	5071-124 5071-424 5071-524	N/A	N/A	N/A	Page 56
KM71-545	N/A	5071-145 5071-445 5071-545	N/A	N/A	N/A	Page 57
KM71-551	N/A	5071-151 5071-451 5071-551	N/A	N/A	N/A	Page 57
KM71-552	N/A	5071-152 5071-452 5071-552	N/A	N/A	N/A	Page 57
KM71-553	N/A	5071-153 5071-453 5071-553	N/A	N/A	N/A	Page 57
KM71-556	N/A	5071-156 5071-456 5071-556	N/A	N/A	N/A	Page 58
KM71-561	N/A	5071-161 5071-461 5071-561	N/A	N/A	N/A	Page 58
KM71-563	N/A	5071-163 5071-463 5071-563	N/A	N/A	N/A	Page 58
KM71-567	N/A	5071-167 5071-467 5071-567	N/A	N/A	N/A	Page 59
KM71-572	N/A	5071-172 5071-472 5071-572	N/A	N/A	N/A	Page 59
KM71-580	N/A	5071-180 5071-480 5071-580	N/A	N/A	N/A	Page 59
KM71-581	N/A	5071-181 5071-481 5071-581	N/A	N/A	N/A	Page 59
KM71-585	N/A	5071-185 5071-485 5071-585	N/A	N/A	N/A	Page 60
KM71-586	N/A	5071-186 5071-486 5071-586	N/A	N/A	N/A	Page 60
ZB12-0761-00U-R1	N/A	N/A	N/A	12-0761-00U-R1	N/A	Page 71
ZB12-0901-19U-R1	N/A	N/A	N/A	12-0901-19U-R1	N/A	Page 71

Cross Reference

Small Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM-060	5301-35-20	5053-20 5053-020 5053-120 5053-220 Z320	5608-28 5608-41	N/A	ZR-120	Page 90
KM-062	1675-133	5053-341 S433	5606-55	N/A	SR-007	Page 82
KM-070	5301-35-21	5053-21 5053-021 5053-121 5053-221 Z321	5608-29	N/A	ZR-021	Page 90
KM-21K	5305-25-121	K121	5610-36	N/A	MS-121	Page 92
KM-22K	5301-25-122 5305-25-122	K122	5610-37	N/A	MS-122	Page 92
KM-23K	5305-25-123	K123	5610-38	N/A	MS-123	Page 92
KM-24K	5305-25-124	K124	5610-39	N/A	MS-124	Page 92
KM-3000	5301-30-62 5302-30-62	5023-62 5023-062 5023-162 5023-262 5053-62 5053-062 5053-162 5053-262 Z362	5607-98	N/A	ZO-062 ZS-362	Page 87
KM-3001	5301-30-64 5302-30-64	5023-64 5023-064 5023-164 5023-264 5053-64 5053-064 5053-164 5053-264 Z364	MO-3001	N/A	ZO-064 ZS-364	Page 87
KM-3002	5301-30-73 5302-30-73	5023-69 5023-069 5023-169 5023-269 5053-73 5053-073 5053-173 5053-273 Z373	MO-3002	N/A	ZO-073 ZS-369	Page 87
KM-3003	N/A	N/A	MO-3003	N/A	N/A	Page 87
KM-3004	N/A	Z365	5608-02	N/A	ZO-065	Page 87
KM-3005	5301-30-63 5302-30-63	5023-63 5023-063 5023-163 5023-263 5053-63 5053-063 5053-163 5053-263	MO-3005	N/A	ZO-063 ZS-363	Page 87

Cross Reference

Small Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM-3006	5301-30-60 5302-30-60	5023-60 5023-060 5023-160 5023-260 5053-60 5053-060 5053-160 5053-260 Z360	5607-96 MO-3006	N/A	ZO-060 ZS-360	Page 88
KM-3008	5301-30-66 5302-30-66	5023-66 5023-066 5023-166 5023-266 5053-66 5053-066 5053-166 5053-266 Z366	N/A	N/A	ZO-066 ZS-366	Page 88
KM-3100C	5301-25-33 5302-25-33	5023-33 5023-033 5023-133 5023-233 5053-33 5053-033 5053-133 5053-233	MS-3100C	N/A	ZS-333	Page 84
KM3-101	2296-3-101 2296-003-101	N/A	N/A	N/A	SP-101	Page 72
KM-3101	5301-25-38 5302-25-38	5023-38 5023-038 5023-138 5023-238 5053-38 5053-038 5053-138 5053-238 Z338	5608-55 MS-3101	N/A	ZS-038 ZS-338	Page 84
KM-3101C	5301-25-43 5302-25-43	5023-43 5023-043 5023-143 5023-243 5053-43 5053-043 5053-143 5053-243	MS-3101C	N/A	ZS-343	Page 84
KM3-102	2296-3-102 2296-003-102	N/A	N/A	N/A	SP-102	Page 72
KM-3102	5301-25-39 5302-25-39	5023-39 5023-039 5023-139 5023-239 5053-39 5053-039 5053-139 5053-239 Z339	5608-56 MS-3102	N/A	ZS-039 ZS-339	Page 84
KM3-103	2296-3-103 2296-003-103	N/A	N/A	N/A	SP-103	Page 72
KM3-104	2296-3-104 2296-003-104	N/A	N/A	N/A	SP-104	Page 72
KM-3104	N/A	5044-145 Z332	MS-3104	N/A	ZS-032	Page 84

Cross Reference

Small Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM3-105	2296-3-105 2296-003-105	N/A	N/A	N/A	SP-105 ZS-105	Page 72
KM-3105	N/A	5023-044 5023-144 5023-244 5053-044 5053-144 5053-244	N/A	N/A	ZS-015	Page 85
KM3-106	2296-3-106 2296-003-106	N/A	N/A	N/A	SP-106	Page 73
KM-3106	5301-25-32 5302-25-32	5023-32 5023-032 5023-132 5023-232 5053-32 5053-032 5053-132 5053-232	MS-3106	N/A	ZS-032 ZS-332	Page 85
KM-3107	5301-25-41 5302-25-41	5023-41 5023-041 5023-141 5023-241 5053-41 5053-041 5053-141 5053-241 Z341	N/A	N/A	ZS-041 ZS-341	Page 85
KM3-108	2296-3-108 2296-003-108	N/A	N/A	N/A	SP-108	Page 73
KM3-111	2296-3-111 2296-003-111	N/A	N/A	N/A	SP-111	Page 73
KM3-114	2296-3-114 2296-003-114	N/A	N/A	N/A	SP-114 ZS-114	Page 73
KM3-115	2296-3-115 2296-003-115	N/A	N/A	N/A	SP-115 ZS-115	Page 73
KM3-125	2296-3-125 2296-003-125	N/A	N/A	N/A	SP-125 ZS-125	Page 73
KM-3170	5301-30-170	5023-70 5023-170 5023-270 5053-70 5053-170 5053-270 Z370	N/A	N/A	ZO-070 ZS-370	Page 89
KM-3171	5301-30-171	5023-71 5023-171 5023-271 5053-71 5053-171 5053-271 Z371	N/A	N/A	ZO-071 ZS-371	Page 89
KM-3172	5301-30-172	5023-72 5023-172 5023-272 5053-72 5053-172 5053-272 Z372	N/A	N/A	ZO-072 ZS-372	Page 89

Cross Reference

Small Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM-3176	5301-30-176	5023-76 5023-176 5023-276 5053-76 5053-176 5053-276	N/A	N/A	ZS-373	Page 89
KM-3183	N/A	5023-183	N/A	N/A	N/A	Page 85
KM-31K	5305-25-131	K131	5610-41	N/A	MS-131	Page 93
KM3-206	2296-3-206 2296-003-206	N/A	N/A	N/A	SP-206	Page 74
KM3-212	2296-3-212 2296-003-212	N/A	N/A	N/A	SP-212	Page 74
KM3-225	2296-3-225 2296-003-225	N/A	N/A	N/A	SP-225	Page 74
KM-3231	5301-25-31 5302-25-31	5023-31 5023-031 5023-131 5023-231 5053-31 5053-031 5053-131 5053-231 Z331	MS-3231	N/A	ZS-331	Page 85
KM-3237	5301-25-37 5302-25-37	5023-37 5023-037 5023-137 5023-237 5053-37 5053-037 5053-137 5053-237	N/A	N/A	ZS-037 ZS-337	Page 86
KM-3240	5301-25-40 5302-25-40	5023-40 5023-040 5023-140 5023-240 5053-40 5053-040 5053-140 5053-240 Z340	MS-3240 5608-57	N/A	ZS-340	Page 86
KM3-254	2296-3-254 2296-003-254	N/A	N/A	N/A	N/A	Page 74
KM-3261	5301-30-61 5302-30-61	5023-61 5023-061 5023-161 5023-261 5053-61 5053-061 5053-161 5053-261 Z361	N/A	N/A	ZO-061 ZS-361	Page 88
KM-3265	5301-30-65 5302-30-65	5023-65 5023-065 5023-165 5023-265 5053-65 5053-065 5053-165 5053-265	N/A	N/A	ZO-165	Page 88

Cross Reference

Small Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM33-101	2296-33-101 2296-033-101	N/A	N/A	N/A	SP-101A	Page 76
KM33-102	2296-33-102 2296-033-102	N/A	N/A	N/A	SP-102A	Page 76
KM33-103	2296-33-103 2296-033-103	N/A	N/A	N/A	SP-103A	Page 76
KM33-104	2296-33-104 2296-033-104	N/A	N/A	N/A	SP-104A	Page 76
KM33-105	2296-33-105 2296-033-105	N/A	N/A	N/A	SP-105A	Page 76
KM33-106	2296-33-106 2296-033-106	N/A	N/A	N/A	SP-106A	Page 76
KM33-108	2296-33-108 2296-033-108	N/A	N/A	N/A	SP-108A	Page 77
KM33-111	2296-33-111 2296-033-111	N/A	N/A	N/A	SP-111A	Page 77
KM33-114	2296-33-114 2296-033-114	N/A	N/A	N/A	SP-114A	Page 77
KM33-115	2296-33-115 2296-033-115	N/A	N/A	N/A	SP-115A	Page 77
KM33-125	2296-33-125 2296-033-125	N/A	N/A	N/A	SP-125A	Page 77
KM33-206	2296-33-206 2296-033-206	N/A	N/A	N/A	SP-206A	Page 77
KM33-212	2296-33-212 2296-033-212	N/A	N/A	N/A	SP-212A	Page 78
KM33-225	2296-33-225 2296-033-225	N/A	N/A	N/A	SP-225A	Page 78
KM33-410	2296-33-410 2296-033-410	N/A	N/A	N/A	SP-410A	Page 78
KM33-411	2296-33-411 2296-033-411	N/A	N/A	N/A	SP-411A	Page 78
KM33-412	2296-33-412 2296-033-412	N/A	N/A	N/A	SP-412A	Page 78
KM33-413	2296-33-413 2296-033-413	N/A	N/A	N/A	SP-413A	Page 78
KM33-414	2296-33-414 2296-033-414	N/A	N/A	N/A	SP-414A	Page 79
KM33-521	2296-33-521 2296-033-521	N/A	N/A	N/A	SP-521A	Page 79
KM-33K	5305-25-133	K133	5610-43	N/A	MS-133	Page 93
KM3-410	2296-3-410 2296-003-410	N/A	N/A	N/A	SP-410	Page 74
KM3-411	2296-3-411 2296-003-411	N/A	N/A	N/A	SP-411	Page 75
KM3-412	2296-3-412 2296-003-412	N/A	N/A	N/A	SP-412	Page 75
KM3-413	2296-3-413 2296-003-413	N/A	N/A	N/A	SP-413	Page 75

Cross Reference

Small Bone Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
KM3-414	2296-3-414 2296-003-414	N/A	N/A	N/A	SP-414	Page 75
KM-34K	5305-25-134	K134	5610-44	N/A	MS-134	Page 93
KM3-51	2296-3-51 2296-003-051	N/A	N/A	N/A	SP-051	Page 72
KM3-511	2296-3-511 2296-003-511	N/A	N/A	N/A	SP-511	Page 75
KM3-525	2296-3-525 2296-003-525	N/A	N/A	N/A	SP-525	Page 75
KM-41K	N/A	K141	N/A	N/A	N/A	Page 94
KM-43K	N/A	K143	N/A	N/A	N/A	Page 94
KM-44K	N/A	K144	N/A	N/A	N/A	Page 94
KM-48K	N/A	K148	N/A	N/A	N/A	Page 94
KM5100-137-123	5100-137-123	N/A	N/A	N/A	N/A	Page 80
KM5100-137-133	5100-137-133	N/A	N/A	N/A	N/A	Page 80
KM5100-137-233	5100-137-233	N/A	N/A	N/A	N/A	Page 80
KM5100-37-113	5100-37-113 5100-037-113	N/A	N/A	N/A	SR-51082	Page 81
KM5100-37-114	5100-37-114 5100-037-114	N/A	N/A	N/A	SR-51083	Page 81
KM5100-37-115	5100-37-115 5100-037-115	N/A	N/A	N/A	SR-51080	Page 81
KM5100-37-116	5100-37-116 5100-037-116	N/A	N/A	N/A	SR-51081	Page 81
KM-51K	N/A	K151	N/A	N/A	N/A	Page 94
KM-54	5301-35-54	5053-54 5053-054 5053-154 5053-254 Z354	5608-42	N/A	ZR-054 ZR-154	Page 90
KM-80	N/A	5053-080 5053-80	N/A	N/A	ZR-080	Page 91
KM-800	1675-116	N/A	N/A	N/A	SR-081	Page 83
KM-801	1675-115	S415	N/A	N/A	SR-080	Page 83
KM-802	1675-113	S414	N/A	N/A	SR-082	Page 83
KM-81	N/A	5053-081 5053-81	5603-31 5303-031	N/A	ZR-081	Page 91
KM-82	N/A	5053-082 5053-82	5603-32 5303-032	N/A	ZR-082	Page 91
KM-83	N/A	5053-083 5053-83	N/A	N/A	ZR-083	Page 91

Cross Reference

Sternum Blades

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
BR2108-125	2108-125 2108-125-000	5072-125	N/A	N/A	S06-125 S0-625	Page 97
KM-278	298-97-100S3	N/A	N/A	N/A	TR-032	Page 102
KM-278N	298-97-100S3	N/A	N/A	N/A	TR-032	Page 102
KM-279	5301-10-279	5059-279	5608-13	N/A	ZO-079	Page 99
KM-279PZ	N/A	N/A	N/A	N/A	N/A	Page 99
KM-31	5301-40-33	5059-31 5059-031 5059-531 5059-031-01	N/A	N/A	ZR-031	Page 98
KM-3179SN	N/A	N/A	N/A	N/A	N/A	Page 100
KM-31N	N/A	5059-31 5059-031 5059-531 5059-031-01	N/A	N/A	N/A	Page 98
KM-32	5301-40-32	5059-32 5059-032 5059-532 5059-032-01	5608-30	N/A	ZR-032	Page 98
KM-3231SN	N/A	N/A	N/A	N/A	N/A	Page 100
KM-3231SS	N/A	N/A	N/A	N/A	N/A	Page 100
KM-32N	N/A	5059-32 5059-032 5059-532 5059-032-01	N/A	N/A	ZR-032X	Page 98
KM-7010	5310-40-110	N/A	5605-03	N/A	CD-110	Page 101
KM-97	298-97-100 298-097-100	5059-235 S490	5606-58	N/A	SR-100	Page 95
KM-97N	298-97-100 298-097-100	5059-235 S490	5606-58	N/A	SR-100	Page 95
KM-98	298-97-101 298-097-101	N/A	N/A	N/A	SR-101	Page 96
KM-98N	298-97-101 298-097-101	N/A	N/A	N/A	SR-101	Page 96
ST296-97-102	296-91-102 296-097-102	N/A	N/A	N/A	N/A	Page 95
ST296-97-102N	296-91-102 296-097-102	N/A	N/A	N/A	N/A	Page 95
ST6207-97-101	6207-97-101 6207-097-101	N/A	N/A	N/A	N/A	Page 95
ST6207-97-101N	6207-97-101 6207-097-101	N/A	N/A	N/A	N/A	Page 95

Cross Reference

High Speed Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Medtronic® Midas Rex®	Page Number(s)
BR5820-010-120	5820-010-120	N/A	N/A	N/A	N/A	N/A	Page 104
BR5820-010-130	5820-010-130	N/A	N/A	N/A	N/A	N/A	Page 104
BR5820-010-140	5820-010-140	N/A	N/A	N/A	N/A	N/A	Page 104
BR5820-010-150	5820-010-150	N/A	N/A	N/A	N/A	N/A	Page 104
BR5820-010-160	5820-010-160	N/A	N/A	N/A	N/A	N/A	Page 104
BR5820-010-170	5820-010-170	N/A	N/A	N/A	N/A	N/A	Page 104
BR5820-010-230	5820-010-230	N/A	N/A	N/A	N/A	N/A	Page 105
BR5820-010-240	5820-010-240	N/A	N/A	N/A	N/A	N/A	Page 105
BR5820-010-250	5820-010-250	N/A	N/A	N/A	N/A	N/A	Page 105
BR5820-010-640	5820-010-640	N/A	N/A	N/A	N/A	N/A	Page 106
BR5820-010-650	5820-010-650	N/A	N/A	N/A	N/A	N/A	Page 107
BR5820-010-660	5820-010-660	N/A	N/A	N/A	N/A	N/A	Page 107
BR5820-012-020	5820-012-020	N/A	N/A	N/A	N/A	N/A	Page 107
BR5820-012-030	5820-012-030	N/A	N/A	N/A	N/A	N/A	Page 107
BR5820-012-040	5820-012-040	N/A	N/A	N/A	N/A	N/A	Page 107
BR5820-012-050	5820-012-050	N/A	N/A	N/A	N/A	N/A	Page 107
BR5820-012-060	5820-012-060	N/A	N/A	N/A	N/A	N/A	Page 108
BR5820-012-310	5820-012-310	N/A	N/A	N/A	N/A	N/A	Page 109
BR5820-012-315	5820-012-315	N/A	N/A	N/A	N/A	N/A	Page 109
BR5820-012-320	5820-012-320	N/A	N/A	N/A	N/A	N/A	Page 109
BR5820-013-020	5820-013-020	N/A	N/A	N/A	N/A	N/A	Page 108
BR5820-013-030	5820-013-030	N/A	N/A	N/A	N/A	N/A	Page 108
BR5820-013-040	5820-013-040	N/A	N/A	N/A	N/A	N/A	Page 108
BR5820-013-050	5820-013-050	N/A	N/A	N/A	N/A	N/A	Page 108
BR5820-013-140	5820-013-140	N/A	N/A	N/A	N/A	N/A	Page 108
BR5820-013-150	5820-013-150	N/A	N/A	N/A	N/A	N/A	Page 109
BR5820-015-040	5820-015-040	N/A	N/A	N/A	N/A	N/A	Page 109
BR5820-015-050	5820-015-050	N/A	N/A	N/A	N/A	N/A	Page 110
BR5820-066-012	5820-066-012	N/A	N/A	N/A	N/A	N/A	Page 111
BR5820-066-016	5820-066-016	N/A	N/A	N/A	N/A	N/A	Page 111
BR5820-066-021	5820-066-021	N/A	N/A	N/A	N/A	N/A	Page 111
BR5820-066-216	5820-066-216	N/A	N/A	N/A	N/A	N/A	Page 111
BR5820-070-011	5820-070-011	N/A	N/A	N/A	N/A	N/A	Page 110
BR5820-080-021	5820-080-021	N/A	N/A	N/A	N/A	N/A	Page 110
BR5820-103-010	5820-103-010	N/A	N/A	N/A	N/A	N/A	Page 110
BR5820-103-015	5820-103-015	N/A	N/A	N/A	N/A	N/A	Page 110
BR5820-103-020	5820-103-020	N/A	N/A	N/A	N/A	N/A	Page 110
BR5820-107-020	5820-107-020	N/A	N/A	N/A	N/A	N/A	Page 105

Cross Reference

High Speed Burs

BUSA®	Stryker®	ConMed® Lintratec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Medtronic® Midas Rex®	Page Number(s)
BR5820-107-025	5820-107-025	N/A	N/A	N/A	N/A	N/A	Page 105
BR5820-107-030	5820-107-030	N/A	N/A	N/A	N/A	N/A	Page 105
BR5820-107-040	5820-107-040	N/A	N/A	N/A	N/A	N/A	Page 106
BR5820-107-130	5820-107-130	N/A	N/A	N/A	N/A	N/A	Page 109
BR5820-107-430	5820-107-430	N/A	N/A	N/A	N/A	N/A	Page 106
BR5820-107-430S1	5820-107-430S1	N/A	N/A	N/A	N/A	N/A	Page 106
BR5820-107-530	5820-107-530	N/A	N/A	N/A	N/A	N/A	Page 106
BR5820-107-530C	5820-107-530C	N/A	N/A	N/A	N/A	N/A	Page 106
KM5110-15-40	5110-15-40	N/A	N/A	N/A	N/A	N/A	Page 104
KM5120-10-230	5120-10-230 5220-10-230	N/A	N/A	N/A	N/A	N/A	Page 104
KM5120-10-240	5120-10-240 5220-10-240	N/A	N/A	N/A	N/A	N/A	Page 104
KM5130-10-240	5130-10-240 5230-10-240	N/A	N/A	N/A	N/A	N/A	Page 105
KM5130-10-250	5130-10-250 5230-10-250	N/A	N/A	N/A	N/A	N/A	Page 105
KM-AM32U	5130-10-75 5230-10-75	7021-331	N/A	N/A	N/A	AM-32	Page 106
KM-AM330U	5130-10-30 5230-10-30	7021-326	N/A	N/A	N/A	AM-330	Page 105
KM-AM35U	5130-10-50 5230-10-50 5430-10-50	7021-328	N/A	N/A	MRAM-35	AM-35	Page 105
KM-AM8U	5130-107-30 5230-107-30 5430-107-30	7021-337	5613-12	N/A	MRAM-8	AM-8	Page 103
KM-M330U	5120-10-30 5220-10-30	7021-226	N/A	N/A	N/A	M-330	Page 104
KM-M33U	5120-10-60 5220-10-60	7021-229	N/A	N/A	N/A	M-33	Page 105
KM-M340DU	5120-12-40	7021-264	N/A	N/A	N/A	M-340D	Page 104
KM-M340U	5120-10-40 5220-10-40	7021-227	N/A	N/A	N/A	M-340	Page 104
KM-M35U	5120-10-50 5220-10-50 5420-10-50	7021-228	5613-29	N/A	MRM-35	M-35	Page 105
KM-M51U	5120-130-250	7021-292	N/A	N/A	N/A	M-51	Page 103
KM-M52U	5120-80-30 5420-80-30	7021-294	N/A	N/A	N/A	M-52	Page 103
KM-M8U	5120-107-30 5220-107-30 5420-107-30	7021-237	5613-23	N/A	MRM-8	M-8	Page 103
KM-U295	5120-15-140	7021-295	N/A	N/A	N/A	N/A	Page 103

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
03-C0312	5300-10-925	5091-54 5091-154 5091-254 M40S	5607-28	ZB-154	N/A	Page 125
03-C0901	5300-10-901	5091-36 5091-136 5091-236 M40B	5607-25	N/A	ZB-136	Page 119
03-C0905	5300-10-903	5091-38 5091-138 5091-238 M55B	5607-26	N/A	ZB-138	Page 119
03-D0135	5300-10-307	5091-264 5091-164 5091-264 M40D	5607-30	N/A	ZB-164	Page 126
04-C0101	N/A	5091-20 5091-120 5091-220	5607-17	N/A	ZB-121	Page 120
04-C0103	5300-10-503	5091-22 5091-122 5091-222 M15R	5607-18	N/A	ZB-122	Page 120
04-C0104	5300-10-504	5091-21 5091-121 5091-221 M18R	5607-06	N/A	ZB-75	Page 120
04-C0105	5300-10-505	5091-24 5091-124 5091-224 M20R	5607-19	N/A	ZB-124	Page 120
04-C0106	5300-10-506 5300-10-507	5091-23 5091-123 5091-223 M23R 5091-25 5091-125 5091-225 M24R	5607-01 5607-20	N/A	ZB-070 ZB-125 ZB-70	Page 120
04-C0107	5300-10-509	5091-26 5091-126 5091-226 M30R	5607-21	N/A	ZB-126	Page 120
04-C0108	5300-10-511	5091-28 5091-128 5091-228 M40R	5605-07 5607-22	N/A	ZB-128	Page 121
04-C0109	5300-10-513	5091-30 5091-130 5091-230 M50R	5607-23	N/A	ZB-130	Page 121
04-C0110	N/A	M60R	N/A	N/A	N/A	Page 121
04-C0170	5300-10-516	5091-32 5091-132 5091-232 M65R	5607-24	N/A	ZB-132	Page 121
04-C0302	5300-10-210	5091-01 5091-101 5091-201 M10S	5607-10	N/A	ZB-101	Page 125

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
04-C0306	5300-10-214	5091-04 5091-104 5091-204 M15S	5607-13	N/A	ZB-104	Page 125
04-C0307	5300-10-216	5091-17 5091-117 5091-217 M17S	5607-02	N/A	ZB-071	Page 125
04-C0601	N/A	5091-15 5091-115 5091-215	N/A	N/A	N/A	Page 123
04-C0606	N/A	M12T	N/A	N/A	ZB-72	Page 123
04-C0607	5300-10-710	5091-18 5091-118 5091-218	N/A	N/A	ZB-73	Page 124
04-C0609	5300-10-708	5091-02 5091-102 5091-202 M15T	5607-11	N/A	ZB-102	Page 124
04-C0612	5300-10-719	5091-13 5091-113 5091-213 M18T	5607-05	N/A	ZB-74	Page 124
04-C0615	5300-10-718	N/A	N/A	N/A	N/A	Page 124
04-C0617	5300-10-714	5091-06 5091-106 5091-206 M27T	5607-14	N/A	ZB-106	Page 124
04-C0623	5300-10-704	N/A	N/A	N/A	N/A	Page 123
04-D0133	N/A	5091-62 5091-162 5091-262	N/A	N/A	N/A	Page 126
04-D0139	5300-10-303	5091-60 5091-160 5091-260 M20D	5607-29	N/A	ZB-160	Page 126
04-S1501	N/A	N/A	5607-39	N/A	ZB-192	Page 128
04-S1701	5300-10-59	5091-47 5091-147 5091-247 M10W	5606-66	N/A	ZB-147	Page 129
04-S1703	5300-10-63	5091-48 5091-148 5091-248 M15W	5606-67	N/A	ZB-148	Page 129
04-S1704	5300-10-71	5091-49 5091-149 5091-249 M20W	5608-69	N/A	ZB-149	Page 129
06-C0619	5300-10-716	5091-14 5091-114 5091-214 M20T	5607-16	N/A	ZB-114	Page 124
06-S1601	5300-10-131	5091-42 5091-142 5091-242 M20F	5608-62	N/A	ZB-142	Page 128

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
07-S1602	5300-10-133	5091-43 5091-143 5091-243 M30F	5608-63	N/A	ZB-143	Page 128
07-S1605	5300-10-135	5091-44 5091-144 5091-244 M40F	5608-64	N/A	ZB-144	Page 128
08-S1802	5300-10-41	5091-12 5091-112 5091-212 L17L	5607-15 5607-40	N/A	ZB-112 ZB-200	Page 127
09-C0109	5300-20-513	5092-30 5092-130 5092-230 L50R	5607-50	N/A	ZB-230	Page 121
09-C0110	5300-20-516	L60R	5607-51	N/A	N/A	Page 121
09-C0170	N/A	5092-32 5092-132 5092-232 L65R	N/A	N/A	ZB-232	Page 122
09-C0901	5300-20-901	5092-36 5092-136 5092-236 L40B	5607-52	N/A	ZB-236	Page 119
09-C0905	5300-20-903	5092-38 5092-138 5092-238 L55B	5607-53	N/A	ZB-238	Page 119
09-D0135	5300-20-307	5092-64 5092-164 5092-264 L40D L60D	5607-58	N/A	ZB-264	Page 126
09-D0138	5300-20-311	5092-68 5092-168 5092-268	5607-60	N/A	ZB-268	Page 127
09-D0141	5300-20-309	5092-66 5092-166 5092-266 L50D	5607-59	N/A	ZB-266	Page 127
09-S0115	N/A	N/A	N/A	N/A	N/A	Page 127
09-S0120	N/A	N/A	N/A	N/A	N/A	Page 128
10-C0101	5300-20-501	5092-20 5092-120 5092-220 L10R	5607-45	N/A	ZB-220	Page 122
10-C0103	5300-20-503	5092-22 5092-122 5092-222 L15R	5607-46	N/A	ZB-222	Page 122
10-C0105	5300-20-505	5092-24 5092-124 5092-224 L20R	5607-47	N/A	ZB-224	Page 122

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
10-C0107	5300-20-509	5092-26 5092-126 5092-226 L30R	5607-48	N/A	ZB-226	Page 122
10-C0108	5300-20-511	5092-28 5092-128 5092-228 L40R	5607-49	N/A	ZB-228	Page 122
10-C0306	5300-20-214	5092-04 5092-104 5092-204 L13S L17S	5607-42	N/A	ZB-204	Page 126
10-C0310	5300-20-921	5092-52 5092-152 5092-252 L20S	5607-54	N/A	ZB-252	Page 126
10-C0617	5300-20-714	5092-06 5092-106 5092-206 L27T	5607-43	N/A	ZB-206	Page 125
10-D0133	5300-20-305	5092-62 5092-162 5092-262 L30D	5607-57	N/A	ZB-262	Page 127
10-D0139	5300-20-303	5092-60 5092-160 5092-260 L20D	5607-56	N/A	ZB-260	Page 127
11-C0619	5300-20-716	5092-14 5092-114 5092-214 L20T	5607-44	N/A	ZB-214	Page 125
12-C0105	5300-30-505	5093-24 5093-124 5093-224 X20R	5607-64	N/A	ZB-324	Page 123
12-C0108	5300-30-511	5093-28 5093-128 5093-228 X40R	5607-66	N/A	ZB-328	Page 123
12-C0109	5300-30-513	5093-30 5093-130 5093-230 X50R	5607-67	N/A	ZB-330	Page 123
12-C0901	5300-30-901	5093-36 5093-136 5093-236 X40B	5607-68	N/A	ZB-336	Page 119
12-C0905	5300-30-903	5093-38 5093-138 5093-238 X55B	5607-69	N/A	ZB-338	Page 119
12-S0153	N/A	N/A	N/A	N/A	N/A	Page 128
14-C0104	277-10-418 277-010-418	N/A	5605-79	N/A	K-75	Page 107
14-C0105	277-10-225 277-010-225	N/A	N/A	N/A	N/A	Page 107

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
14-C0106	277-10-423 277-010-423	N/A	5605-74	N/A	K-70	Page 107
14-C0107	277-10-230 277-010-230	N/A	N/A	N/A	N/A	Page 108
14-C0108	277-10-240 277-010-240	N/A	N/A	N/A	N/A	Page 108
14-C0109	277-10-250 277-010-250	N/A	N/A	N/A	N/A	Page 108
14-C0110	N/A	N/A	N/A	N/A	N/A	Page 108
14-C0308	1607-2-103 1607-002-103	N/A	5605-75	N/A	K-71	Page 109
14-C0405	N/A	N/A	N/A	N/A		Page 109
14-C0606	277-10-212 277-010-212	N/A	5605-76	N/A	K-72	Page 109
14-C0608	277-10-216 277-010-216	N/A	5605-77	N/A	K-73	Page 109
14-C0616	277-10-221 277-010-221	N/A	5605-78	N/A	K-74	Page 109
14-C0901	277-10-62 277-010-062	N/A	N/A	N/A	N/A	Page 107
14-C0905	N/A	N/A	N/A	N/A	N/A	Page 107
14-S0201	1607-2-57 1607-002-057	N/A	N/A	N/A	K-29	Page 111
14-S0402	1607-2-1 1607-002-001	N/A	5605-15	N/A	K-01	Page 110
14-S0903	1607-2-35 1607-002-035	N/A	5605-16	N/A	K-02	Page 112
14-S0909	1607-2-3 1607-002-003	N/A	5605-17	N/A	K-03	Page 112
14-S1201	1607-2-5 1607-002-005	N/A	5605-18	N/A	K-04	Page 113
15-S1101	1607-2-47 1607-002-047	N/A	5605-21	N/A	K-08	Page 111
16-C0107	N/A	N/A	N/A	N/A	N/A	Page 108
16-C0108	N/A	N/A	N/A	N/A	N/A	Page 108
16-C0901	N/A	N/A	N/A	N/A	N/A	Page 107
16-D0130	1608-6-87 1608-006-087	N/A	5605-52	N/A	K-42	Page 109
16-D0131	1608-6-89 1608-006-089	N/A	5605-53	N/A	K-43	Page 110
16-D0132	1608-6-91 1608-006-091	N/A	5605-54	N/A	K-44	Page 110
16-D0134	1608-6-93 1608-006-093	N/A	5605-55	N/A	K-45	Page 110
16-D0135	1608-6-95 1608-006-095	N/A	5605-55	N/A	K-45	Page 110
16-D0137	1608-6-99 1608-006-099	N/A	5605-58	N/A	K-48	Page 110

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
16-S0106	1608-6-147 1608-006-147	N/A	5605-59	N/A	K-50	Page 113
16-S0107	1608-6-149 1608-006-149	N/A	5605-60	N/A	K-51	Page 113
16-S0108	1608-6-137 1608-006-137	N/A	5605-62	N/A	K-53	Page 113
16-S0109	1608-6-141 1608-006-141	N/A	5605-65	N/A	K-56	Page 113
16-S0110	1608-6-143 1608-006-143	N/A	5605-66	N/A	K-57	Page 114
16-S0115	1608-2-9 1608-002-009	N/A	5605-23	N/A	K-10	Page 114
16-S0116	1608-6-163 1608-006-163	N/A	5605-72	N/A	K-64	Page 114
16-S0117	1608-2-11 1608-002-011	N/A	5605-24	N/A	K-11	Page 114
16-S0118	1608-6-165 1608-006-165	N/A	5605-73	N/A	K-65	Page 114
16-S0119	1608-6-151 1608-006-151	N/A	5605-61	N/A	K-52	Page 114
16-S0120	1608-2-55 1608-002-055	N/A	5605-25	N/A	K-12	Page 115
16-S0121	1608-6-153 1608-006-153	N/A	5605-63	N/A	K-54	Page 115
16-S0122	1608-6-139 1608-006-139	N/A	5605-64	N/A	K-55	Page 115
16-S0123	1608-2-13 1608-002-013	N/A	5605-26	N/A	K-13	Page 115
16-S0126	1608-6-145 1608-006-145	N/A	5605-67	N/A	K-58	Page 115
16-S0402	1608-2-1 1608-002-001	N/A	5605-44	N/A	K-33	Page 111
16-S0903	1608-2-35 1608-002-035	N/A	5605-49	N/A	K-38	Page 112
16-S0909	1608-2-3 1608-002-003	N/A	5605-45	N/A	K-34	Page 112
16-S1201	1608-2-5 1608-002-005	N/A	5605-46	N/A	K-35	Page 113
16-S1401	1608-2-101 1608-002-101	N/A	5605-43	N/A	K-32	Page 115
16-S1601	1608-2-131 1608-002-131	N/A	5605-39	N/A	K-26	Page 116
16-S1801	1608-2-21 1608-002-021	N/A	5605-34	N/A	K-21	Page 111
16-S1805	N/A	N/A	N/A	N/A	N/A	Page 111
17-S1602	1608-2-133 1608-002-133	N/A	5605-40	N/A	K-27	Page 116
17-S1603	1608-2-135 1608-002-135	N/A	5605-41	N/A	K-28	Page 116

Cross Reference

Micro Burs

BUSA®	Stryker®	ConMed® Linvatec®	DePuy Synthes®	Zimmer Biomet®	MicroAire®	Page Number(s)
18-S1701	1608-2-63 1608-002-063	N/A	5605-29	N/A	K-16	Page 116
18-S1703	1608-2-59 1608-002-059	N/A	5605-31	N/A	K-18	Page 116
19-S1705	1608-2-61 1608-002-061	N/A	5605-32	N/A	K-19	Page 116
19-S1802	1608-2-43 1608-002-043	N/A	5605-35	N/A	K-22	Page 111
22-S1808	N/A	N/A	N/A	N/A	N/A	Page 112
23-S1809	N/A	N/A	N/A	N/A	N/A	Page 112
KM94901N	2296-100-26 2296-100-026	N/A	N/A	N/A	N/A	Page 118
KMS94701N	2296-100-212	N/A	N/A	N/A	N/A	Page 117
KMS94702N	2296-100-216	N/A	N/A	N/A	N/A	Page 117
KMS94703LN	2296-100-321	N/A	N/A	N/A	N/A	Page 118
KMS94703N	2296-100-221	N/A	N/A	N/A	N/A	Page 117
KMS941703LN	2296-100-421	N/A	N/A	N/A	N/A	Page 118
KMS951703LN	2296-101-421	N/A	N/A	N/A	N/A	Page 118
KMS95701N	2296-101-212	N/A	N/A	N/A	N/A	Page 117
KMS95702LN	2296-101-316	N/A	N/A	N/A	N/A	Page 118
KMS95702N	2296-101-216	N/A	N/A	N/A	N/A	Page 117
KMS958N	2296-101-423	N/A	N/A	N/A	N/A	Page 117
KMS95901N	2296-101-26 2296-101-026	N/A	N/A	N/A	N/A	Page 118

Notes

[illegible]

Notes

MANAGE COSTS WITHOUT SACRIFICING QUALITY

- ✓ Compatible With Major OEM Power Systems
- ✓ Seamless & Easy Conversion Process
- ✓ Over 35 Years of Sales and Manufacturing Experience
- ✓ No Minimum Order Quantity
- ✓ Extensive Cross Referencing System by OEM
- ✓ Free Evaluation Samples

Visit our website at BUSAMedical.com

To order contact your distributor or email productinfo@BUSAMedical.com

